## A-400 GGT

Nr kat. **7-424** (PL)

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### LITERATURA

- Szasz G., Weimann G., Suhler F., Wahlefrld A.W., Persijn J.P.: Z. Klin. Chem. Klin. Biochem., 12 (1974): 228.
- Persijn J.P., van der Silk W.: J. Clin. Chem. Clin. Biochem. 14 (1976): 421-427
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, (1994), 850-1.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA (1995), WB Saunders: 286
- Kaplan L.A., Pesce A.J.: Clinical Chemistry. Theory, analysis and correlation 3<sup>rd</sup> Ed. (1996), the C.V. Mosby Company, St. Louis: 1072.
- Shaw LM, Strømme JH, London JL, Theodorsen L. International Federation Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section: IFCC method for gamma-glutamyltransferase J Clin Chem Clin Biochem. Vol 21(1983): 633-46
- Schumann G et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of gamma-glutamyltransferase. Clin Chem Lab Med. (2002), 40(7): 734-8.
- Schumann G, Klauke R. New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clin Chim Acta 327 (2003): 69-79

Data wydania: 10. 2020.

## A-400 GGT

Cat. No **7-424** (EN)

### INTENDED USE

Diagnostic kit for determination of  $\gamma$ -glutamyltransferase activity intended to use in automatic analyzers: BS-400 and BS-480.


The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

$\gamma$ -Glutamyltransferase (GGT, GGTP) is a membrane localized enzyme that catalyzes the transfer of glutamyl groups from glutathione to amino acids or peptides. Large GGT amounts are present in secretory organs: kidney, liver, bile duct, pancreas. Although the GGT activity is highest in renal tissue, serum GGT is generally elevated as a result of liver disease. Since alcohol induces GGT production, measurement of GGT activity is used for monitoring of abstinence in withdrawal treatment.

### METHOD PRINCIPLE

IFCC method. Kinetic method with L- $\gamma$ -glutamyl-3-carboxy-4-nitroanilide.

The rate of 5-amino-2-nitrobenzoate creation measured colorimetrically is directly proportional to  $\gamma$ -glutamyltransferase activity.

### REAGENTS

**Package**  
 1-Reagent 2 x 32 ml  
 2-Reagent 1 x 17.5 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 11 weeks.

### Concentrations in the test

Tris (pH 8.25) 100 mmol/l  
 glycylglycine 100 mmol/l  
 L- $\gamma$ -glutamyl-3-carboxy-4-nitroanilide 4 mmol/l

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- 1-Reagent and 2-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Ingredients:

1-REAGENT and 2-REAGENT contain reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-2H-isothiazol-3-one, mixture (3:1).

### Warning

H317 May cause an allergic skin reaction.  
 P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

P333+P313 If skin irritation or rash occurs: Get medical advice.

P363 Wash contaminated clothing before reuse.

### SPECIMEN

Serum, EDTA plasma free from haemolysis.

Do not use citrate, oxalate and fluoride as anticoagulants because of GGT activity inhibition!

Heparin causes turbidity in the reaction mixture!

GGT activity remains stable in specimen up to 2 days at 15-25°C or 1 week at 2-8°C or 1 month at -25°C. Freezing of sample causes a loss of enzyme activity. Frozen specimens should be thawed and kept at room temperature for 18 to 24 hours before measurement to achieve full enzyme reactivation.

Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

Deionised water is recommended as a reagent blank.

### Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: GGT – URINE PROTEINS II GEN. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_002\_BS-400\_CARRYOVER.

### REFERENCE VALUES <sup>7</sup>

serum / plasma	37°C
female	< 38 U/l (0.633 $\mu$ kat/l)
male	< 55 U/l (0.917 $\mu$ kat/l)

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionized water should be used as a calibrator 0.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

### Sensitivity

8 U/l (0.133  $\mu$ kat/l) – BS-400  
 2 U/l (0.033  $\mu$ kat/l) – BS-480

### Linearity

up to 1000 U/l (16.667  $\mu$ kat/l) – BS-400  
 up to 980 U/l (16.333  $\mu$ kat/l) – BS-480

### Specificity / Interferences

Haemoglobin up to 0.16 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
<b>BS-400</b> (n = 20)	level 1	19.48	0.31	1.60
	level 2	210.52	1.18	0.56
<b>BS-480</b> (n = 10)	level 1	52.94	1.13	2.14
	level 2	141.67	1.41	1.00
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
<b>BS-400</b> (n = 14)	level 1	54.34	0.78	1.43
	level 2	159.46	1.33	0.83
<b>BS-480</b> (n = 10)	level 1	51.19	0.99	1.93
	level 2	139.72	1.04	0.75

### Method comparison

A comparison between GGT values determined at **BS-400** (y) and at **Advia 1650** (x) using 80 samples gave following results:  
 $y = 1.0607x - 7.0066$  U/l;  
 $R = 0.996$  (R – correlation coefficient)

A comparison between GGT values determined at **BS-480** (y) and at **Cobas Integra 400 Plus** (x) using 41 samples gave following results:  
 $y = 0.9867x + 0.6235$  U/l;  
 $R = 1.000$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Szasz G., Weimann G., Suhler F., Wahlefeld A.W., Persijn J.P.: Z. Klin. Chem. Klin. Biochem., 12 (1974): 228.
- Persijn J.P., van der Silk W.: J. Clin. Chem. Clin. Biochem. 14 (1976): 421-427
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, (1994), 850-1.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA (1995), WB Saunders: 286
- Kaplan L.A., Pesce A.J.: Clinical Chemistry. Theory, analysis and correlation 3<sup>rd</sup> Ed. (1996), the C.V. Mosby Company, St. Louis: 1072.
- Shaw LM, Strömme JH, London JL, Theodorsen L. International Federation Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section: IFCC method for gamma-glutamyltransferase J Clin Chem Clin Biochem. Vol 21 (1983): 633-46
- Schumann G et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of gamma-glutamyltransferase. Clin Chem Lab Med. (2002), 40(7): 734-8.
- Schumann G, Klauke R. New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clin Chim Acta 327 (2003): 69-79.

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- Szasz G., Weimann G., Suhler F., Wahlefrld A.W., Persijn J.P.: Z. Klin. Chem. Klin. Biochem., 12 (1974): 228.
- Persijn J.P., van der Silk W.: J. Clin. Chem. Clin. Biochem. 14 (1976): 421-427
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 2nd ed. Philadelphia, PA: WB Saunders, (1994), 850-1.
- Tietz N.W., ed. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia, PA (1995), WB Saunders: 286
- Kaplan L.A., Pesce A.L.: Clinical Chemistry. Theory, analysis and correlation 3<sup>rd</sup> Ed. (1996), the C.V. Mosby Company, St. Louis: 1072.
- Shaw LM, Strømme JH, London JL, Theodorsen L. International Federation Clinical Chemistry, (IFCC), Scientific Committee, Analytical Section: IFCC method for gamma-glutamyltransferase J Clin ChemClin Biochem. Vol 21(1983): 633-46
- Schumann GJ et al. International Federation of Clinical Chemistry and Laboratory Medicine. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 degrees C. International Federation of Clinical Chemistry and Laboratory Medicine. Part 6. Reference procedure for the measurement of gamma-glutamyltransferase. Clin Chem Lab Med. (2002), 40(7): 734-8.
- Schumann GJ, Klauke R. New IFCC reference procedures for the determination of catalytic activity concentrations of five enzymes in serum: preliminary upper reference limits obtained in hospitalized subjects. Clin Chim Acta 327 (2003): 69-79

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**BS-480**

Chem	GGT	No.	023	Sample Type	SERUM
Chemistry	GGT	Print name		GGT	
Reaction Type		Reaction Direction		Increase	
Pri Wave		Sec Wave		700	
Unit		Decimal		0.1	
Blank Time		Reaction Time		58 78	
Standard		Reagent Vol		Diluent	
Decreased		R1 160 μL		μL	
Increased		R2 40 μL		μL	
		R3		μL	
		R4		μL	
Linearity Range (Standard)		Linearity Limit		0.2	
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs		-33000 33000	
R1 Blank Abs		Uncapping Time		77 Day(s)	
Blank Response		Reagent Alarm Limit			
Twin Chemistry		Enzyme Linear Extension			
Prozone Check		Rate Check		Antigen Addition	
Q1 0		Q2 0		Q3 0 Q4 0	
PC 0		ABS 0			
<b>Calibration Settings</b>		<b>Auto Calibration</b>			
Math Model		Bottle Changed			
Factor		Lot Changed			
Replicates		Cal Time			
<b>Acceptance Limits</b>					
Cal Time		SD			
Slope Diff		Repeatability			
Sensitivity					
Deter Coeff					

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