A-400 CRP ULTRA

Nr kat. 7-425 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of c-reactive protein concentration, used in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

CRP (C-reactive protein) is an acute phase protein whose concentration is seen to increase as a result of the inflammatory process, most notably in response to pneumococcal (bacterial) infectious, histolytic disease and a variety of disease states. CRP to be used as a marker or general diagnostic indicator of infections and inflammation, in addition to serving as a monitor of patient response to therapy and surgery. Furthermore, regular measurements of CRP in infants can be a useful aid in the early diagnosis of infectious disease.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between CRP in a sample and anti-CRP antibody which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of CRP in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package

1-Reagent 1 x 31 ml
2-Reagent 1 x 31 ml

The reagents when stored at 2-10°C are stable up to expiry date printed on the package. The reagents are stable for 9 weeks on board the analyser at 2-10°C.

Concentrations in the test

suspension of latex particles sensitized with anti-CRP antibodies (rabbit) (pH 7.3) 0.20 w/v%
glycine buffer solution (pH 7.0)
preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C.
- Care should be taken not to interchange the caps of reagent bottles.
- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citrate).

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
For reagent blank 0.9% NaCl is recommended.

REFERENCE VALUES^{3,4}

serum, plasma	
adults	< 0.5 mg/dl (< 5 mg/l)
children (2 months – 15 years)	0.01 – 0.28 mg/dl (0.1 – 2.8 mg/l)
newborns (0 – 3 weeks)	0.01 – 0.41 mg/dl (0.1 – 4.1 mg/l)

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY CRP ULTRA CALIBRATORS kit (Cat. No 4-276) is recommended. As a 0 calibrator 0.9% NaCl should be used.

The calibration curve should be prepared every 9 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analysers BS-400 and/or BS-480. Results may vary if a different instrument is used.

- Sensitivity:** 0.08 mg/dl (BS-480)
- LOQ (Limit of Quantitation):** 0.03 mg/dl (BS-400)

- Linearity:** up to 32 mg/dl

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- Specificity / Interferences**

Haemoglobin up to 0.5 g/dl, bilirubin up to 30 mg/dl, triglycerides up to 500 mg/dl and RF up to 560 IU/ml do not interfere with the test.

- Precision**

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 20)	level 1	0.90	0.01	0.57
	level 2	2.57	0.01	0.25
BS-480 (n = 10)	level 1	0.91	0.02	2.51
	level 2	2.51	0.04	1.54

Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 20)	level 1	1.0	0.07	6.8
	level 2	2.7	0.06	2.1
BS-480 (n = 10)	level 1	0.88	0.02	2.01
	level 2	2.45	0.04	1.47

- Method comparison**

A comparison between CRP values determined at **BS-400** (y) and at **ADVIA 1800** (x) using 43 serum samples gave following results:

$$y = 1.0272 x - 0.1005 \text{ mg/dl};$$

$$R = 0.993 \quad (R - \text{correlation coefficient})$$

A comparison between CRP values determined at **BS-480** (y) and at **Advia 1650** (x) using 46 samples gave following results:

$$y = 0.9707 x + 0.2987 \text{ mg/dl};$$

$$R = 0.985 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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• Calibration		Judgment Criteria		Auto QC	
Calibration					
Rule	Spline	Sensitivity		Blank Abs.	
Replicate	2	Factor Diff.		Error Limit	
K		SD		Corr. Coeff.	
• QC					
Rules					
Westgard Multi-rule		Cum. Sum Check		Interval	
<input type="checkbox"/> 1-2S	<input type="checkbox"/> R-4S	<input type="checkbox"/> 1.0 - 2.7			
<input type="checkbox"/> 1-3S	<input type="checkbox"/> 4-1S	<input checked="" type="checkbox"/> 1.0 - 3.0			
<input type="checkbox"/> 2-2S	<input type="checkbox"/> 10-X	<input type="checkbox"/> 0.5 - 5.1			

• **BS-480**

Chem	CRP U	No.	038	Sample Type	SERUM
Chemistry	CRP ULTRA	Print name	CRP U		
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	570	Sec Wave	800		
Unit	mg/dl	Decimal	0.01		
Blank Time	48 49	Reaction Time	81 82		
Standard	Sample Vol 3 μL	Aspirated	Diluent	Reagent Vol	Diluent
Decreased	3 μL	20 μL	180 μL	R1 140 μL	μL
Increased	μL	μL	μL	R2 140 μL	μL
	μL	μL	μL	R3 μL	μL
	μL	μL	μL	R4 μL	μL
	Sample Blank <input type="checkbox"/>	Auto Retun	<input checked="" type="checkbox"/>		
Linearity Range (Standard)	0.08 32	Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	63 Day(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		Enzyme Linear Extension	<input type="checkbox"/>		
<input type="checkbox"/> Prozone Check	<input type="checkbox"/> Rate Check	Antigen Addition	<input checked="" type="checkbox"/>		
Q1 0	Q2 0	Q3 0	Q4 0		
PC 0	ABS 0				

Calibration Settings		Auto Calibration	
Math Model	Spline	<input type="checkbox"/> Bottle Changed	
Factor		<input type="checkbox"/> Lot Changed	
Replicates	3	<input type="checkbox"/> Cal Time	
Acceptance Limits			
Cal Time	1512 Hour	SD	
Slope Diff		Repeatability	
Sensitivity			
Deter Coeff			

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