A-400 CK-MB

Nr kat. 7-427 (PL)

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(EN)

INTENDED USE

Diagnostic kit for determination of CK-MB fraction activity, used in automatic analysers: BS-400 and BS-480.




The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatine kinase (CK) catalyzes the transfer of phosphate group between creatine phosphate and adenosine diphosphate (ADP). The product of this reaction is adenosine triphosphate (ATP) – molecular source of energy. CK is a dimer, composed of two different subunits called M and B. Three different isoenzymes formed from these subunits are found in brain and smooth muscle (BB), skeletal muscle (MM) and cardiac muscle (MM and MB). Increased CK-MB serum level is a strong marker of myocardial infarction.

METHOD PRINCIPLE

Optimized kinetic method according to International Federation of Clinical Chemistry (IFCC) with use of antibodies against CK-M fraction. Specific antibodies against CK-M inhibit the complete CK-MM activity (which is the main part of total CK activity) and the CK-M subunit of CK-MB. Only CK-B activity is measured.

The rate of absorbance changes at $\lambda=340$ nm is directly proportional to half of CK-MB activity (B subunit activity).

REAGENTS

Package

1-Reagent 2 x 25.5 ml
2-Reagent 2 x 7.5 ml

The reagents, stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 12 weeks.

Concentrations in the test

1-Reagent
imidazole buffer pH 6.7 100 mmol/l
glucose 20 mmol/l
N-acetylcysteine 20 mmol/l
magnesium acetate 10 mmol/l
EDTA 2 mmol/l
NADP 2 mmol/l
ADP 2 mmol/l
AMP 5 mmol/l
HK > 2.5 U/ml
polyclonal antibodies against CK-M; inhibiting capacity 8000 U/l

2-Reagent

diadenosinepentaphosphate 10 μmol/l
glucose-6-phosphate-dehydrogenase > 1.5 U/ml (G6P-DH)
creatine phosphate 30 mmol/l
preservatives

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- Do not use reagents past the expiry date.
- Do not interchange caps among reagents.
- Results CK-MB can be falsely high in case of prostate, kidney, ovary, breast and bladder cancer when isoenzyme CK-BB appears in the blood.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-Reagent contains imidazole.

Danger

H360: May damage fertility or the unborn child.
P201: Obtain special instructions before use.
P202: Do not handle until all safety precautions have been read and understood.
P308+P313: IF exposed or concerned: Get

medical advice/attention.

P405: Store locked up.

P501: Dispose of the contents/containers in accordance with the current legislation on waste treatment.

SPECIMEN

Serum, free from hemolysis.
CK activity is unstable and is rapidly lost during storage. Probes should be stored tightly closed and protected from light. Specimens can be stored up to 4-8 hours at 15-25°C or 1-2 days at 2-8°C or 1 month at -20°C.
Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent and 2-Reagent are ready to use.
Deionised water is recommended as a reagent blank.

Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: CK-MB – CK, CK-MB – CREATININE. To avoid this effect follow the recommendations contained in the instruction 51_03_24_002_BS-400_CARRYOVER.

REFERENCE VALUES ⁹

serum	37°C
adults	up to 24 U/l up to 0.401 μkat/l

The probability that cardiac infarction has occurred is high when CK-MB and total CK activities are above normal values and CK-MB activity is between 6 and 25% of the total CK activity.

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY CK-MB CONTROL N (Cat. No 5-183) and CORMAY CK-MB CONTROL P (Cat. No 5-184) with each batch of samples.

For the calibration the CORMAY CK-MB CALIBRATOR (Cat. No 5-182) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

Sensitivity:

11 U/l (0.18 μkat/l) – BS-400

8 U/l (0.13 μkat/l) – BS-400

Linearity

up to 2000 U/l (33.4 μkat/l) – BS-400

up to 4500 U/l (75.0 μkat/l) – BS-480

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin interfere even in small amounts, bilirubin up to 20 mg/dl, ascorbate up to 62 mg/l and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
BS-400 (n = 10)	level 1	42.6	0.53	1.25
	level 2	100.9	0.65	0.64
BS-480 (n = 10)	level 1	40.20	0.81	2.01
	level 2	126.60	0.38	0.30
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
BS-400 (n = 10)	level 1	41.90	0.41	0.99
	level 2	157.53	1.14	0.73
BS-480 (n = 10)	level 1	42.10	0.53	1.26
	level 2	153.46	1.08	0.70

Method comparison

A comparison between CK-MB values determined at **BS-400** (y) and at **Cobas Integra 400** (x) using 40 samples gave following results:

$$y = 1.084x - 2.7 \text{ U/l;}$$

$$R = 0.962 \quad (R - \text{correlation coefficient})$$

A comparison between CK-MB values determined at **BS-480** (y) and at **BS-800** (x) using 49 samples gave following results:

$$y = 0.9686x - 0.0311 \text{ U/l;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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BS-480

Chem	CK-MB	No.	022	Sample Type	SERUM
Chemistry	CK-MB	Print name	CK-MB		
Reaction Type	Kinetic	Reaction Direction	Increase		
Pri Wave	340	Sec Wave	412		
Unit	U/L	Decimal	0.1		
Blank Time	0	Reaction Time	72	82	
Standard	Sample Vol	Aspirated	Diluent	Reagent Vol	Diluent
Decreased	12 μL	20 μL	180 μL	R1 200 μL	
Increased	12 μL	20 μL	180 μL	R2 40 μL	
				R3	
				R4	
	Sample Blank	Auto Retun			
Linearity Range (Standard)	8.0	4500	Linearity Limit	0.2	
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-33000 33000	
R1 Blank Abs	-33000	33000	Uncapping Time	77 Day(s)	
Blank Response	-33000	33000	Reagent Alarm Limit		
Twin Chemistry			Enzyme Linear Extension		
	Prozone Check	Rate Check	Antigen Addition		
Q1	0	Q2	0	Q3	0
PC	0	ABS	0	Q4	0
Calibration Settings		Auto Calibration			
Math Model	Two-point Linear		Bottle Changed		
Factor		Replicates	3	Lot Changed	
				Cal Time	
Acceptance Limits					
Cal Time	1848	Hour	SD		
Slope Diff			Repeatability		
Sensitivity					
Deter Coeff					

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