A-400 MG

Nr kat. 7-429 (PL)

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Cat. No 7-429 (EN)

INTENDED USE

Diagnostic kit for determination of magnesium concentration, used in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Magnesium in human organism occurs mainly in bone (about 50%) but is present also intracellularly in other tissues. Magnesium serves as a cofactor for multiple enzymatic reactions involved in nucleic acids synthesis, transport and production of energy. Magnesium is important in neuromuscular conduction and activation. Reduced magnesium level generates: concentration disturbances, fatigue, muscle tremor, anxiety state.

METHOD PRINCIPLE

Magnesium forms a purple coloured complex in alkaline solution. In the presence of EGTA, the reaction is specific. The intensity of the purple colour is proportional to the magnesium concentration.

REAGENTS

Package

1-Reagent 2 x 32.5 ml

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. The reagents are stable for 3 weeks on board the analyser at 2-10°C.

Concentrations in the test

xylylidyl blue	≤ 0.18 mmol/l
EGTA	≤ 0.12 mmol/l
CAPS	≤ 0.06 mol/l
potassium hydroxide	≤ 0.07 mol/l
buffer (pH 11.5)	≤ 0.14 mol/l
detergent	
preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- The reagent is air sensitive, to extend reagents stability it is recommended to keep reagent's bottles recapped on the board of analyser
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-Reagent contains potassium hydroxide.

Danger.

H314 Causes severe skin burns and eye damage.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305 +P351 +P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER or doctor.

SPECIMEN ^{6,7}

Serum, heparinized plasma free from hemolysis, 24-hours urine.

Recommended anticoagulants: heparine lithium, sodium or ammonium salt.

Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain approximately 3 times the magnesium concentration found in normal serum.

Urine preparation: Acidify urine with some drops of concentrated hydrochloride acid to pH 1.0. Then dilute 1 part of acidified urine with 4 parts of distilled water. Multiply the result by 5. Mix well samples before analysis.

Serum and plasma can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C. 24-hours urine samples can be stored up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

Avoid foaming.

For reagent blank deionized water is recommended.

Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: ALP-MG, AMYLASE EPS-MG, CHOL-MG, CK-MG, GLUCOSE HEX-MG, dTIBC-MG. To avoid this effect, it is recommended to **precede the MG test by wash cycle** with A-400 ACID WASHING SOLUTION (Cat. No 3-111) according to the guidelines contained in the instruction: 51_03_24_002_BS-400_CARRYOVER.

If possible, we recommend moving the MG test to the top of the Test Order list.

REFERENCE VALUES ⁶

serum / plasma	mg/dl	mmol/l
newborn 2 – 4 d	1.5 – 2.2	0.62 – 0.91
children 5 mo – 6 y	1.7 – 2.3	0.70 – 0.95
6 – 12 y	1.7 – 2.1	0.70 – 0.86
12 – 20 y	1.7 – 2.2	0.70 – 0.91
adults	1.6 – 2.6	0.66 – 1.07

24-hours urine:	mg/24h	mmol/24h
	72.9 – 145.8	3 – 5

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the following controls: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for **determination in serum**, CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for **determination in urine**. For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177). Deionised water should be used as a calibrator 0. The calibration curve should be prepared every week, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

Sensitivity:

0.13 mg/dl (0.053 mmol/l) – BS-400

0.06 mg/dl (0.025 mmol/l) – BS-480

Linearity:

up to 5.0 mg/dl (2.05 mmol/l) – BS-400

up to 5.0 mg/dl (2.05 mmol/l) – BS-480

Specificity / Interferences

Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 15 mg/dl, triglycerides up to 1000 mg/dl and calcium up to 20 mg/dl do not interfere with the test.

Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n = 10	level 1	1.98	0.01	0.56
	level 2	4.27	0.03	0.76
BS-480 n = 10	level 1	2.00	0.00	0.24
	level 2	4.30	0.02	0.53
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n = 20	level 1	2.01	0.02	1.23
	level 2	4.20	0.07	1.67
BS-480 n = 20	level 1	2.01	0.05	2.24
	level 2	4.19	0.08	1.82

Method comparison

A comparison between magnesium values determined at **BS-400** (y) and at **ADVIA 1650** (x) using 144 samples gave following results:

$$y = 0.9311x + 0.0314 \text{ mg/dl};$$

$$R = 0.980 \quad (R - \text{correlation coefficient})$$

A comparison between magnesium values determined at **BS-480** (y) and at **ADVIA 1650** (x) using 56 samples gave following results:

$$y = 0.9342x + 0.1362 \text{ mg/dl};$$

$$R = 0.922 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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• **Calibration**

Calibration		Judgment Criteria	
Rule	Multi-point Linear	Sensitivity	
Replicate	3	Factor Diff.	
K		SD	
• QC			
Rules		Auto QC	
Westgard Multi-rule		Cum. Sum Check	
v 1-2S	v R-4S	1.0 - 2.7	Interval
v 1-3S	v 4-1S	• 1.0 - 3.0	
v 2-2S	v 10-X	0.5 - 5.1	

• **BS-480**

Chem	MG	No.	014	Sample Type	SERUM/URINE
Chemistry	Magnesium	Print name	MG		
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	505	Sec Wave	700		
Unit	mg/dl	Decimal	0.01		
Blank Time	11 12	Reaction Time	20 21		
Standard	Sample Vol 2 μL	Aspirated	μL	Diluent	μL
Decreased	2 μL	30 μL	120 μL	R1	200 μL
Increased	μL	μL	μL	R2	μL
	μL	μL	μL	R3	μL
	μL	μL	μL	R4	μL
	Sample Blank	<input type="checkbox"/>	Auto Retun		
Linearity Range (Standard)	0.06 5	Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	21 Day(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check	<input type="checkbox"/> Rate Check	• Antigen Addition			
O1 0	O2 0	O3 0	O4 0		
PC 0	ABS 0				

Calibration Settings		Auto Calibration	
Math Model	Multi-point Linear	<input type="checkbox"/> Bottle Changed	
Factor		<input type="checkbox"/> Lot Changed	
Replicates	3	<input type="checkbox"/> Cal Time	
Acceptance Limits			
Cal Time	168 Hour	SD	
Slope Diff		Repeatability	
Sensitivity			
Deter Coeff			

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