## A-400 ASO

Nr kat. **7-440** (PL)

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Cat. No **7-440** (EN)

### INTENDED USE

Diagnostic kit for determination of anti-streptolysin O levels, intended to use in automatic analyzers: BS-400 and BS-480. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Most people infected with hemolytic streptococcus produce anti-streptolysin O (ASO), antibodies against streptolysin O (SLO), an exotoxin of Streptococcus. Measuring the level of ASO is effective for diagnosing, judging the progress of medical treatment, and assessing recovery from diseases caused by hemolytic streptococcus such as rheumatic fever, acute glomerulonephritis, scarlatina and tonsillitis.

### METHOD PRINCIPLE

When an antigen-antibody reaction occurs between ASO in a sample and SLO which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of ASO in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

### REAGENTS

**Package**  
 1-Reagent 1 x 40 ml  
 2-Reagent 1 x 28.5 ml

The reagents, stored at 2-10°C are stable up to expiry date printed on the package. The reagents stored on board of the analyser at 2-10°C are stable for 11 weeks.

### Concentrations in the test

suspension of latex particles sensitized with SLO (pH 8.2) 0.17 w/v%  
 glycine buffer solution (pH 8.3)  
 preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Reagents with different lot numbers should not be interchanged or mixed.
- Do not interchange the caps of reagent bottles.
- Reagents should be mixed before use by gentle inverting the bottle several times
- Immediately after use, recap the bottles and store at 2-10°C.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

### SPECIMEN

Serum or plasma collected on heparine (sodium or lithium salt), EDTA (sodium or potassium salt) or citric acid. After the blood is completely clotted, the sample should be centrifuged and separated from the cells and fibrinogen. If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Avoid repeated freezing and thawing. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
 0.9% NaCl is recommended as a reagent blank

### REFERENCE VALUES <sup>3</sup>

serum, plasma	< 160 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples. For the calibration of automatic analysers systems the CORMAY ASO CALIBRATOR kit (Cat. No 4-278) is recommended. Calibrator and 0.9% NaCl should be used for calibration.

The calibration curve should be prepared every 11 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and/or BS-480. Results may vary if a different instrument is used.

- Sensitivity:** 18 IU/ml
- Linearity:** up to 630 IU/ml

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 20 mg/dl and triglycerides up to 500 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [IU/ml]	SD [IU/ml]	CV [%]
<b>BS-400</b> (n = 10)	level 1	149.05	2.54	1.70
	level 2	291.79	1.36	0.47
<b>BS-480</b> (n = 10)	level 1	151.66	4.92	3.24
	level 2	291.93	4.43	1.52

Reproducibility (day to day)		Mean [IU/ml]	SD [IU/ml]	CV [%]
<b>BS-400</b> (n = 10)	level 1	149.59	2.29	1.53
	level 2	293.90	4.82	1.64
<b>BS-480</b> (n = 10)	level 1	133.43	6.46	4.84
	level 2	277.08	3.69	1.33

### Method comparison

A comparison between ASO values determined at **BS-400** (y) and at **Olympus AU400** (x) using 30 samples gave following results:

$$y = 1.0527x - 2.4214 \text{ IU/ml};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

A comparison between ASO values determined at **BS-400** (y) and at **Olympus AU400** (x) using 46 samples gave following results:

$$y = 1.0527x - 2.4214 \text{ IU/ml};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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• **BS-480**

Chem	ASO	No.	033	Sample Type	SERUM
Chemistry	ASO	Print name	ASO		
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	570	Sec Wave	800		
Unit	IU/mL	Decimal	0.1		
Blank Time	48	Reaction Time	80	82	
Standard	Sample Vol 4 μL	Aspirated 49 μL	Diluent μL	Reagent Vol 160 μL	Diluent μL
Decreased	4 μL	20 μL	180 μL	R2 110 μL	
Increased				R3	
	Sample Blank	v	Auto Return	R4	
Linearity Range (Standard)	18.0	630	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-33000 33000	
R1 Blank Abs	-33000	33000	Uncapping Time	77 Day(s)	
Blank Response	-33000	33000	Reagent Alarm Limit		
Twin Chemistry			<input type="checkbox"/> Enzyme Linear Extension		
	<input type="checkbox"/> Prozone Check	<input type="checkbox"/> Rate Check	<input type="checkbox"/> Antigen Addition		
Q1	0	Q2	0	Q3	0
PC	0	ABS	0	Q4	0
<b>Calibration Settings</b>		<b>Auto Calibration</b>			
Math Model	Two-point Linear	<input type="checkbox"/> Bottle Changed			
Factor		<input type="checkbox"/> Lot Changed			
Replicates	3	<input type="checkbox"/> Cal Time			
<b>Acceptance Limits</b>					
Cal Time	1848	Hour	SD		
Slope Diff		Repeatability			
Sensitivity					
Deter Coeff					

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