## A-400 BIL TOTAL

Nr kat. **7-445** (PL)

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### LITERATURA

- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2):116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

Data wydania: 10.2020.

## A-400 BIL TOTAL

Cat. No **7-445** (EN)

### INTENDED USE

Diagnostic kit for determination of total bilirubin concentration intended to use in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions

### INTRODUCTION

Bilirubin is a yellow pigment – product of heme degradation. For clinical purposes, bilirubin is expressed as two fractions: conjugated and unconjugated. In hepatocytes bilirubin is enzymatically conjugated with glucuronic acid residues. This form is called direct or conjugated. Bilirubin without glucuronic acid modification is bound to albumin and is termed unconjugated or indirect. Indirect bilirubin is calculated as the difference between total and direct bilirubin.

Serum bilirubin measurement is widely used as a screening test for liver functions. Hiperbilirubinemia is usually the result of jaundice (mechanical, hemolytic), Dubin-Jonson syndrome, Gilbert's syndrome, Crigler-Najjar syndrome, bile ducts disease.

### METHOD PRINCIPLE

Method is based on chemical oxidation, utilizing vanadate as an oxidizing agent.

In the presence of detergent and vanadate in an acidic solution, total bilirubin (both conjugated – direct, and unconjugated bilirubin) is oxidized to produce biliverdin.

This oxidation reaction causes change of the yellow colour, which is specific to bilirubin to the green colour typical for biliverdin. Therefore, the total bilirubin concentration in the sample can be obtained by measuring the absorbance before and after the vanadate oxidation.

### REAGENTS

#### Package

1-Reagent 2 x 31 ml  
2-Reagent 1 x 17 ml

The reagents when stored at 10-25°C are stable up to expiry date printed on the package. The reagents are stable for 11 weeks on board the analyser at 2-10°C.

### Concentrations in the test

**1-Reagent**  
citrate buffer (pH 2.8) 90 mmol/l  
detergent

**2-Reagent**  
phosphate buffer (pH 7.0) 4.6 mmol/l  
sodium metavanadate 3.0 mmol/l

### Warnings and notes

- Do not freeze reagents.
- Protect from direct sunlight and avoid contamination!
- Do not use after expiry date.
- Do not interchange caps.
- Reagent bottles should be shaken before use by gently inverting several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lack of significant changes in the color of the reaction mixture at the samples with low bilirubin concentration does not indicate the assay malfunction.
- 1-Reagent meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

### Warning

H319 Causes serious eye irritation.  
P280 Wear protective gloves/protective clothing/eye protection/face protection.  
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

### SPECIMEN

Serum free from haemolysis.  
Serum should be separated from red blood cells as soon as possible after blood collection.  
Lipemic specimens may show falsely decreased bilirubin concentration thus fasting specimen is recommended. It is recommended to follow CLSI procedures regarding specimen collecting and handling.  
Because bilirubin is photooxidized when exposed to light, specimen should be protected from direct exposure to either artificial light or sunlight. Therefore it is essential to store specimens in the dark at 2-8°C, at the most 3 days. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.  
For reagent blank deionized water is recommended.

#### Actions required:

When performing assays at analyzer **BS-400**, there is a probability of **cross-contamination** affecting the tests results: CREATININE – BIL TOTAL II GEN. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_002\_BS-400\_CARRYOVER.

### REFERENCE VALUES <sup>4</sup>

Age	mg/dl	µmol/l
0 - 1 days	< 8	< 137
1 - 2 days	< 12	< 205
3 - 5 days	< 16	< 274
5 days - 60 years	0.3 - 1.2	5 - 21
60 - 90 years	0.2 - 1.1	3 - 19
> 90 years	0.2 - 0.9	3 - 15

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers systems are recommended the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177). As a 0 calibrator deionized water should be used.

The calibration curve should be prepared every 9 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity:

0.06 mg/dl (1.03 µmol/l) – BS-400  
0.07 mg/dl (1.20 µmol/l) – BS-480

#### Linearity:

up to 69 mg/dl (1180 µmol/l) – BS-400  
up to 48.7 mg/dl (832.9 µmol/l) – BS-480

#### Specificity / Interferences

Haemoglobin up to 0.25 g/dl, ascorbic acid up to 500 mg/l and intralipid up to 250 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
<b>BS-400</b> (n = 20)	level 1	1.00	0.02	1.73
	level 2	4.20	0.02	0.49
<b>BS-480</b> (n = 10)	level 1	1.17	0.01	0.49
	level 2	4.33	0.01	0.13
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
<b>BS-400</b> (n = 19)	level 1	1.03	0.08	7.26
	level 2	4.20	0.14	3.28
<b>BS-480</b> (n = 10)	level 1	1.20	0.01	1.20
	level 2	4.46	0.05	1.21

#### Method comparison

A comparison between CORMAY reagent (y) and another commercially available assay (x) using 41 samples gave following results:

$y = 0.8872x + 0.0143$  mg/dl;  
 $R = 0.999$  (R – correlation coefficient)

A comparison between total bilirubin values determined at **BS-480** (y) and at **Advia 1650** (x) using 82 samples gave following results:

$y = 0.9187x + 0.0385$  U/l;  
 $R = 0.999$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2);116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999; 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996; 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

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- Tokuda K. Tanimoto K. New method of measuring serum bilirubin using vanadic acid. Jpn J Clin. Chem. 1993;22(2);116-122.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders; 1999: 1803.
- Tietz NW. Fundamentals of Clinical Chemistry, 4th ed. Edited by Burtis CA. and Ashwood ER. WB Saunders Company; 1996: 547.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 172 (2006).

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• **Calibration**

<b>Calibration</b>	<b>Judgment Criteria</b>
Rule <input type="text" value="Multi-point Linear"/>	Sensitivity <input type="text"/>
Replicate <input type="text" value="3"/>	Factor Diff. <input type="text"/>
K <input type="text"/>	SD <input type="text"/>
<b>• QC</b>	Blank Abs. <input type="text"/>
<b>Rules</b>	Error Limit <input type="text"/>
Westgard Multi-rule	Corr. Coeff. <input type="text"/>
<input type="checkbox"/> 1-2S	<input type="checkbox"/> R-4S
<input type="checkbox"/> 1-3S	<input type="checkbox"/> 4-1S
<input type="checkbox"/> 2-2S	<input type="checkbox"/> 10-X
Cum. Sum Check <input type="text"/>	Auto QC Interval <input type="text"/>
<input type="checkbox"/> 1.0 - 2.7	
<input type="checkbox"/> 1.0 - 3.0	
<input type="checkbox"/> 0.5 - 5.1	

• **BS-480**

Chem <input type="text" value="BIL T"/>	No. <input type="text" value="005"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="Bilirubin Total"/>	Print name <input type="text" value="BIL T"/>	
Reaction Type <input type="text" value="Endpoint"/>	Reaction Direction <input type="text" value="Decrease"/>	
Pri Wave <input type="text" value="450"/>	Sec Wave <input type="text" value="546"/>	
Unit <input type="text" value="mg/dl"/>	Decimal <input type="text" value="0.01"/>	
Blank Time <input type="text" value="46"/> <input type="text" value="47"/>	Reaction Time <input type="text" value="81"/> <input type="text" value="82"/>	
Standard <input type="text" value="11"/> <input type="text" value="11"/> <input type="text" value="11"/>	Reagent Vol <input type="text" value="240"/> <input type="text" value="60"/>	Diluent <input type="text"/>
Decreased <input type="text" value="11"/> <input type="text" value="20"/> <input type="text" value="180"/>	R1 <input type="text" value="240"/> <input type="text" value="60"/>	R2 <input type="text" value="60"/>
Increased <input type="text" value="11"/> <input type="text" value="20"/> <input type="text" value="180"/>	R3 <input type="text" value="60"/>	R4 <input type="text" value="60"/>
<input type="checkbox"/> Sample Blank <input type="checkbox"/> Auto Retun		
Linearity Range (Standard) <input type="text" value="0.07"/> <input type="text" value="48.7"/>	Linearity Limit <input type="text"/>	
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>	
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	
R1 Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	Uncapping Time <input type="text" value="84"/> Day(s)	
Blank Response <input type="text" value="-33000"/> <input type="text" value="33000"/>	Reagent Alarm Limit <input type="text"/>	
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check <input type="checkbox"/> Rate Check	• Antigen Addition	
Q1 <input type="text" value="0"/> <input type="text" value="0"/>	Q3 <input type="text" value="0"/> <input type="text" value="0"/>	Q4 <input type="text" value="0"/>
PC <input type="text" value="0"/> <input type="text" value="0"/>	ABS <input type="text" value="0"/>	

<b>Calibration Settings</b>	<b>Auto Calibration</b>
Math Model <input type="text" value="Multi-point Linear"/>	<input type="checkbox"/> Bottle Changed
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed
Replicates <input type="text" value="3"/>	<input type="checkbox"/> Cal Time
<b>Acceptance Limits</b>	
Cal Time <input type="text" value="2016"/> Hour	
Slope Diff <input type="text"/>	SD <input type="text"/>
Sensitivity <input type="text"/>	Repeatability <input type="text"/>
Deter Coeff <input type="text"/>	

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