## A-400 CALCIUM ARSENAZO

Nr kat. 7-447 (PL)

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Cat. No 7-447 (EN)

### INTENDED USE

Diagnostic kit for determination of calcium concentration, used in automatic analysers: BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Calcium plays an essential role in many cell functions: intracellularly in muscle contraction and glycogen metabolism, extracellularly, in bone mineralization, in blood coagulation and in transmission of nerve impulses. Calcium is present in plasma in three forms: free, bound to proteins or complexed with anions as phosphate, citrate and bicarbonate. Decreased total calcium levels can be associated with diseases of the bone apparatus (especially osteoporosis), kidney diseases (especially under dialysis), defective intestinal absorption and hypoparathyroidism. Increased total calcium can be measured in hyperparathyroidism, malignant diseases with metastases and sarcoidosis. Calcium measurements also help in monitoring of calcium supplementation mainly in the prevention of osteoporosis.

### METHOD PRINCIPLE

Photometric test using arsenazo III.

Calcium with arsenazo III at neutral pH yields a blue colored complex, whose intensity is proportional to the calcium concentration. Interference by magnesium is eliminated by addition of 8-hydroxyquinoline-5-sulfonic acid.

### REAGENTS

#### PACKAGE

1-Reagent 4 x 35.5 ml

The reagent, stored at 2-8°C is stable up to expiry date printed on the package

### Concentrations in the test

phosphate buffer (pH 7.5) ≤ 80 mmol/l  
8-hydroxyquinoline-5-sulfonic acid ≤ 6 mmol/l  
arsenazo III ≤ 120 μmol/l  
detergents, preservative

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.

### SPECIMEN

Serum. Random or 24-hours urine.

Serum can be stored up to 7 days at 20-25°C or up to 3 weeks at 4-8°C. Samples frozen at -20°C can be stored up to 8 months. Discard contaminated specimens.

**24-hours urine preparation:** To prevent calcium salt precipitation specimens should be collected in 10 ml of 6M HCl. In case of presence of precipitants they can be solved by lowering pH of the urine to below 2.0. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent is ready to use.

Deionized water is recommended as a reagent blank.

#### Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: LIPASE II GEN – CALCIUM ARSENAZO, RF – CALCIUM ARSENAZO, CALCIUM ARSENAZO – PHOSPHORUS. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_002\_BS-400\_CARRYOVER

### REFERENCE VALUES <sup>4</sup>

serum	mg/dl	mmol/l
adult	8.6 – 10.3	2.15 – 2.57
random urine	mg/dl	mmol/l
male	0.9 – 37.9	0.225 – 9.47
female	0.5 – 35.7	0.125 – 8.92
24-hours urine	mg/24h	mmol/24h
adult	100 – 300	2.5 – 7.5

It is recommended for each laboratory to establish its own reference ranges for local population.

### Calcium concentration 24-hours urine – calculation

calcium concentration in 24-hours urine [mg/24h]	=	calcium concentration in sample of 24-hours urine [mg/dl]	×	urine volume of 24-hours urine [dl/24h]
--	---	---	---	---

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples the following controls: CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum; CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers BS-400 and/or BS-480. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**  
0.08 mg/dl (0.02 mmol/l) - BS-400, BS-480
- LoD (Limit of Detection):**  
0.13 mg/dl (0.033 mmol/l) - BS-400, BS-480
- LoQ (Limit of Quantitation):**  
0.61 mg/dl (0.153 mmol/l) - BS-400  
0.31 mg/dl (0.078 mmol/l) - BS-480
- Linearity:**  
up to 20.3 mg/dl (5.075 mmol/l) – BS-400  
up to 20.1 mg/dl (5.025 mmol/l) – BS-480

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

#### Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl and magnesium up to 20 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=20	level 1	10.1	0.12	1.18
	level 2	12.8	0.07	0.56
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 n=80	level 1	9.9	0.16	1.6
	level 2	12.8	0.14	1.1

#### Method comparison

A comparison between calcium values determined at **BS-400** (y) and at **BECKMAN COULTER AU680** (x) using 78 serum samples gave following results:

$$y = 1.0253x + 0.099 \text{ mg/dl};$$

$$R = 0.959 \quad (R - \text{correlation coefficient})$$

A comparison between calcium values determined at **BS-400** (y) and at **BECKMAN COULTER AU680** (x) using 67 urine samples gave following results:

$$y = 0.9577x + 0.4428 \text{ mg/dl};$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between calcium values determined at **BS-480** (y) and at **BECKMAN COULTER AU680** (x) using 30 serum samples gave following results:

$$y = 0.9163x + 1.1433 \text{ mg/dl};$$

$$R = 0.987 \quad (R - \text{correlation coefficient})$$

A comparison between calcium values determined at **BS-480** (y) and at **BECKMAN COULTER AU680** (x) using 33 urine samples gave following results:

$$y = 0.8907x + 1.09 \text{ mg/dl};$$

$$R = 0.993 \quad (R - \text{correlation coefficient})$$

### WASTE MANAGEMENT

Please refer to local legal requirements.

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#### BS-480

Chem Ca Ars		No. 060		Sample Type SERUM/URINE	
Chemistry Calcium Arsenazo		Print name Ca Ars			
Reaction Type Endpoint		Reaction Direction Increase			
Pri Wave 660		Sec Wave 700			
Unit mg/dL		Decimal 0.01			
Blank Time 11 12		Reaction Time 26 27			
Standard 2.5 μL Aspirated μL Diluent μL		Reagent Vol 250 μL Diluent μL			
Decreased 2.5 μL 20 μL 180 μL		R1 250 μL			
Increased μL μL μL		R2 μL			
Sample Blank v Auto Retun		R3 μL			
		R4 μL			
Linearity Range (Standard) 0.31 20.1		Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs -33000 33000			
R1 Blank Abs -33000 33000		Uncapping Time Day(s)			
Blank Response -33000 33000		Reagent Alarm Limit			
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension			
<input type="checkbox"/> Prozone Check <input type="checkbox"/> Rate Check		• Antigen Addition			
Q1 Q2 Q3 Q4					
PC ABS					
<b>Calibration Settings</b>		<b>Auto Calibration</b>			
Math Model Multi-point Linear		<input type="checkbox"/> Bottle Changed			
Factor Replicates 2		<input type="checkbox"/> Lot Changed			
		<input type="checkbox"/> Cal Time			
<b>Acceptance Limits</b>					
Cal Time Hour		SD			
Slope Diff		Repeatability			
Sensitivity					
Deter Coeff					

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