A-400 BIL DIRECT

Nr kat. **7-448** (PL)

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
(EN)

REAGENT STABILITY

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 10 weeks on board the analyser at 2-10°C (BS-400).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not freeze reagents.
- Do not interchange caps.
- Reagent bottles should be shaken before use by gently inverting several times.
- The appearance of turbidity or control sera values outside the manufacturer's acceptable range may indicate of the reagents instability.
- Lack of significant changes in the color of the reaction mixture at the samples with low bilirubin concentration does not indicate the assay malfunction.
- Do not use the reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-REAGENT contains hydroxylammonium chloride. May produce an allergic reaction (EUH208).

SPECIMEN ^{6,7,8}

Serum free from hemolysis.

Serum should be separated from red blood cells as soon as possible after blood collection.

Lipemic specimens may show falsely decreased bilirubin concentration thus fasting specimen is recommended.

Because bilirubin is photooxidized when exposed to light, specimen should be protected from direct exposure to either artificial light or sunlight. Therefore it is essential to store specimens in the dark! Serum can be stored at 4-8°C or at room temperature (20-25°C) up to 2 days.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.

Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

ASSAY PROCEDURE

Reagents are ready to use.

For reagent blank deionized water is recommended.

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of automatic analysers the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176)

INTENDED USE

BIL DIRECT reagent is intended to determine quantitatively direct bilirubin level in serum. It is intended to diagnosis, monitoring and as an aid to diagnosis of clinical conditions connected with abnormal bilirubin level. BIL DIRECT reagent is intended to use on automatic analysers BS-400 and BS-480. It is only for *in vitro* diagnostics, for healthcare professional users.

SUMMARY ^{1,2,3}

Determinations of direct bilirubin level in serum are mainly used in diagnosis of liver functions. Measurement of direct bilirubin concentration is used as an aid in diagnosis of clinical conditions associated with hyperbilirubinemia. Increased direct bilirubin levels can be observed in mechanical jaundice, Dubin-Johnson syndrome, bile ducts diseases, gallbladder diseases. Direct bilirubin is also the marker of perforation in acute sigmoid diverticulitis.

METHOD PRINCIPLE ^{4,5}

Method is based on chemical oxidation, utilizing vanadate as an oxidizing agent.

In the presence of detergent and vanadate in an acidic solution, conjugated (direct) bilirubin is oxidized to produce biliverdin.

This oxidation reaction causes change of the yellow colour, which is specific to bilirubin to the green colour typical for biliverdin. Therefore, the direct bilirubin concentration in the sample can be obtained by measuring the absorbance before and after the vanadate oxidation.

REAGENTS

Package

1-REAGENT	1 x 40 ml
2-REAGENT	1 x 12.5 ml

Number of tests:

BS-400	200
BS-480	200

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-REAGENT

citrate buffer	
hydroxylammonium chloride	
stabilizer	
detergent	0.34 %
pH adjuster	
preservative	
2-REAGENT	
sodium metavanadate	4.1 mmol/l
phosphate buffer	
pH adjuster	
stabilizer	

and LEVEL 2 (Cat. No 5-175; 5-177) are recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 6 weeks (BS-400).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ³

serum (adults)	0.1 – 0.3 mg/dl 1.7 – 5.1 µmol/l
----------------	-------------------------------------

It is recommended for each laboratory to establish its own reference ranges for local population and method principle.

The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and/or BS-480. Results may vary if a different instrument is used.

- **LoB (Limit of Blank):**
0.00 mg/dl (0.00 µmol/l) – BS-400
- **LoD (Limit of Detection):**
0.01 mg/dl (0.171 µmol/l) – BS-400
- **LoQ (Limit of Quantitation):**
0.09 mg/dl (1.539 µmol/l) – BS-400
0.06 mg/dl (1.026 µmol/l) – BS-480
- **Linearity:**
up to 40 mg/dl (684 µmol/l) – BS-400, BS-480
- **Measurement range:**
0.09 mg/dl (1.539 µmol/l) – 40 mg/dl (684 µmol/l) – BS-400
0.06 mg/dl (1.026 µmol/l) – 40 mg/dl (684 µmol/l) – BS-480

- **Specificity / Interferences**

Ascorbic acid up to 62 mg/l and triglycerides up to 550 mg/dl (in samples containing low concentration of direct bilirubin) and up to 900 mg/dl (in samples containing high concentration of direct bilirubin) do not interfere with the test. Haemoglobin interferes even in small amount with the determination.

- **Precision**

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 20)	poziom 1	0.36	0.00	0.8
	poziom 2	1.14	0.02	1.9
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
BS-400 (n = 80)	poziom 1	0.37	0.01	3.2
	poziom 2	1.15	0.06	5.2

- **Method comparison**

A comparison between direct bilirubin values determined at **BS-400** (y) and at **BECKMAN COULTER AU680** (x) using 119 serum samples gave following results:
 $y = 1.0769x - 0.0197$ mg/dl
 $R = 0.999$ (R – correlation coefficient)

A comparison between direct bilirubin values determined at **BS-480** (y) and at **BECKMAN COULTER AU680** (x) using 30 serum samples gave following results:
 $y = 1.1953x - 0.0526$ mg/dl
 $R = 0.997$ (R – correlation coefficient)

WASTE MANAGEMENT ⁹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents:
1-REAGENT – 18 01 07
2-REAGENT – 18 01 07
- Empty packages:
1-REAGENT – 15 01 02
2-REAGENT – 15 01 02
- Wastewater from the analyzer:
1-REAGENT – 18 01 03*
2-REAGENT – 18 01 03*

INCIDENTS ¹⁰

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

1. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 767-771 (2018).
2. Käser S.A., Fankhauser G., Glauser P.M. et al. Diagnostic Value of Inflammation Markers in Predicting Perforation in Acute Sigmoid Diverticulitis. World J Surg 34, 2717–2722 (2010).
3. Pagana K, Pagana T, Pagana T. Mosby's Diagnostic & Laboratory Test Reference 14th ed., St. Louis, Missouri, Elsevier, 137-139 (2019).

4. Tokuda K., Tanimoto K. A new method of measuring bilirubin in serum by vanadic acid. Jpn J Clin Chem (in Japanese with English abstract) 22: 116–22 (1993).
5. Dhungana N., Morris C., Krasowski M.D. Operational impact of using a vanadate oxidase method for direct bilirubin measurements at an academic medical center clinical laboratory. Practical laboratory medicine 8: 77-85 (2017).
6. Tietz N.W., ed. Clinical Guide to Laboratory Test. 3rd ed. Philadelphia: WB Saunders 518-525 (1995).
7. Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 1013 (1989).
8. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64.
9. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
10. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 03	Current version: 04
<i>Sections updated: INTENDED USE; SUMMARY; CONCENTRATION OF THE ACTIVE INGREDIENTS; WARNING AND NOTES; SPECIMEN; QUALITY CONTROL; REFERENCE VALUES; PERFORMANCE CHARACTERISTICS; WASTE MANAGEMENT; LITERATURE; APPLICATION</i>	
<i>Sections added: REAGENT STABILITY, INCIDENTS</i>	

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- Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 767-771 (2018).
- Käser S.A., Fankhauser G., Glauser P.M. et al. Diagnostic Value of Inflammation Markers in Predicting Perforation in Acute Sigmoid Diverticulitis. World J Surg 34, 2717-2722 (2010).
- Pagana K, Pagana T, Pagana T. Mosby's Diagnostic & Laboratory Test Reference 14th ed., St. Louis, Missouri, Elsevier, 137-139 (2019).
- Tokuda K., Tanimoto K. A new method of measuring bilirubin in serum by vanadic acid. Jpn J Clin Chem (in Japanese with English abstract) 22: 116-22 (1993).
- Dhungana N., Morris C., Krasowski M.D. Operational impact of using a vanadate oxidase method for direct bilirubin measurements at an academic medical center clinical laboratory. Practical laboratory medicine 8: 77-85 (2017).
- Tietz N.W., ed. Clinical Guide to Laboratory Test. 3rd ed. Philadelphia: WB Saunders 518-525 (1995).
- Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 1013 (1989).
- World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64.
- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.

10. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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• BS-480

Chem	BIL D	No.	006	Sample Type	SERUM
Chemistry	Bilirubin Direct			Print name	BIL D
Reaction Type	Endpoint	Reaction Direction	Decrease	Sec Wave	546
Pri Wave	450	Decimal	0.01	Reaction Time	80
Unit	mg/dL	Blank Time	48		82
Standard	5.5	Aspirated		Diluent	
Decreased	5.5	20	180	R1	160
Increased				R2	40
				R3	
				R4	
		Sample Blank	V	Auto Retun	
Linearity Range (Standard)	0.06	40	Linearity Limit		
Linearity Range (Decreased)			Substrate Depletion		
Linearity Range (Increased)			Mixed Blank Abs	-33000	33000
R1 Blank Abs	-33000	33000	Uncapping Time	70	Day(s)
Blank Response	-33000	33000	Reagent Alarm Limit		
Twin Chemistry			Enzyme Linear Extension		
	Prozone Check	Rate Check	Antigen Addition		
O1	0	O2	0	O3	0
PC	0	ABS	0	O4	0

Calibration Settings		Auto Calibration	
Math Model	Multi-point Linear		Bottle Changed
Factor			Lot Changed
Replicates	2		Cal Time
Acceptance Limits			
Cal Time	1008	Hour	
Slope Diff		SD	
Sensitivity		Repeatability	
Deter Coeff			

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PZ CORMAY S.A.