## A-400 AMYLASE

Nr kat. 7-455 (PL)

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Cat. No **7-455** (EN)

### INTENDED USE

Diagnostic kit for determination of  $\alpha$ -amylase activity intended to use in automatic analyzers: BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

$\alpha$ -Amylase is a digestive enzyme secreted by salivary glands and pancreas. Low level of amylase is also found in skeletal muscle, adipose tissue and fallopian tubes.  $\alpha$ -Amylase is measured generally in pancreas diseases. Elevation of amylase activity is observed also due to inflammation of abdominal cavity or salivary glands.

### METHOD PRINCIPLE

2-Chloro-4-nitrophenyl- $\alpha$ -maltotriose (CNP-G3) is a direct substrate for determination of  $\alpha$ -amylase activity, which does not require the presence of ancillary enzymes.

10 CNP-G3  $\alpha$ -amylase  $\rightarrow$  9 CNP + CNP-G2 + 9 maltotriose + glucose

The rate of 2-chloro-4-nitrophenol formation can be monitored at 405 nm and is proportional to the  $\alpha$ -amylase activity.

### REAGENTS

**Package**  
1-Reagent 2 x 38 ml

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. The reagents are stable for 8 weeks on board the analyzer at 2-10°C.

### Concentrations in the test

MES buffer < 120 mmol/l  
calcium acetate < 7 mmol/l  
potassium hydroxide < 40 mmol/l  
potassium thiocyanate < 1100 mmol/l  
2-chloro-4-nitrophenyl- $\alpha$ -maltotriose < 2mmol/l  
preservative, stabilizer

### Warnings and notes

- Protect from direct sunlight.
- Prevent the reagent from microbiological contamination and from saliva and sweat  $\alpha$ -amylase! Saliva and sweat contain  $\alpha$ -amylase. Do not pipette by the mouth, avoid skin contact with reagent, specimens, tips, cuvettes. Ensure to use automatic pipettes and laboratory gloves.

### SPECIMEN

Serum or plasma collected on heparin, free from hemolysis, urine.

Do not use anticoagulants: EDTA, citrates and oxalates as they inhibit amylase activity.

Serum / plasma can be stored for 7 days at 15-25°C or for one month at 2-8°C.<sup>7</sup>

A-400 AMYLASE

Urine can be stored for 2 days at 15-25°C or for 10 days at 2-8°C.<sup>9</sup> Amylase is very unstable in acid urine. Adjust pH to approximately 7.0 before storage. Nevertheless it is recommended to perform the assay with freshly collected samples!

### PROCEDURE

1-Reagent is ready to use.  
For reagent blank deionized water is recommended.

### REFERENCE VALUES <sup>8</sup>

serum / plasma	U/l	$\mu$ kat/l
	20 – 104	0.34 – 1.77
urine	U/l	$\mu$ kat/l
	32 – 641	0.54 – 10.90

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples the following controls:

CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum  
CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

**For the calibration of automatic analyzer system BS-400** the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

**For the calibration of automatic analyzer system BS-480** the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 5 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyzers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity:

2.6 U/l (0,043  $\mu$ kat/l) – BS-400  
3.23 U/l (0,054  $\mu$ kat/l) – BS-480

#### Linearity:

up to 1500 U/l (25  $\mu$ kat/l) – BS-400  
up to 1630 U/l (27  $\mu$ kat/l) – BS-480

#### Specificity / Interferences

Haemoglobin up to 2.5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

### Precision

Repeatability (run to run)		Mean [U/l]	SD [U/l]	CV [%]
<b>BS-400</b> (n = 20)	level 1	36.79	0.33	0.91
	level 2	504.70	1.87	0.37
<b>BS-480</b> (n = 10)	level 1	70.07	1.05	1.50
	level 2	399.23	1.57	0.39
Reproducibility (day to day)		Mean [U/l]	SD [U/l]	CV [%]
<b>BS-400</b> (n = 14)	level 1	77.86	3.41	4.38
	level 2	254.01	14.59	5.75
<b>BS-480</b> (n = 10)	level 1	70.07	1.05	1.50
	level 2	399.23	1.57	0.39

#### Method comparison

A comparison between amylase values determined at **BS-400** (y) and at **OLYMPUS AU640** (x) using 42 samples gave following results:

$y = 0.966 x - 5.834$  U/l;  
 $R = 0.999$  (R – correlation coefficient)

A comparison between amylase values determined at **BS-480** (y) and at **BS-800** (x) using 90 samples gave following results:

$y = 1.0013 x + 0.1359$  U/l;  
 $R = 0.999$  (R – correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

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<b>• Calibration</b>		<b>Judgment Criteria</b>	
<b>Calibration</b>		Sensitivity	
Rule	Two-point Linear	Blank Abs.	
Replicate	3	Error Limit	
K		Corr. Coeff.	
<b>• QC</b>			
<b>Rules</b>		<b>Auto QC</b>	
Westgard Multi-rule		Cum. Sum Check	
v	1-2S	v	R-4S
v	1-3S	v	4-1S
v	2-2S	v	10-X
		•	1.0 - 2.7
			1.0 - 3.0
			0.5 - 5.1
		Interval	

• **BS-480**

Chem	AMYL	No.	019	Sample Type	SERUM/URINE
Chemistry	Amylase	Print name	AMYL		
Reaction Type	Kinetic	Reaction Direction	Increase		
Pri Wave	412	Sec Wave	700		
Unit	U/L	Decimal	0.1		
Blank Time	0 0	Reaction Time	51 71		
Standard	Sample Vol 3 μL Aspirated 20 μL Diluent 180 μL	Reagent Vol	R1 200 μL	Diluent	
Decreased	3 μL 20 μL 180 μL	R2			
Increased		R3			
	Sample Blank v Auto Retun	R4			
Linearity Range	3.23 1630	Linearity Limit	0.2		
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	84 Day(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		Enzyme Linear Extension			
Prozone Check		Rate Check		• Antigen Addition	
O1 0	O2 0	O3 0	O4 0		
PC 0	ABS 0				

<b>Calibration Settings</b>		<b>Auto Calibration</b>	
Math Model	Multi-point Linear	<input type="checkbox"/> Bottle Changed	
Factor		<input type="checkbox"/> Lot Changed	
Replicates	3	<input type="checkbox"/> Cal Time	
<b>Acceptance Limits</b>			
Cal Time	840 Hour		
Slope Diff		SD	
Sensitivity		Repeatability	
Deter Coeff			

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