A-400 dTIBC

Nr kat. 7-465

(PL)

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sodium bicarbonate 325 mmol/l
buffer, preservatives, conservants

INTENDED USE

Diagnostic kit for determination of total iron binding capacity, intended to use in automatic analyzers: BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Iron is transported as Fe (III) bound to the apotransferrin. The apotransferrin-Fe (III) complex is called transferrin. Normally only about one third of the iron binding sites of transferrin are occupied by Fe (III). TIBC (Total iron-binding capacity) is a measurement for the maximum iron concentration that transferrin can bind.

Serum TIBC levels vary in disorders of iron metabolism, e.g. elevated TIBC level when iron levels are low is useful in the diagnosis and monitoring of iron deficiency anaemia and also anaemia in late pregnancy. Decreased level of TIBC may indicate disorders such as anaemia associated with chronic inflammatory disorders, haemochromatosis, other anaemias, thalassemia, cirrhosis and malignancies.

METHOD PRINCIPLE

Colorimetric method, direct, with chromazurol B.

The assay consists of 2 reaction steps:

1. The sample is added to 1-Reagent acidic buffer which contains excess iron and iron-binding dye (Chromazurol B).

The low pH causes transferrin to release all iron atoms which then are bound by the iron-binding dye along with the excess iron, producing a coloured dye-iron complex.

$Fe^{3+} + \text{Chromazurol B} + \text{Fe-transferrin} \xrightarrow{H^+} \text{Fe-Chromazurol B} + \text{transferrin}$

2. Neutral buffer 2-Reagent addition results in pH change and causes a large increase in affinity of transferrin for iron. Then serum transferrin rapidly binding iron ions separating it from the dye-iron complex. The observed decrease in the absorbance of a coloured dye-iron complex is directly proportional to the total iron-binding capacity TIBC in the serum sample.

$Fe\text{-Chromazurol B} + \text{transferrin} \xrightarrow{pH \text{ neutral}} \text{Fe-transferrin} + \text{Fe-Chromazurol B}$

REAGENTS

Package

1-Reagent 1 x 27 ml
2-Reagent 1 x 9 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 11 weeks on board the analyzer at 2-10°C.

Concentrations in the test

1-Reagent

sodium acetate 200 mmol/l
iron chloride 0.02 mmol/l

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyzers BS-400 and BS-480. Results may vary if a different instrument is used.

▪ Sensitivity

100 µg/dl (17.90 µmol/l) – BS-400
100 µg/dl (17.90 µmol/l) – BS-480

▪ Linearity

up to 600 µg/dl (107.40 µmol/l) – BS-400
up to 600 µg/dl (107.40 µmol/l) – BS-480

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

▪ Specificity / Interferences

Haemoglobin up to 1.25 g/dl, ascorbate up to 31 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl, copper up to 3 mg/dl, zinc up to 250 µg/dl, iron up to 357.5 µg/dl, Desferal up to 11.5 µg/ml, Cuprimine up to 250 µg/dl and Iron Dextran (Imferon) up to 1430 µg/ml do not interfere with the test.

▪ Precision

Repeatability (run to run)		Mean [µg/dl]	SD [µg/dl]	CV [%]
BS-400 (n = 10)	level 1	322	1.2	0.4
	level 2	223	0.8	0.4
BS-480 (n = 10)	level 1	337	3.0	0.9
	level 2	217	1.4	0.6
Reproducibility (day to day)		Mean [µg/dl]	SD [µg/dl]	CV [%]
BS-400 (n = 10)	level 1	357	5.0	1.4
	level 2	245	2.4	1.0
BS-480 (n = 10)	level 1	355	3.8	1.1
	level 2	228	3.6	1.6

▪ Method comparison

A comparison between TIBC values determined at **BS-400** (y) and at **Advia 1650** (x) using 149 samples gave following results:

$y = 0.9628 x + 6.373 \mu\text{g/dl}$;
 $R = 0.960$ (R – correlation coefficient)

A comparison between TIBC values determined at **BS-480** (y) and at **Advia 1650** (x) using 54 samples gave following results:

$y = 0.9988 x - 5.1433 \mu\text{g/dl}$;
 $R = 0.953$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

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• BS-480

Chem	dTIBC	No.	065	Sample Type	SERUM
Chemistry	Total Iron Binding Capacity	Print name	dTIBC		
Reaction Type	Endpoint	Reaction Direction	Decrease		
Pri Wave	660	Sec Wave			
Unit	ug/dl	Decimal	0		
Blank Time	48 49	Reaction Time	80 82		
Standard	Sample Vol 13.5 uL Aspirated 20 uL Diluent 180 uL	Reagent Vol	R1 180 uL R2 45 uL R3 uL R4 uL	Diluent	uL uL uL uL
Decreased	13.5 uL				
Increased	uL				
	<input type="checkbox"/> Sample Blank <input checked="" type="checkbox"/> Auto Retun				
Linearity Range (Standard)	100 600	Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-33000 33000		
R1 Blank Abs	-33000 33000	Uncapping Time	77 Dav(s)		
Blank Response	-33000 33000	Reagent Alarm Limit			
Twin Chemistry		<input type="checkbox"/> Enzyme Linear Extension			
	<input type="checkbox"/> Prozone Check <input type="checkbox"/> Rate Check	• Antigen Addition			
O1	0 O2 0	O3	0 O4 0		
PC	0 ABS 0				

Calibration Settings		Auto Calibration	
Math Model	Two-point Linear	<input type="checkbox"/> Bottle Changed	
Factor		<input type="checkbox"/> Lot Changed	
Replicates	3	<input type="checkbox"/> Cal Time	
Acceptance Limits			
Cal Time	1848 Hour	SD	
Slope Diff		Repeatability	
Sensitivity			
Deter Coeff			

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