A-400 ETHANOL

Nr kat. **7-478** (PL)

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LITERATURA

- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, Toronto: WB Saunders Co., 1171-1175 (1994).
- William H., Porter Ph.D. Clinical Toxicology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 922-923.
- Bonnichsen, R.K.,Theorell, H., An Enzymatic Method for the Micro Determination of Ethanol, Scand. J. Clin. Lab. Invest. 3, 58 (1951).
- Lundquist, F., The Determination of Ethyl Alcohol in Blood and Tissues., Methods of Biochemical Analysis, Vol III, Editor, D. Glick, Interscience Publishers, New York, 1957 p. 217-251.
- CLSI Guidelines and Standards, Clinical Laboratory Tests, AACC Press, Washington, 1990.
- Preeti Dalwari: [medscape.com/article/2090019-overview#a3](https://www.medscape.com/article/2090019-overview#a3).
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 923, (1999).
- Krajowa Izba Diagnostów Laboratoryjnych, Oznaczanie etanolu w krwi i moczu pacjentów w laboratorium medycznym, Warszawa 2016. p.18-19.

Data wydania: 10.2020.

A-400 ETHANOL

Cat. No **7-478** (EN)

INTENDED USE


Diagnostic kit for determination of ethyl alcohol concentration intended to use in automatic analyzers: BS-400 and BS-480. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Ethanol after consumption is absorbed in alimentary tract and transported with blood to whole organism. A little amount of ethanol is excreted with urine or exhaled with air, but majority is metabolized by liver to acetic aldehyde, then to acetate, finally to carbon dioxide and water. Ethanol consumption is often cause of various types of accidents, and ethanol poisoning may be cause of decease.

METHOD PRINCIPLE

Enzymatic method with alcohol dehydrogenase.

The rate of absorbance changes measured at $\lambda=340$ nm or $\lambda=380$ nm is proportional to the quantity of ethyl alcohol in the sample.

REAGENTS

Package
 1-Reagent 2 x 25 ml

The reagent, stored tightly capped at 2-8°C is stable up to expiry date printed on the package. The reagent, stored on board of the analyser at 2-10°C is stable for 10 weeks.

Concentrations in the test

buffer (pH 7,3) 1.35 mmol/l
 NAD analog >203 kU/l
 alcohol dehydrogenase
 stabilizers, preservatives

Warnings and notes

- Protect from direct sunlight and avoid contamination !
- Avoid contact with skin and mucous membranes.
- Concentration of ethyl alcohol in serum or plasma is approximately 1.2 times higher than in whole blood.⁶
- Turbidity would indicate contamination or deterioration of reagent.

SPECIMEN

Serum or plasma (heparin, EDTA or citrate) free from hemolysis; urine.

Samples may be stored tightly closed in temp. 2-8°C up to 3 days. Urine samples should be stored with as little dead air space as possible.

Nevertheless it is recommended to perform the assay with freshly collected samples!

- Do not use ethanol to swab venipuncture site and to clean or sterilize glassware or other equipment used to collect the specimen or perform the assay.
- Samples should be at room temperature for testing.
- Before opening mix samples by inverting the container.

PROCEDURE

1-Reagent is ready to use.
 Deionised water or 0.9% NaCl is recommended as a reagent blank.

REFERENCE VALUES 7

Ethanol consumption cause following symptoms:

blood ethanol concentration (mg / dl)	clinical symptoms*
30-120	euphoria; sociability; talkativeness; diminution of attention, judgment and control; some sensory-motor impairment
90-250	emotional instability; loss of critical judgment; impairment of perception, memory and comprehension; increased reaction time; sensory-motor incoordination; drowsiness;
180-300	disorientation; mental confusion; dizziness; staggering gait; slurred speech; disturbances of vision; increased pain threshold; apathy; lethargy;
250-400	general inertia; markedly decreased response to stimuli; vomiting; incontinence of urine and feces; sleep or stupor;
350-500	complete unconsciousness; coma; depressed or abolished reflexes; impairment of circulation and respiration; possible death.

*Effects are dependent on individual use, chronicity, tolerance and other factors.⁸

To calculate the concentration of ethyl alcohol in blood the result received for serum or plasma should be divided by 1.2 factor. Otherwise, doctor should be informed that the result is overstated and requires a proper interpretation.⁸

Please refer to the local regulations concerning legal ethanol limits for drivers.

Conversion factor of the blood ethanol concentration:
 100 mg/dl = 1 ‰.

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY AMMONIA/ETHANOL CONTROLS (Cat. No 5-163). For the calibration of automatic analysers the CORMAY ETHANOL CALIBRATORS (Cat. No 5-105) is recommended. Two levels of calibrator and deionised water or 0.9% NaCl (as a calibrator 0) should be used for calibration. The calibration curve should be prepared every 4 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analysers: BS-480 and/or BS-400 and Beckman CX. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
 0.2 mg/dl (0.043 mmol/l) – BS-400
 1.9 mg/dl (0.41 mmol/l) – BS-480
- LoD (Limit of Detection):**
 4.0 mg/dl (0.87 mmol/l) – BS-400
 3.1 mg/dl (0.67 mmol/l) – BS-480
- LoQ (Limit of Quantitation):**
 12.4 mg/dl (2.69 mmol/l) - BS-400
 7.0 mg/dl (1.52 mmol/l) - BS-480
- Linearity:**
 up to 800 mg/dl (173.6 mmol/l) – BS-400
 up to 800 mg/dl (173.6 mmol/l) – BS-480

Specificity / Interferences

Haemoglobin up to 0.6 g/dl, bilirubin up to 40 mg/dl, intralipid up to 1000 mg/dl, LDH up to 10 000 U/l (in the presence of 200 mg/dl lactate) do not interfere with the test.

Acetone up to 2000 mg/dl, butanol up to 200 mg/dl, ethylene glycol up to 2000 mg/dl, methanol up to 2000 mg/dl, izopropanol up to 2000 mg/dl do not interfere with the test, which confirms that assay is specific to ethanol.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	42.2	1.02	2.43
level 2	105.7	1.05	1.00
level 3	258.2	1.73	0.67
Reproducibility (day to day) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	42.32	2.44	5.76
level 2	101.06	4.42	4.37
level 3	252.63	8.69	3.44

Method comparison

A comparison between CORMAY reagent (y) and commercially available assay (x) using 40 serum samples, gave following results:
 $y = 1.064x - 2.33$ mg/dl;
 $R = 0.9994$ (R – correlation coefficient)

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, Toronto: WB Saunders Co., 1171-1175 (1994).
- William H., Porter Ph.D. Clinical Toxicology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 922-923.
- Bonnichsen, R.K.,Theorell, H., An Enzymatic Method for the Micro Determination of Ethanol, Scand. J. Clin. Lab. Invest. 3, 58 (1951).
- Lundquist, F., The Determination of Ethyl Alcohol in Blood and Tissues., Methods of Biochemical Analysis, Vol III, Editor, D. Glick, Interscience Publishers, New York, 1957 p. 217-251.
- CLSI Guidelines and Standards, Clinical Laboratory Tests, AACC Press, Washington, 1990.
- Preeti Dalawari: emedicine.medscape.com/article/2090019-overview#a3.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 923, (1999).
- Polish National Chamber of Laboratory Diagnosticians, Oznaczanie etanolu w krwi i moczu pacjentów w laboratorium medycznym, Warszawa 2016. p.18-19.

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- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, Toronto: WB Saunders Co., 1171-1175 (1994).
- William H., Porter Ph.D. Clinical Toxicology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 922-923.
- Bonnichsen, R.K.,Theorell, H., An Enzymatic Method for the Micro Determination of Ethanol, Scand. J. Clin. Lab. Invest. 3, 58 (1951).
- Lundquist, F., The Determination of Ethyl Alcohol in Blood and Tissues., Methods of Biochemical Analysis, Vol III, Editor, D. Glick, Interscience Publishers, New York, 1957 p. 217-251.
- CLSI Guidelines and Standards, Clinical Laboratory Tests, AACC Press, Washington, 1990.
- Preeti Dalawari: emedicine.medscape.com/article/2090019-overview#3.
- Burtis C.A., Ashwood E.R., ed. Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA: WB Saunders, 923, (1999).
- Polish National Chamber of Laboratory Diagnosticians, Oznaczenie etanolu w krwi i moczu pacjentów w laboratorium medycznym, Warszawa 2016. p.18-19.

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• Calibration	
Calibration	Judgment Criteria
Rule <input type="text" value="Multi-point Linear"/>	Sensitivity <input type="text"/>
Replicate <input type="text" value="2"/>	Factor Diff. <input type="text"/>
K <input type="text"/>	SD <input type="text"/>
• QC	Blank Abs. <input type="text"/>
Rules	Error Limit <input type="text"/>
Westgard Multi-rule	Corr. Coeff. <input type="text"/>
v 1-2S <input type="text"/>	
v 1-3S <input type="text"/>	
v 2-2S <input type="text"/>	
v R-4S <input type="text"/>	
v 4-1S <input type="text"/>	
v 10-X <input type="text"/>	
Cum. Sum Check <input type="text"/>	Auto QC
• 1.0 - 2.7	Interval <input type="text"/>
• 1.0 - 3.0	
• 0.5 - 5.1	

BS-480

Chem <input type="text" value="ETH"/>	No. <input type="text" value="071"/>	Sample Type <input type="text" value="SERUM"/>
Chemistry <input type="text" value="ETHANOL"/>	Print name <input type="text" value="ETH"/>	
Reaction Type <input type="text" value="Kinetic"/>	Reaction Direction <input type="text" value="Increase"/>	
Pri Wave <input type="text" value="380"/>	Sec Wave <input type="text" value="505"/>	
Unit <input type="text" value="mg/dL"/>	Decimal <input type="text" value="0.1"/>	
Blank Time <input type="text"/>	Reaction Time <input type="text" value="15"/> <input type="text" value="20"/>	
Standard <input type="text" value="3"/> <input type="text" value="μL"/>	Reagent Vol <input type="text" value="200"/> <input type="text" value="μL"/>	Diluent <input type="text" value="μL"/>
Aspirated <input type="text" value="20"/> <input type="text" value="μL"/>	R1 <input type="text" value="μL"/>	Diluent <input type="text" value="μL"/>
Diluent <input type="text" value="180"/> <input type="text" value="μL"/>	R2 <input type="text" value="μL"/>	
Decreased <input type="text" value="μL"/>	R3 <input type="text" value="μL"/>	
Increased <input type="text" value="μL"/>	R4 <input type="text" value="μL"/>	
Sample Blank <input type="text"/>	Auto Retun <input type="text"/>	
Linearity Range (Standard) <input type="text" value="7"/> <input type="text" value="800"/>	Linearity Limit <input type="text" value="0.2"/>	
Linearity Range (Decreased) <input type="text"/>	Substrate Depletion <input type="text"/>	
Linearity Range (Increased) <input type="text"/>	Mixed Blank Abs <input type="text"/>	
R1 Blank Abs <input type="text" value="-33000"/> <input type="text" value="33000"/>	Uncapping Time <input type="text" value="70"/> Day(s)	
Blank Response <input type="text" value="-33000"/> <input type="text" value="33000"/>	Reagent Alarm Limit <input type="text"/>	
Twin Chemistry <input type="text"/>	<input type="checkbox"/> Enzyme Linear Extension	
<input type="checkbox"/> Prozone Check <input type="checkbox"/> Rate Check	• Antigen Addition	
Q1 <input type="text" value="0"/> Q2 <input type="text" value="0"/>	Q3 <input type="text" value="0"/> Q4 <input type="text" value="0"/>	
PC <input type="text" value="0"/> ABS <input type="text" value="0"/>		
Calibration Settings	Auto Calibration	
Math Model <input type="text" value="Multi-point Linear"/>	<input type="checkbox"/> Bottle Changed	
Factor <input type="text"/>	<input type="checkbox"/> Lot Changed	
Replicates <input type="text" value="2"/>	<input type="checkbox"/> Cal Time	
Acceptance Limits		
Cal Time <input type="text" value="672"/> Hour		
Slope Diff <input type="text"/>	SD <input type="text"/>	
Sensitivity <input type="text"/>	Repeatability <input type="text"/>	
Deter Coeff <input type="text"/>		

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