## A-400 HDL DIRECT

Nr kat. 7-479 (PL)

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### LITERATURA

- Gotto, A.M. Lipoprotein metabolism and the etiology of hyperlipidemia. Hospital Practice 1988; 23 Suppl:1 4-13.
- Badimon, J. J., Badimon, L., Fuester V. Regression of Atherosclerotic Lesions by High Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit. J Clin Invest 1990; 85:1234-41.
- Warnick, G. Russell, Wood, Peter D. National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol: Executive Summary. Clin Chem 1995; 41(10):1427-1433.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 564 (2006).
- Campois, J. Altered Composition of Lipoproteins in Liver Cirrhosis Compromises Three Homogeneous Methods for HDL-Cholesterol. Clinical Chemistry, 1999; 45:685-688

Data wydania: 10.2020

## A-400 HDL DIRECT

Cat. No **7-479** (EN)

### INTENDED USE

Diagnostic kit for determination of HDL-cholesterol concentration (direct method), used in automatic analysers BS-400 and BS-480.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

### INTRODUCTION

Plasma lipoproteins are spherical particles containing varying amounts of cholesterol, triglycerides, phospholipids and proteins. The relative protein and lipid determine the density of these lipoproteins and provide the basis on which to begin their classification. The classes are: chylomicron, very-low-density lipoprotein (VLDL), low-density-lipoprotein (LDL) and high-density lipoprotein (HDL). The principle role of HDL in lipid metabolism is the uptake and transport of cholesterol from peripheral tissues to the liver. Low HDL cholesterol (HDL-C) levels are strongly associated with an increased risk of coronary artery disease. The HDL-C determination is used to identify high-risk patients.


### METHOD PRINCIPLE

The assay is a homogeneous method for directly measuring HDL-cholesterol concentration in serum or plasma, without any off-line pretreatment or centrifugation steps.

Accelerator selective detergent methodology.

During the first phase, LDL, VLDL particles and Chylomicrons generate free non-HDL cholesterol, which through an enzymatic reaction, produce hydrogen peroxide. The generated peroxide is consumed by a peroxidase reaction with DSBmT yielding a colourless product.

During the second phase, specific detergent solubilises HDL-Cholesterol. In conjunction with cholesterol oxidase (CO) and cholesterol esterase (CE) action, peroxidase and 4-AAP develop a coloured reaction which is proportional to HDL-Cholesterol concentration.

### REAGENTS

#### Package

1-Reagent 3 x 29 ml  
2-Reagent 3 x 11.6 ml

The reagents are stable up to the kit expiry date printed on the package when stored at 2-8°C. The reagents are stable for 12 weeks on board the analyser at 2-10°C.

### Concentrations in the test

#### 1-Reagent

Buffer  
Cholesterol oxidase (*E.coli*) < 1000 U/l  
Peroxidase (horseradish) < 1300 ppg U/l  
N,N-bis(sulfobutyl)-toluidine, disodium (DSBmT) < 1 mM  
Accelerator < 1 mM  
Preservative < 0.06 %  
Ascorbic acid oxidase (*Cureubita sp.*) < 3000 U/l

#### 2-Reagent

Buffer  
Cholesterol esterase (*Pseudomonas sp.*) < 1500 U/l  
4-aminoantipyrine (4-AAP) < 1 mM  
Detergent < 2 %  
Preservative < 0.06 %

### Warnings and notes

- Protect from direct sunlight and avoid contamination!

### SPECIMEN

Serum, heparinized or EDTA plasma.

Anticoagulants containing citrate should not be used.

Blood should be collected only if the patient has been fasting for 12 -14 hours.

Serum: Collect whole blood by venepuncture and allow to clot. Centrifuge and remove the serum as soon as possible after collection (within 3 hours).

Plasma: Specimens may be collected in EDTA or lithium or sodium heparin. Centrifuge and remove the plasma as soon as possible after collection (within 3 hours).

Serum and plasma should not remain at 15-30°C longer than 14 hours. If assays are not completed within 14 hours, serum or plasma should be stored at 2-8°C for up to 1 week. If specimens need to be stored for more than 1 week, they may be preserved at less than -20°C for up to 3 months. Samples may be frozen once.

Nevertheless it is recommended to perform the assay with freshly collected samples.

### PROCEDURE

1-Reagent and 2-Reagent are ready to use.

For reagent blank deionized water is recommended.

#### Actions required:

When performing assays at analyser **BS-400**, there is a probability of **cross-contamination** affecting the tests results: HDL DIRECT II GEN – TG, HDL DIRECT II GEN – TG mono, HDL DIRECT II GEN – UA, HDL DIRECT II GEN – URINE PROTEINS II GEN, HDL DIRECT II GEN – UA PLUS. To avoid this effect follow the recommendations contained in the instruction 51\_03\_24\_002\_BS-400\_CARRYOVER.

### REFERENCE VALUES <sup>4</sup>

serum / plasma	40 – 60 mg/dl	1.04 – 1.55 mmol/l
----------------	---------------	--------------------

As HDL cholesterol is affected by a number of factors such as smoking, exercise, hormones, age and sex, each laboratory should establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples CORMAY LIPID CONTROL 1 (Cat. No 5-179) and CORMAY LIPID CONTROL 2 (Cat. No 5-180) or CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of automatic analysers the CORMAY HDL/LDL CALIBRATOR (Cat. No 5-178) is recommended. Deionised water should be used as a calibrator 0.

The calibration curve should be prepared every 12 weeks with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers BS-400 and BS-480. Results may vary if a different instrument or a manual procedure is used.

#### Sensitivity:

0.4 mg/dl (0.010 mmol/l) – BS-400

0.4 mg/dl (0.010 mmol/l) – BS-480

#### Linearity:

up to 140 mg/dl (3.63 mmol/l) – BS-400

up to 170 mg/dl (4.40 mmol/l) – BS-480

For higher concentrations, dilute sample with 0.9% NaCl and repeat the assay. Multiply the result by the dilution factor.

#### Specificity / Interferences

Bilirubin conjugated up to 60 mg/dl, bilirubin total up to 60 mg/dl, haemoglobin up to 1 g/dl, ascorbic acid up to 100 mg/dl, intralipid up to 1800 mg/dl, triglycerides up to 2000 mg/dl and gamma-globulins up to 5000 mg/dl do not interfere with the test.

#### Precision

Repeatability (run to run)		Mean [mg/dl]	SD [mg/dl]	CV [%]
<b>BS-400</b> n=10	level 1	44.26	1.19	2.68
	level 2	57.58	1.53	2.65
<b>BS-480</b> n=10	level 1	46.26	0.23	0.49
	level 2	60.34	0.41	0.68
Reproducibility (day to day)		Mean [mg/dl]	SD [mg/dl]	CV [%]
<b>BS-400</b> n=20	level 1	44.18	1.35	3.06
	level 2	57.58	1.63	2.84
<b>BS-480</b> n=20	level 1	32.65	0.62	1.90
	level 2	66.27	1.16	1.76

#### Method comparison

A comparison between HDL cholesterol values determined at **BS-400** (y) and at **Hitachi 912** (x) using 23 samples gave following results:

y = 0.9603 x – 1.2892 mg/dl

R = 0.981

(R - correlation coefficient)

A comparison between HDL cholesterol values determined at **BS-480** (y) and at **COBAS INTEGRA 400 Plus** (x) using 38 samples gave following results:

y = 0.8454 x + 6.6173 mg/dl

R = 0.978

(R - correlation coefficient)

### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- Gotto, A.M. Lipoprotein metabolism and the etiology of hyperlipidemia. Hospital Practice 1988; 23 Suppl:1 4-13.
- Badimon, J. J., Badimon, L., Fuester V. Regression of Atherosclerotic Lesions by High Density Lipoprotein Plasma Fraction in the Cholesterol-Fed Rabbit. J Clin Invest 1990; 85:1234-41.
- Warnick, G. Russell, Wood, Peter D. National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol: Executive Summary. Clin Chem 1995; 41(10):1427-1433.
- Alan H.B. Wu: Tietz Clinical Guide to Laboratory Tests, 4th ed. WB Saunders, 564 (2006).
- Camps, J, Altered Composition of Lipoproteins in Liver Cirrhosis Compromises Three Homogeneous Methods for HDL-Cholesterol, Clinical Chemistry, 1999; 45:685-688.

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• **BS-480**

Chem	HDL D	No.	011	Sample Type	SERUM
Chemistry	HDL Direct	Print name	HDL D		
Reaction Type	Endpoint	Reaction Direction	Increase		
Pri Wave	605	Sec Wave	800		
Unit	mg/dL	Decimal	0.1		
Blank Time	48   49	Reaction Time	80   82		
Standard	Sample Vol   Aspirated   Diluent	Reagent Vol	Diluent		
Decreased	3   20   180	R1	180		
Increased		R2	60		
		R3			
		R4			
Linearity Range (Standard)	0.4   170	Linearity Limit			
Linearity Range (Decreased)		Substrate Depletion			
Linearity Range (Increased)		Mixed Blank Abs	-33000   33000		
R1 Blank Abs	-33000   33000	Uncapping Time	84 Day(s)		
Blank Response	-33000   33000	Reagent Alarm Limit			
Twin Chemistry		Enzyme Linear Extension			
	Prozone Check   Rate Check	Antigen Addition			
Q1	0   0	Q3	0   0		
PC	0   0	ABS	0   0		

<b>Calibration Settings</b>		<b>Auto Calibration</b>	
Math Model	Two-point Linear	<input type="checkbox"/>	Bottle Changed
Factor		<input type="checkbox"/>	Lot Changed
Replicates	3	<input type="checkbox"/>	Cal Time
<b>Acceptance Limits</b>			
Cal Time	2016 Hour	SD	
Slope Diff		Repeatability	
Sensitivity			
Deter Coeff			

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