HAEM 8 CONTROL / HAEM 12 CONTROL

Nazwa zestawu	Objętość zestawu	(PL) Nr kat.
HAEM 8 CONTROL L	1 x 2,5 ml	8-601
HAEM 8 CONTROL N	1 x 2,5 ml	8-602
HAEM 8 CONTROL H	1 x 2,5 ml	8-603
HAEM 12 CONTROL L	1 x 2,5 ml	8-604
HAEM 12 CONTROL N	1 x 2,5 ml	8-605
HAEM 12 CONTROL H	1 x 2,5 ml	8-606

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INTENDED USE

Control Blood is assigned for quality control of automated, semi-automated and manual procedures that measure whole blood parameters. Control material will provide values within the expected range indicated on the assay sheet when assayed routinely (like a patient sample) on a properly calibrated and functioning instrument.

Daily determination of these controls confirms the precision and accuracy of hematological analyzers.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

PRODUCT DESCRIPTION

Control Blood is the mixture composed of human RBC's, leukocytes and animal derived platelets in a plasma like fluid with preservatives.

INSTRUCTION FOR USE

- Remove the vial from the refrigerator and allow to warm at room temperature (15-30°C) for 15 minutes.
- Place the control on a mechanical mixer for 20 minutes or mix the control manually according step 3. Do not use a Vortex-mixer.
- Roll the vial forth and back for 30 seconds, holding the vial horizontally between the palms of the hands. Next, gently invert the vial 10 times. Do not shake. Continue mixing in this manner until the cells are completely suspended. Vials stored for a long time may require extra mixing.
- After mixing allow the vial to stand about 15 seconds for small air bubbles dispersion. Gently invert the vial 10 times immediately before sampling. Analyze the control blood as a patient sample.
- After sampling wipe the vial ring and the inner side of cap with lint-free swab and replace the cap.
- Place the vial back in the refrigerator within 30 minutes after measuring the controls. Store in upright position.

STORAGE AND STABILITY

Control Blood when unopened remains stable at 2-8°C until expiry date given on the product label. **Protect vials from overheating and freezing.**

Opened vials, when stored at 2-8°C and handled according to the instruction, remain stable for:

HAEM 8 CONTROL L, N, H	2 weeks
HAEM 12 CONTROL L, N, H	2 weeks

EXPECTED VALUES

The expected values and standard deviation for each control parameter are calculated from results of multiple determinations on analysers calibrated with whole blood and by manual reference methods. The values obtained on Control Blood should be within the expected range. The expected ranges listed on the assay sheet

are useful for routine control of analysers. It is not recommended to use control material for analysers calibration.

Assigned values are presented as a Mean and Range. The Mean is derived from replicate testing on instruments operated and maintained according to the manufacturer's instructions. The Range is an estimate of variation between laboratories and also takes into account inherent imprecision of the method and expected biological variability of the control material.

Verify that the lot number on the tube/vial matches the lot number on the table of assay values.

It is recommended for each laboratory to establish its own reference ranges for local population.

PRECAUTIONS

- Control material should be treated as potentially infectious.
- Human blood donors used in preparation of Control Blood were tested for antibodies to human immunodeficiency virus (anti-HIV 1,2), hepatitis C virus (anti-HCV), antibodies to Trypanosoma cruzi (T cruzi, the causative agent of Chagas disease), hepatitis B surface antigen (HBsAg), nucleic acid tests (NAT) for HCV ribonucleic acid (RNA), HIV-1 RNA, HBV deoxyribonucleic acid (DNA), West Nile virus (WNV) RNA, serologic test for syphilis and found to be non-reactive.
- No known test methods can provide complete assurance that human blood derived products will not transmit infectious diseases. Keep special precautions when handling or disposing of vials, as with patients samples.
- Do not interchange caps, the color of the cap is helpful on identifying the level of the control blood.

INDICATIONS OF DETERIORATION

After mixing, product should be similar in appearance to fresh whole blood. In unmixed vials, the supernatant may appear cloudy and reddish; this is normal and does not indicate deterioration.

DO NOT USE THE PRODUCT IF DETERIORATION IS SUSPECTED.

WASTE MANAGEMENT

Please refer to local legal requirements.

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