## LYSING REAGENT KX 21 CN FREE

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### INTENDED USE

LYSING REAGENT KX 21 CN FREE is intended to lyse of erythrocytes for hemoglobin release to allow further hemoglobin's assay on hematological analyzers: SYSMEX KX-21 and SYSMEX KX-21N. It prepares leukocytes for subpopulations differentiation process. It is only or *in vitro* diagnostics, for healthcare professional users.

### OPERATION PRINCIPLE

LYSING REAGENT KX 21 CN FREE is located on board the analyzer.  
Due to the lysing properties of the reagent, it is possible to disrupt the cell membrane of red blood cells and release cellular contents in preparation for further analysis. LYSING REAGENT KX 21 CN FREE contains surfactants that lyses the red blood cells and leukoprotective agents that protects the state of the leukocytes to permit the differentiation in subpopulations.

### CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

quaternary ammonium salts 0,87 %  
imidazole  
sodium chloride  
stabilizers

### REAGENT STABILITY

The reagent must be used before the expiry date indicated on the label.  
Once opened, the reagent is stable for 90 days at room temperature.

### WARNINGS AND NOTES

- Protect from direct sunlight!
- Do not use reagent beyond expiry date printed on the label.
- Do not combine the remaining reagent with the reagent from a new package.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Avoid discharge to sewage system or to environment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- LYSING REAGENT KX 21 CN FREE meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

H412 Harmful to aquatic life with long lasting effects.  
P273 Avoid release to the environment.

### DEVICES INTENDED TO BE USED IN COMBINATION WITH THE REAGENT

LYSING REAGENT KX 21 CN FREE is intended to be used in combination with the diluent, FLUSH / FLUSH SET cleaning solution and appropriate hematology analyzer.

Name	Cat. No
DILUENT PLUS	8-875
FLUSH (1 l)	8-832
FLUSH (500 ml)	8-891
FLUSH SET	8-896
SYSMEX KX21	N/A
SYSMEX KX21N	N/A

### ASSAY PROCEDURE

LYSING REAGENT KX 21 CN FREE is ready for use.  
Follow the analyzer's user manual.

Prior the first use, unopened reagent must be stored for minimum 24 hours at room temperature.

### WASTE MANAGEMENT<sup>1</sup>

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 06\*
- Empty packages: 15 01 10\*
- Wastewater from the analyzer: 18 01 03\*

### INCIDENTS<sup>2</sup>

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (websites address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

**Serious incident** means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

### LITERATURE

- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices.

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