

AccuDiag™
Cysticercosis IgG
(Taenia Solium)
ELISA Kit

Cat # 8105-35



Test	Cysticercosis IgG (T. solium)ELISA
Method	Enzyme Linked Immunosorbent Assay
Principle	Sandwich Complex
Detection Range	Semi-quantitative : Positive, Negative
Sample	5 µL
Total Time	~ 20 min.
Shelf Life	12 Months from the manufacturing date
Specificity	100%
Sensitivity	98.4%

INTENDED USE

The **Cysticercosis ELISA Test** is an enzyme linked immunosorbent assay (ELISA) for the qualitative identification of serum antibodies (IgG) to *Taenia solium* antigen.

SUMMARY AND EXPLANATION

The causative agents of *Taenia solium* infection are generally associated with contaminated food, infected workers, gastric reflux in tapeworm carriers, or contaminated water. These contamination routes are generally created by ingestion of *T. solium* eggs. Cysticercosis is known as an infection of the larval form (cysticerci) of *Taenia* in organs or tissues. Areas most affected by the infection revolve around the central nervous system. In addition, Cysticercosis can cause seizures and cranial pressure in the brain.

This Cysticercosis ELISA test should not be the sole method of diagnosis. Because of possible cross reactions with Echinococcosis, other testing procedures should be done so that the Echinococcus infection can be ruled out. Even though cyst vesicular antigen aids in sensitivity and specificity, other non-serological tests are recommended to verify accurate results of the *T. solium* infection.

TEST PRINCIPLE

The principle of the Cysticercosis ELISA test involves three incubation steps. Before the first incubation, the microwells are coated with *T. solium* cyst antigen. The patient's sera is added, and if there are any antibodies present, they will bind to the wells during the first incubation. Next, the wells must be washed of any test sample, and added to the wells at this point is the Enzyme Conjugate. During the second incubation, the Enzyme Conjugate will bind to any antibodies present. Before the third incubation step, more washings are necessary. Then a chromogen (tetramethylbenzidine or TMB) is added. With the presence of the Enzyme

Conjugate and the peroxidase causing the consumption of peroxide, the chromogen changes to a blue color. The blue color turns to a bright yellow color after the addition of the stop solution, which ends the reaction. ELISA readers can be used to obtain results or the results can be read visually.

SPECIMEN COLLECTION AND PREPARATION

Serological specimens should be collected under aseptic conditions. Coagulate blood and remove serum. Hemolysis is avoided through prompt separation of the serum from the clot. Serum should be stored at 2-8 °C if it is to be analyzed within a few days. Serum may be held for 3 to 6 months by storage at -20 °C or lower. Lipemic and strongly hemolytic serum should be avoided. Do not heat inactive serum and avoid repeated freezing and thawing of samples. Freeze sample at -20 °C or lower if not used immediately.

Do not heat inactivate serum and avoid repeated freezing and thawing of samples. Test samples: Make a 1:64 dilution of patients' sera using the dilution buffer (e.g. 5 µl sera and 315 µl dilution buffer).

MATERIALS AND COMPONENTS

Materials provided with the test kits

- Plate:** Microwells containing *T. solium* antigens - 96 test wells in a test strip holder.
- Enzyme Conjugate:** One (1) bottle containing 11 ml of Protein A conjugated to peroxidase.
- Positive Control:** One (1) vial containing 1 ml of diluted positive rabbit serum.
- Negative Control:** One (1) vial containing 1 ml of diluted negative human serum.
- TMB Substrate Solution:** One (1) bottle containing 11 ml of the chromogen tetramethylbenzidine (TMB).
- Wash Concentrate 20X:** One (1) bottle containing 25 ml of concentrated buffer and surfactant.
- Dilution Buffer:** One (1) or Two (2) bottles containing 30 ml of buffered protein solution.
- Stop Solution:** One (1) bottle containing 11 ml of 0.73 M phosphoric acid.

Materials required but not provided

- Pipettes
- Squeeze bottle for washing strips (narrow tip is recommended)
- Reagent grade water and graduate cylinder
- Tubes for sample dilution
- Absorbent paper
- ELISA plate reader with a 450 nm and a 620-650 nm filter(optional if results are read visually).

Preparation

Wash Buffer - Remove cap and add contents of bottle to 475 ml of reagent grade water. Place diluted wash buffer into a squeeze bottle with a narrow tip opening. **Note:** Washings consist of filling to the top of each well, shaking out the contents and refilling. Avoid generating bubbles in the wells during the washing steps.

ASSAY PROCEDURE

- Break off number of wells needed (two for controls plus number of samples) and place in strip holder.
- Add 100 µl (or two drops) of the negative control to well #1, 100 µl of the positive control to well #2 and 100 µl of the diluted (1:64) test samples to the remaining wells.
 Note: Negative and positive controls are supplied prediluted. Do not dilute further.
- Incubate at room temperature (15 to 25 °C) for 10 minutes.
- Shake out contents and wash 3 times with the diluted wash buffer.
- Add 100 µl of Enzyme Conjugate to each well.
- Incubate at room temperature for 5 minutes.

7. Shake out contents and wash 3 times with wash buffer. Slap wells against paper towels to remove excess moisture.
8. Add 100 µl of the Chromogen to every well.
9. Incubate at room temperature for 5 minutes.
10. Add 100 µl of the Stop Solution and mix by tapping strip holder.

RESULTS

Visually: Look at each well against a white background (e.g. paper towel) and record as clear or +, ++ or +++ reaction.

ELISA Reader: Zero reader on air. Set for bichromatic readings at 450/620-650 nm.

Troubleshooting

Negative control has excessive color after development.

Reason: inadequate washings

Correction: wash more vigorously. Remove excessive liquid from the wells by tapping against an Absorbent towel. Do not allow test wells to dry out.

Interpretation of the Test

Zero ELISA reader on air. Read all wells at 450/650 to 620 nm.

Positive –Absorbance reading equal to or greater than 0.3 OD units.

Negative –Absorbance reading less than 0.3 OD units.

Compare results to the controls. A sample should be interpreted as positive if the degree of color development is obvious and significant.

A positive OD reading indicates that the patient may be infected by *T. solium* or a closely related organism (e.g. *Echinococcus*).

A negative OD reading indicates that the patient has no detectable level of antibodies. This may be due to lack of infection or poor immune response by the patient.

QUALITY CONTROL

The use of controls allows validation of kit stability. The kit should not be used if any of the controls are out of range. Expected values for the controls are:

Positive – 0.5 OD units and above.

Negative - Absorbance reading 0-0.3 OD units.

EXPECTED VALUES

The number of individuals showing positive results can vary significantly between populations and geographic regions. If possible, each laboratory should establish an expected range for its patient population.

PERFORMANCE CHARACTERISTICS

		Reference Method *	
		+	-
Diagnostic Automation, Inc.	+	15	1
	-	0	65

Positive Agreement: 100% (15/15)

Negative Agreement: 98.4% (65/66)

*Reference Method refers to a commercially available ELISA.

LIMITATIONS OF PROCEDURE

Serologic results are an aid in diagnosis but cannot be used as the sole method of diagnosis. Significant cross reactions with *Echinococcus* infections will occur in this assay. If *Echinococcus* infection cannot be ruled out in the differential diagnosis, a positive sample should be confirmed by other means (i.e. immunoblot offered by the CDC) or by other non-serological means.

PRECAUTIONS

1. Do not use solutions if they precipitate or become cloudy. Wash concentrate may show crystallization upon storage at 2 – 8 °C. Crystallization will disappear after dilution to working strength.
2. Do not use serum that may have supported microbial growth, or is cloudy due to high lipid content. Samples high in lipids should be clarified before use.
3. Treat all sera as if capable of being infectious. Negative control has been tested and found negative for Hepatitis B surface antigen and for the antibody to HIV by required test methods. This product should be used under appropriate safety conditions that would be used for any potentially infectious agent.
4. Do not add azides to the samples or any of the reagents.

STORAGE



Reagents, strips and bottled components:

Store between 2 – 8 °C.

Squeeze bottle containing diluted wash buffer may be stored at room temperature.

REFERENCES

1. Flisser, A. and Larralde, C., Cysticercosis. Immunodiagnosis of Parasitic Diseases, Vol. 1, *Helminthic Diseases*. Ed. Walls and Schantz. Academic Press, 1986. pp. 109-149
2. Evans, C. et. Al., Controversies in the Management of Cysticercosis, *Emerging Infectious Diseases*, Vol. 3, No. 3, July-September 1997, pp.403-405
3. Larralde, C. et. al., Reliable serology of Taenia solium cysticercosis with antigens from cyst vesicular fluid: ELISA and hemagglutination tests. *Am J Trop Med Hyg.* Vol.35#5, 1986. pp. 965-973.
4. Rhoads, M. and Murrell, D., Taeniasis and Cysticercosis, *Laboratory Diagnosis of Infectious Diseases, Principles and Practice*, Vol.1, 1988, pp.987-992
5. Del Brutto, O. and Sotelo, J., Neurocysticercosis: An Update, *Reviews of Infectious Diseases*, Vol. 10, No. 6, November-December 1988, pp. 1075-1087

<p>ISO 13485 ISO 9001</p> 			
 <p>Diagnostic Automation/Cortez Diagnostics, Inc. 23961 Craftsman Road, Suite E/F, Calabasas, California 91302 USA</p>			
Date Adopted	Cat # 8105-35		
2013-09-16	<p>AccuDiag™- Cysticercosis IgG (Taenia Solium) ELISA -2013</p> <p>CEpartner4U, Esdoornlaan 13, 3951DB Maarn. The Netherlands. www.cepartner4u.eu</p>		
<table border="1" style="display: inline-table;"> <tr> <td>EC</td> <td>REP</td> </tr> </table>	EC	REP	
EC	REP		
Revision C Date: 12-11-2013			