

**AccuDiag™**  
**E. histolytica IgG (Amebiasis)**  
**ELISA Kit**

**Cat # 8201-35**



<b>Test</b>	<b>E. histolytica IgG (Amebiasis) ELISA</b>
<b>Method</b>	<b>Enzyme Linked Immunosorbent Assay</b>
<b>Principle</b>	<b>Sandwich Complex</b>
<b>Detection Range</b>	<b>Qualitative : Positive, Negative Control</b>
<b>Sample</b>	<b>5 µL serum</b>
<b>Total Time</b>	<b>~ 20 min.</b>
<b>Shelf Life</b>	<b>12 Months from the manufacturing date</b>
<b>Specificity</b>	<b>100 %</b>
<b>Sensitivity</b>	<b>100%</b>

**INTENDED USE**

The E. Histolytica (Amebiasis) ELISA Test is an enzyme linked immunosorbent assay (ELISA) for the qualitative identification of serum IgG antibodies to Entamoeba histolytica.

**SUMMARY AND EXPLANATION**

Generally, the disease of Amebiasis is found in a number of tropical regions where living conditions and poor sanitation cause significant health problems. Transmission of the disease centers around native populations and tourists traveling from these areas. A protozoan parasite called Entamoeba histolytica is the causing agent of Amebiasis, and the disease usually shows up as intestinal problems. Symptoms are generally mild, but in some cases the organism becomes extra-intestinal and can lead to abscesses, primarily affecting the liver.

Serological tests are recommended for extra-intestinal diagnosis, making sure to isolate the disease from other diseases of the liver, or ulcerative colitis, for example. This E. histolytica (Amebiasis) ELISA Test should not be used for diagnosing intestinal infections. Intestinal infections are conventionally established through an Ova and Parasite (O&P) test, or an E. histolytic fecal antigen assay. A positive result may not automatically be evidence of an active infection, and a negative outcome at least assures exclusion of a suspected E. histolytica tissue invasion.

**TEST PRINCIPLE**

The principle of the E. histolytica (Amebiasis) ELISA test is a three-incubation process whereby the first incubation involves the coating of the wells with E. histolytic antigen. During this step, any antibodies that are reactive with the antigen, will bind to the wells. Next, the wells must be washed to remove test sample. At this point Enzyme Conjugate is added. During this second incubation, the Enzyme Conjugate will bind to any antibodies present. Before the third incubation step, more washings are necessary. Then a chromogen (tetramethylbenzidine or TMB) is

added. With the presence of Enzyme Conjugate and the peroxidase causing the consumption of peroxide, the chromogen changes to a blue color. The blue color turns to a bright yellow color after the addition of the stop solution, which ends the reaction. ELISA readers can be used to obtain results, or the results can be read visually.

**SPECIMEN COLLECTION AND PREPARATION**

Coagulate blood and remove serum. Freeze sample at -20 °C or lower if not used immediately.

Do not heat inactivate serum and avoid repeated freezing and thawing of samples. Test samples: Make a 1:64 dilution of patients' sera using the dilution buffer (e.g. 5 µl sera and 315 µl dilution buffer).

**MATERIALS AND COMPONENTS**

**Materials provided with the test kits**

- Plate:** Microwells containing E. histolytica strain NIH-200 antigens - 96 test wells in a test strip holder.
- Enzyme Conjugate:** One (1) bottle containing 11ml of Protein-A conjugated to peroxidase.
- Positive Control:** One (1) vial containing 1 ml of diluted positive rabbit serum.
- Negative Control:** One (1) vial containing 1 ml of diluted human serum.
- TMB substrate Solution:** One (1) bottle containing 11 ml of the chromogen tetramethylbenzidine (TMB).
- Wash Concentrate 20X:** One (1) bottle containing 25 ml of concentrated buffer and surfactant.
- Dilution Buffer:** One (1) or Two (2) bottles containing 30 ml of buffered protein solution.
- Stop Solution:** One (1) bottle containing 11 ml of 0.73 M phosphoric acid.

**Materials required but not provided**

- Pipettes
- Squeeze bottle for washing strips (narrow tip is recommended)
- Reagent grade water and graduated cylinder
- Tubes for sample dilution
- Absorbent paper
- ELISA plate reader with a 450 nm and a 620-650 nm filter (optional if results are read visually).

**Preparation**

Wash Buffer - Remove cap and add contents of bottle to 475 ml of reagent grade water. Place diluted wash buffer into a squeeze bottle with a narrow tip opening.

Note: Washings consist of filling to the top of each well, shaking out the contents and refilling.

Avoid generating bubbles in the wells during the washing steps.

**ASSAY PROCEDURE**

- Break off number of wells needed (two for controls plus number of samples) and place in strip holder.
- Add 100 µl (or two drops) of the negative control to well #1, 100 µl of the positive control to well #2 and 100 µl of the diluted (1:64) test samples to the remaining wells.  
**Note:** Negative and positive controls are supplied prediluted. Do not dilute further.
- Incubate at room temperature (15 to 25 °C) for 10 minutes.
- Shake out contents and wash 3 times with the diluted wash buffer.
- Add 100 µl of Enzyme Conjugate to each well.
- Incubate at room temperature for 5 minutes.
- Shake out contents and wash 3 times with wash buffer. Slap wells against paper towels to remove excess moisture.
- Add 100 µl of the Chromogen to every well.
- Incubate at room temperature for 5 minutes.
- Add 100 µl of the Stop Solution and mix by tapping strip holder.



## RESULTS

**Visually:** Look at each well against a white background (e.g. paper towel) and record as clear or +, ++ or +++ reaction.

**ELISA Reader:** Zero reader on air. Set for bichromatic readings at 450/620-650 nm.

### Troubleshooting

Negative control has excessive color after development.

**Reason:** inadequate washings

**Correction:** wash more vigorously. Remove excessive liquid from the wells by tapping against an Absorbent towel. Do not allow test wells to dry out.

### Interpretation of the Results- ELISA Reader

Zero ELISA reader on air. Read all wells at 450/650-620 nm.

**Positive** - Absorbance reading greater than 0.4 OD units.

**Negative** - Absorbance reading less than 0.4 OD units.

A positive OD reading indicates that the patient may be infected by *E. histolytica*.

A negative OD reading indicates that the patient has no detectable level of antibodies. This may be due to lack of infection or poor immune response by the patient.

### Interpretation of Results -Visual

Compare results to the controls. A sample should be interpreted as positive if the degree of color is significant and obvious.

## QUALITY CONTROL

The use of controls allows validation of kit stability. The kit should not be used if any of the controls are out of range.

Expected values for the controls are:

**Negative** - 0.0 to 0.3 OD units

**Positive** - 0.5 OD units and above

## PERFORMANCE CHARACTERISTICS

		Reference Method *	
		+	-
Diagnostic Automation, Inc.	+	14	0
	-	0	65

Positive Agreement: 100% (14/14)

Negative Agreement: 100% (65/65)

\*Reference Method refers to a commercially available ELISA.

## LIMITATIONS OF PROCEDURE

Serologic results are an aid in diagnosis but cannot be used as the sole method of diagnosis.

## EXPECTED VALUES

The number of individuals showing positive results can vary significantly between populations and geographic regions. If possible, each laboratory should establish an expected range for its patient population.

## PRECAUTIONS

- Do not use solutions if they precipitate or become cloudy. Wash concentrate may show crystallization upon storage at 2-8 °C. Crystallization will disappear after dilution to working strength.

- Do not use serum that may have supported microbial growth, or is cloudy due to high lipid content. Samples high in lipids should be clarified before use.
- Treat all sera as if capable of being infectious. Negative control has been tested and found negative for Hepatitis B surface antigen and for the antibody to HIV by required test methods.

This product should be used under appropriate safety conditions that would be used for any potentially infectious agent.



Do not add azides to the samples or any of the reagents.

## STORAGE

- Reagents, strips and bottled components should be stored at 2-8 °C.
- Squeeze bottle containing diluted wash buffer may be stored at room temperature).

## REFERENCES

- Patterson, M. et. al. Serologic Testing for Amebiasis. *Gastroenterology*. 78:136, 1980
- Healy, G. Laboratory Diagnosis of Amebiasis. *Bull NY Acad Med*. 47:478, 1971
- Healy, G. Immunologic Tools in the Diagnosis of Amebiasis: Epidemiology in the United States. *Rev Infect Diseases*. Vol.8,#2:239, 1986
- Walsh, J. Problems in Recognition and Diagnosis of Amebiasis: Estimation of the Global Magnitude of Morbidity and Mortality. *Rev Infect Diseases*. Vol.8,#2:228, 1986

 <b>ISO 13485</b> <b>ISO 9001</b>	
 <b>Diagnostic Automation/ Cortez Diagnostics, Inc.</b> 23961 Craftsman Road, Suite E/F, Calabasas, California 91302 USA	
<b>Date Adopted</b>	<b>Cat # 8201-35</b>
<b>2013-10-18</b>	<b>AccuDiag™- E. histolytica IgG (Amebiasis) ELISA -2013</b>
<b>EC REP</b>	<b>CEpartner4U, Esdoornlaan 13, 3951DB Maarn. The Netherlands. <a href="http://www.cepartner4u.eu">www.cepartner4u.eu</a></b>
<b>Revision B Date: 12-27-2013</b>	