TOTAL PROTEIN

Nazwa zestawu	(PL) Nr kat.	HC- TOTAL PROTEIN	OS- TOTAL PROTEIN	B50- TOTAL PROTEIN
Liquick Cor-TOTAL PROTEIN 30	2-240			
Liquick Cor-TOTAL PROTEIN 60	2-236	1-REAGENT	6 x 96,5 ml	3 x 54 ml
Liquick Cor-TOTAL PROTEIN 120	2-237			4 x 58 ml
HC-TOTAL PROTEIN	4-536	2-STANDARD jest wzorcowym roztworem albuminy: 4 g/dl (40 g/l).		
OS-TOTAL PROTEIN	9-404			
B50-TOTAL PROTEIN	5-501			

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str. / page / crp. 4/13

TOTAL PROTEIN

Kit name	(EN) Cat. No
Liquick Cor-TOTAL PROTEIN 30	2-240
Liquick Cor-TOTAL PROTEIN 60	2-236
Liquick Cor-TOTAL PROTEIN 120	2-237
HC-TOTAL PROTEIN	4-536
OS-TOTAL PROTEIN	9-404
B50-TOTAL PROTEIN	5-501

INTENDED USE

TOTAL PROTEIN reagent is intended to determine quantitatively total protein level in plasma and serum. It is intended to monitoring and as an aid in the diagnosis of clinical conditions associated with abnormal protein concentration (hypo-, hyperproteinemia). TOTAL PROTEIN reagent is intended to use on automatic analysers and manually. 2-STANDARD is intended to calibrate manual determinations of total protein, performed with Liquick Cor - TOTAL PROTEIN 30 and Liquick Cor - TOTAL PROTEIN 60 reagent variants.

It is only for *in vitro* diagnostics, for healthcare professional users.

SUMMARY 1-5

Determinations of total protein level in serum and plasma are used to diagnosis and monitoring of clinical conditions associated with abnormal protein concentration. Increased total protein level can be observed in multiple myeloma, Waldenstrom macroglobulinemia, dehydration, hyperimmunoglobulinemia caused by polyclonal or monoclonal gammopathies, acute phase response. Decreased values can be observed in malnutrition, malabsorption, liver diseases (cirrhosis, fulminant failure), kidney diseases, protein-losing gastroenteropathies and acute burns.

METHOD PRINCIPLE 4

The method is based on the biuret reaction. Protein forms the coloured complex with cupric ions in alkaline medium. The colour intensity is proportional to the protein concentration.

REAGENTS

Package

	Liquick Cor- TOTAL PROTEIN 30	Liquick Cor- TOTAL PROTEIN 60	Liquick Cor- TOTAL PROTEIN 120
1-TOTAL PROTEIN	6 x 30 ml	6 x 60 ml	6 x 120 ml
2-STANDARD	1 x 2 ml	1 x 2 ml	-
	HC- TOTAL PROTEIN	OS- TOTAL PROTEIN	B50- TOTAL PROTEIN
1-REAGENT	6 x 96.5 ml	3 x 54 ml	4 x 58 ml

2-STANDARD is albumin standard solution: 4.0 g/dl (40 g/l).

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

1-TOTAL PROTEIN / 1-REAGENT	
sodium-potassium tartrate	30 mmol/l
potassium iodide	30 mmol/l
copper sulfate	12 mmol/l
sodium hydroxide	
2-STANDARD	
bovine albumin stabilizer	4 g/dl
preservative	

REAGENT STABILITY

The reagent when stored at 2-8°C is stable up to expiry date printed on the package. The reagent is stable for 11 weeks on board the analyser at 2-10°C (Biolis 30i).

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-TOTAL PROTEIN / 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-TOTAL PROTEIN / 1-REAGENT contains sodium hydroxide and copper (II) sulfate pentahydrate.

Danger

H314 Causes severe skin burns and eye damage.
H318 Causes serious eye damage.
H411 Toxic to aquatic life with long lasting effects.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P301+P330+P331 IF SWALLOWED: rinse mouth. Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P273 Avoid release to the environment.

EQUIPMENT REQUIRED BUT NOT PROVIDED

- automatic analyzer or photometer able to read at 546 nm;
- thermostat at 37°C;
- general laboratory equipment.

SPECIMEN 2, 6, 7

Serum or plasma. Nonhemolyzed, nonlipaemic, fresh serum is recommended.

Recommended anticoagulants: EDTA, heparine.

Serum should be separated from red blood cells as soon as possible after blood collection.

Serum and plasma can be stored up to 3 days at 2-8°C or 6 months at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

Follow tube manufacturers' instructions carefully when using collection tubes.


Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

ASSAY PROCEDURE

The reagent is ready to use.

Carry-over effect may occur when the reagent is used on automated chemistry analyzers. It is based on interference between specific test combinations and results in either underestimating or inflating patient sample results. To minimize this effect, the preventive actions should be taken, such as ordering an extra washing program or performing tests in a separate order.

General procedure for automatic determinations:

Main Wavelength: 546 nm

2nd wavelength: 700 nm

Method: End Point

Calibration Type: Linear

Direction: Increasing

R1 Volume: 250 µl

Sample Volume: 6 µl

Temperature: 37°C

Application guides for analyzers are available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

Manual procedure

wavelength	546 nm
temperature	37°C
cuvette	1 cm

Pipette into the cuvettes:

	reagent blank (RB)	test (T)	standard (S)
1-TOTAL PROTEIN / 1-REAGENT	1000 µl	1000 µl	1000 µl

Bring up to the temperature of determination. Then add:

standard / calibrator	-	-	20 µl
sample	-	20 µl	-

Mix well, incubate for 5 min. at the temperature of determination. Read the absorbance of the test A(T) and standard A(S) against reagent blank (RB). The intensity of colour is stable for 30 min.

Calculation

$$\text{protein concentration} = \frac{A(T)}{A(S)} \times \text{standard / calibrator concentration}$$

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173).

For the calibration of manual assay the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) or 2-STANDARD.

For calibration of the automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175, 5-177) is recommended. Deionised water should be used as a calibrator 0.

Calibration stability depends on the type of analyser used for analysis. The calibration curve should be prepared every 2 weeks (Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,
- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES 4,8

serum	g/dl	g/l
premature 1 day	3.6 – 6.0	36 – 60
children 1 day – 2 weeks	5.3 – 8.3	53 – 83
15 days – 12 months	4.4 – 7.1	44 – 71
1 – 3 years	6.1 – 7.5	61 – 75
3 – 5 years	6.3 – 8.1	63 – 81
6 – 19 years	6.8 – 8.2	68 – 82
adults 20 – 29 years	6.5 – 8.3	65 – 83
30 – 79 years	6.5 – 7.8	65 – 78

Results obtained from plasma analysis might be slightly elevated due to fibrinogen presence:

Sample	Protein concentration increase from serum to plasma
Blood donors	+ 0.25 g/dl
Nonhospitalized patients	+ 0.36 g/dl
Hospitalized patients	+ 0.46 g/dl
Hospitalized patients with CRP<50 mg/dl	+ 0.66 g/dl

It is recommended for each laboratory to establish its own reference ranges for local population and method principle. The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 30i. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):**
0.03 g/dl (0.3 g/l)
- **LoD ((Limit of Detection):**
0.05 g/dl (0.5 g/l)
- **LoQ (Limit of Quantitation):**
0.15 g/dl (1.5 g/l)
- **Linearity:**
up to 19 g/dl (190 g/l)

For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

- **Measurement range:**
0.15 g/dl (1.5 g/l) - 19 g/dl (190 g/l)

Specificity / Interferences

Haemoglobin up to 0.31 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 950 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [g/dl]	SD [g/dl]	CV [%]
level 1	4.54	0.03	0.59
level 2	6.88	0.06	0.87
Reproducibility (day to day) n = 80	Mean [g/dl]	SD [g/dl]	CV [%]
level 1	6.7	0.10	1.5
level 2	4.3	0.07	1.7

Method comparison

A comparison between total protein values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 63 serum samples gave following results:

$$y = 0.9766x + 0.258 \text{ g/dl};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

TRACEABILITY

2-STANDARD is traceable to the SRM 927D/1950 reference material.

WASTE MANAGEMENT ⁹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Residual reagents
1-TOTAL PROTEIN / 1-REAGENT - 18 01 06*
2-STANDARD - 18 01 03*
- Empty packages:
1-TOTAL PROTEIN / 1-REAGENT - 15 01 10*
2-STANDARD - 18 01 03*
- Wastewater from the analyzer:
1-TOTAL PROTEIN / 1-REAGENT - 18 01 03*
2-STANDARD/ 1-REAGENT - 18 01 03*

INCIDENTS ¹⁰

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (website address: incidents@cornay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

1. Szczeklik A., Gajewski P., ed. Interna Szczeklika. Mały podręcznik 2018/2019 10th ed., Kraków, Medycyna Praktyczna, (2018)
2. Alan H.B., Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B Saunders Company, 916 (2006)
3. Chowdhury, M., Akhter, N., Haque, M., Aziz, R., Nahar, N., Serum Total Protein and Albumin Levels in Different Grades of Protein Energy Malnutrition, *Journal of Bangladesh Society of Physiologist*, 3, 58-60 (2009).
4. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 396, 1780 (2018)
5. McPherson R.A., Pincus M.R., ed. Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed., ELSEVIER, Philadelphia, (2011).
6. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2:Jan 2002.
7. NCCLS - Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition, NCCLS document H18-A3 (ISBN 1-56238-555-0), NCCLS, 940 West Valley Road, Suite 1400 Wayne, Pennsylvania 19087-1898, USA (2004)
8. Thomas L.:Clinical Laboratory Diagnostics. 1st ed. Frankfurt: THBooks Verlagsgesellschaft; 644-647 (1998).
9. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
10. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 04	Current version: 05
Sections updated: <i>INTENDED USE; SUMMARY; CONCENTRATION OF THE ACTIVE INGREDIENTS; WARNING AND NOTES; EQUIPMENT REQUIRED, BUT NOT PROVIDED; SPECIMEN; ASSAY PROCEDURE; QUALITY CONTROL; REFERENCE VALUES; PERFORMANCE CHARACTERISTICS; TRACEABILITY; WASTE MANAGEMENT; LITERATURE</i>	
Sections added: <i>REAGENT STABILITY, INCIDENTS</i>	

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1. Szczeklik A., Gajewski P., ed. Interna Szczeklika. Mały podręcznik 2018/2019 10th ed., Kraków, Medycyna Praktyczna, (2018)
2. Alan H.B., Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B Saunders Company, 916 (2006)
3. Chowdhury, M., Akhter, N., Haque, M., Aziz, R., Nahar, N., Serum Total Protein and Albumin Levels in Different Grades of Protein Energy Malnutrition, *Journal of Bangladesh Society of Physiologist*, 3, 58-60 (2009).
4. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., St. Louis, Missouri, Elsevier, 396, 1780 (2018)
5. McPherson R.A., Pincus M.R., ed. Henry's Clinical Diagnosis and Management by Laboratory Methods, 22nd ed., ELSEVIER, Philadelphia, (2011).
6. WHO Publication: Use of anticoagulants in diagnostic laboratory investigations, WHO/DIL/LAB/99.1 Rev.2: Jan 2002.
7. NCCLS - Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition, NCCLS document H18-A3 (ISBN 1-56238-555-0), NCCLS, 940 West Valley Road, Suite 1400 Wayne, Pennsylvania 19087-1898, USA (2004)
8. Thomas L.:Clinical Laboratory Diagnostics. 1st ed. Frankfurt: THBooks Verlagsgesellschaft; 644-647 (1998).
9. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
10. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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