PHOSPHORUS

| | |
|----------------------------|---------------|
| Nazwa zestawu | (PL) |
| Liquick Cor-PHOSPHORUS 30 | Nr kat. 3-243 |
| Liquick Cor-PHOSPHORUS 120 | 3-338 |
| HC – PHOSPHORUS | 4-543 |
| OS – PHOSPHORUS | 9-415 |
| B50 – PHOSPHORUS | 5-534 |

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INTENDED USE

Diagnostic kit for determination of inorganic phosphorus concentration used both for manual assay and in several automatic analysers

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions

INTRODUCTION

Phosphorus is present in all body cells as a component of nucleic acids, phospholipids and phosphoproteins. Phosphorus is essential for intracellular storage and conversion of energy (ATP, creatine phosphate) and participates in carbohydrates metabolism. In the blood phosphorus is present as a mixture of inorganic phosphates HPO_4^{2-} and $H_2PO_4^-$. Besides phosphorus and calcium constitute mineral portion of bone. Continuous flux of phosphorus in organism is controlled by parathyroid hormone (PTH), vitamin D and calcitonin. Phosphorus serum level abnormalities are caused usually by disorders of vitamin D metabolism or parathyroid and kidney diseases.

METHOD PRINCIPLE

Direct phosphomolybdate reaction without deproteinization. Phosphate ions form with molybdate ions in acid solution proportional amounts of unreduced phosphomolybdate complex. The concentration of the complex formed is determined by measuring its absorbance.

REAGENTS

| | | | |
|----------------|----------------------------------|-----------------------------------|-------------------------|
| Package | Liquick Cor-PHOSPHORUS 30 | Liquick Cor-PHOSPHORUS 120 | |
| 1-PHOSPHORUS | 6 x 30 ml | 6 x 100 ml | |
| 2-STANDARD | 1 x 2 ml | - | |
| | HC – PHOSPHORUS | OS – PHOSPHORUS | B50 – PHOSPHORUS |
| 1-REAGENT | 6 x 97 ml | 3 x 53 ml | 2 x 58.5 ml |

2-STANDARD is phosphate ions standard solution: 1.6 mmol/l (5.0 mg/dl).

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. On board stability of the reagents depends on type of analyser used for analysis.

Concentrations in the test

| | |
|--------------------|---------------|
| ammonium molybdate | ≤ 0,48 mmol/l |
| sulphuric acid | ≤ 180 mmol/l |
| hydrochloric acid | ≤ 120 mmol/l |
| detergents | |

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Contaminated glassware is the greatest source of error. Disposable plastic ware is recommended for the test.
- The reagent is usable when its absorbance is less than 0.200 (read against distilled water, wavelength $\lambda=340$ nm, cuvette $l=1$ cm, at temp. 37°C).

- 1-PHOSPHORUS meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

1-PHOSPHORUS contains sulfuric acid (VI) and hydrochloric acid.

Danger

 H314 Causes severe skin burns and eye damage.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P301+P330+P331 IF SWALLOWED: rinse mouth.
Do NOT induce vomiting.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing.

P305 +P351 +P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P310 - Immediately call a POISON CENTER or doctor.

ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 340 nm;
- thermostat at 37°C;
- general laboratory equipment;

SPECIMEN

Serum, heparinized plasma (recommended: heparine lithium, sodium or ammonium salt) free from hemolysis, 24-hours urine.

Serum is the preferred specimen. Level of inorganic phosphate in heparinized plasma is about 0.2 to 0.3 mg/dl (0.06 – 0.10 mmol/l) lower than in serum.

Serum should be separated from red blood cells as soon as possible after blood collection, because erythrocytes contain several times higher phosphate concentration than normal serum.

Serum and plasma can be stored up to 7 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Urine preparation: to prevent phosphate precipitation in urine, specimens should be collected in HCl, 20-30 ml of 6 mol/l for 24-h specimen. Then dilute 1 part of acidified urine with 9 parts of distilled water. Multiply the result by the dilution factor.

24-hours urine samples can be stored up to 7 days at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

The reagent is ready to use.

Applications for analyzers are available on request.

Manual procedure

| | |
|-------------|--------|
| wavelength | 340 nm |
| temperature | 37°C |
| cuvette | 1 cm |

Pipette into the cuvette:

| | reagent blank (RB) | test (T) | standard (S) |
|--------------|--------------------|--------------|--------------|
| 1-PHOSPHORUS | 1000 μ l | 1000 μ l | 1000 μ l |

Bring up to the temperature of determination. Then add:

| | | | |
|----------|---|------------|------------|
| standard | - | - | 10 μ l |
| sample | - | 10 μ l | - |

Mix well, incubate for 5 min. at the determination temperature. Read the absorbance of test A(T) and standard A(S) against reagent blank (RB). The absorbance is stable within 30 minutes.

Calculation

$$\text{phosphorus concentration} = \frac{A(T)}{A(S)} \times \text{standard concentration}$$

REFERENCE VALUES

| serum / plasma | mg/dl | mmol/l |
|----------------|-----------|-------------|
| age: | | |
| 0 – 10 d | 4.5 – 9.0 | 1.45 – 2.91 |
| 10 d – 24 mo | 4.5 – 6.7 | 1.45 – 2.16 |
| 24 mo – 12 y | 4.5 – 5.5 | 1.45 – 1.78 |
| 12 – 60 y | 2.7 – 4.5 | 0.87 – 1.45 |
| > 60 y males | 2.3 – 3.7 | 0.74 – 1.20 |
| > 60 y females | 2.8 – 4.1 | 0.90 – 1.32 |
| 24-hours urine | g/24h | mmol/24h |
| | 0.4 – 1.3 | 12.9 – 42.0 |

It is recommended for each laboratory to establish its own reference ranges for local population.

Phosphorus concentration in 24-hours urine – calculation

| | | | | | |
|--|--|---|----------|--------|------|
| phosphorus concentration in 24-hours urine [g/24h] | phosphorus concentration in sample of 24-hours urine [mg/dl] | urine volume of 24-hours urine [dl/24h] | \times | \div | 1000 |
|--|--|---|----------|--------|------|

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum, CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine. For the calibration of manual assay the PHOSPHORUS STANDARD (Cat. No 5-120) is recommended.

For calibration of the automatic analysers systems CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

Calibration stability depends on type of analyser used for analysis. The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using a Multi + analyser for manual assay and/or an automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):**
0.02 mg/dl (0.006 mmol/l) - Biolis 24i Premium
- LoD (Limit of Detection):**
0.05 mg/dl (0.016 mmol/l) - Biolis 24i Premium
- LoQ (Limit of Quantitation):**
0.4 mg/dl (0.129 mmol/l) - Multi +
0.12 mg/dl (0.039 mmol/l) - Biolis 24i Premium
- Linearity:**
up to 18 mg/dl (5.8 mmol/l) – Multi +
up to 20 mg/dl (6.46 mmol/l) - Biolis 24i Premium

For higher concentration dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 0.12 g/dl, ascorbate up to 62 mg/l, bilirubin up to 20 mg/dl and triglycerides up to 1000 mg/dl do not interfere with the test.

Precision (Multi +)

| Repeatability (run to run) n = 20 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|-----------------------------------|--------------|------------|--------|
| level 1 | 3.24 | 0.11 | 3.4 |
| level 2 | 7.18 | 0.20 | 2.8 |

Precision (Biolis 24i Premium)

| Repeatability (run to run) n = 20 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
|-------------------------------------|--------------|------------|--------|
| level 1 | 3.46 | 0.04 | 1.1 |
| level 2 | 7.31 | 0.07 | 1.0 |
| Reproducibility (day to day) n = 80 | Mean [mg/dl] | SD [mg/dl] | CV [%] |
| level 1 | 3.37 | 0.06 | 1.9 |
| level 2 | 7.19 | 0.12 | 1.7 |

Method comparison

A comparison between phosphorus values determined at Multi + (y) and at BECKMAN COULTER AU680 (x) using 30 serum samples gave following results:

$$y = 1.03 x - 0.0731 \text{ mg/dl;}$$

$$R = 0.995 \quad (R - \text{correlation coefficient})$$

A comparison between phosphorus values determined at Biolis 24i Premium (y) and at BECKMAN COULTER AU680 (x) using 155 serum samples gave following results:

$$y = 0.9883 x + 0.0965 \text{ mg/dl;}$$

$$R = 0.997 \quad (R - \text{correlation coefficient})$$

A comparison between phosphorus values determined at Biolis 24i Premium (y) and at BECKMAN COULTER AU680 (x) using 57 urine samples gave following results:

$$y = 0.9613 x + 0.7032 \text{ mg/dl;}$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

TRACEABILITY

PHOSPHORUS STANDARD is traceable to the spectrophotometric reference method.

WASTE MANAGEMENT

Please refer to local legal requirements.

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