ASAT

Nazwa zestawu	(PL) Nr kat.	HC-ASAT	OS-ASAT	B50-ASAT
Liquick Cor-ASAT 30	1-222	1-REAGENT 6 x 76 ml	4 x 53,5 ml	2 x 58,5 ml
Liquick Cor-ASAT 60	1-214	2-REAGENT 6 x 19,5 ml	4 x 16 ml	2 x 18,4 ml
Liquick Cor-ASAT 120	1-215			
HC-ASAT	4-514			
OS-ASAT	9-418			
B50-ASAT	5-522			

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INTENDED USE


ASAT reagent is intended to determine quantitatively aspartate aminotransferase (ASAT) activity in serum. It is used for monitoring and as an aid to diagnosis of clinical conditions associated with abnormal aspartate aminotransferase activity. ASAT reagent is intended to use on automatic analysers and manually. It is only for *in vitro* diagnostics, for healthcare professional users.

INTRODUCTION^{1,2}

Determinations of elevated aspartate aminotransferase activity (ASAT, AST, GOT) in serum are mainly used to aid in the diagnosis and monitor of liver diseases e.g. hepatitis, hepatic cirrhosis, drug-induced liver injury, hepatic metastasis, mononucleosis. Also skeletal muscle diseases such as stroke trauma, surgery, burns, muscular dystrophy, heat stroke and other diseases such as acute hemolytic anemia, acute pancreatitis and patients with higher risk of COVID-19 adverse outcome. Decreased levels of ASAT can be observed in acute renal disease, beriberi, diabetic ketoacidosis, pregnancy, chronic renal dialysis.

METHOD PRINCIPLE^{3,4}

Optimized, modified method according to International Federation of Clinical Chemistry (IFCC), without pyridoxal phosphate.

The rate of absorbance changing at $\lambda=340$ nm is directly proportional to aspartate aminotransferase activity.

REAGENTS

Package	Liquick Cor - ASAT 30	Liquick Cor - ASAT 60	Liquick Cor - ASAT 120
1-ASAT	5 x 24 ml	5 x 48 ml	5 x 96 ml
2-ASAT	1 x 30 ml	1 x 60 ml	1 x 120 ml
	HC-ASAT	OS-ASAT	B50-ASAT
1-REAGENT	6 x 76 ml	4 x 53.5 ml	2 x 58.5 ml
2-REAGENT	6 x 19.5 ml	4 x 16 ml	2 x 18.4 ml

CONCENTRATION OF THE ACTIVE INGREDIENTS

1-ASAT / 1-REAGENT	
L-aspartate	300 mmol/l
LDH	1.95 U/ml
MDH	1.15 U/ml
buffer Tris	
pH adjuster	
stabilizers	
preservative	
2-ASAT / 2-REAGENT	
2-oxoglutarate	62 mmol/l
NADH	1.4 mmol/l
buffer	
pH adjuster	
preservatives	

REAGENT STABILITY

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents stored on board of the analyzer at 2-10°C are stable for 12 weeks (Biolis 30i).

Working reagent preparation and stability

Assay can be performed with use of separate 1-ASAT and 2-ASAT reagents (1-REAGENT and 2-REAGENT) or with use of working reagent. For working reagent preparation mix gently 4 parts of 1-ASAT (1-REAGENT) with 1 part of 2-ASAT (2-REAGENT). Avoid foaming!

Stability of working reagent: 4 weeks at 2-8°C
5 days at 15-25°C

WARNINGS AND NOTES

- Protect from direct sunlight and avoid contamination!
- The reagents are usable when the absorbance of the working reagent is higher than 1.400 (read against distilled water, wavelength $\lambda=340$ nm, cuvette l = 1 cm, at temp. 25°C).
- Do not use reagent beyond expiry date printed on the package.
- Do not mix reagents from different kits or lots.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- 1-ASAT/ 1-REAGENT meeting the criteria for classification in accordance with Regulation (EC) No 1272/2008

Warning

H315 Causes skin irritation.
H319 Causes serious eye irritation.
P280 Wear protective gloves/protective clothing/eye protection/face protection.

P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

EQUIPMENT REQUIRED BUT NOT PROVIDED

- automatic analyzer or photometer able to read at 340 nm (Hg 334nm, 365nm);
- thermostat at 37°C;
- general laboratory equipment;

SPECIMEN^{5,6,7,8}

Serum free from hemolysis.
Hemolysis should be avoided, since ASAT activity in erythrocytes is 15 times higher than in normal serum. ASAT activity remains stable in specimen up to 4 days at 15-25°C, up to 7 days at 2-8°C or up to 3 months at 2-8°C. Nevertheless it is recommended to perform the assay with freshly collected samples!
Follow tube manufacturers' instructions carefully when using collection tubes.


Human-origin material should be handled as potentially infectious. Standard precautions in normal laboratory work are required.

ASSAY PROCEDURE

Reagents are ready to use.

Carry-over effect may occur when the reagent is used on automated chemistry analyzers. It is based on interference between specific test combinations and results in either underestimating or inflating patient sample results. To minimize this effect, the preventive actions should be taken, such as ordering an extra washing program or performing tests in a separate order.

General procedure for automatic determinations:

Main Wavelength: 340 nm

2nd wavelength: 450 nm

Method: Kinetic

Calibration Type: Linear

Direction: Decreasing

R1 Volume: 200 µl

Sample Volume: 18 µl

R2 Volume: 50 µl

Application guides for analysers are available on request.

The application guide should be verified and validated by the user prior to testing patient samples.

Manual procedure

wavelength	340 nm (Hg 334 nm, 365 nm)
temperature	37°C
cuvette	1 cm

Sample Start method

Pipette into the cuvette:

working reagent	1000 µl
sample	100 µl

Bring up to the temperature of determination 37°C for 10 minutes. Then add:

sample	100 µl
--------	--------

Mix and incubate at adequate temperature. After about 1 min. read the absorbance against air or water. Repeat the reading after exactly 1, 2 and 3 minutes. Calculate the mean absorbance change per minute ($\Delta A/\text{min.}$).

Calculation

ASAT activity [U/l] = $\Delta A/\text{min.}$ x calibrator concentration

or

ASAT activity [U/l] = $\Delta A/\text{min.}$ x F

F value depends on the used wavelength:

λ	334 nm	340 nm	365 nm
F	1979	1939	4800

Reagent Start method

The determination can be also performed with use of separate 1-ASAT and 2-ASAT reagents (1-REAGENT and 2-REAGENT).

Pipette into the cuvette:

1-ASAT/ 1-REAGENT	1000 µl
sample	100 µl

Bring up to the temperature of determination. Then add:

2-ASAT/ 2-REAGENT	250 µl
-------------------	--------

Mix well, incubate for 5 min. Then add:

2-ASAT/ 2-REAGENT	250 µl
-------------------	--------

Mix well; perform measurement as described for Sample Start method.

Calculation

ASAT activity [U/l] = $\Delta A/\text{min.}$ x calibrator concentration

or

ASAT activity [U/l] = $\Delta A/\text{min.}$ x F

F value depends on the used wavelength:

λ	334 nm	340 nm	365 nm
F	2670	2441	5922

QUALITY CONTROL

For internal quality control it is recommended to use with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173). For the calibration of manual assay the CORMAY MULTICALIBRATOR LEVEL 2 (Cat. No 5-175; 5-177) or appropriate assigned factor (F) is recommended.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) is recommended. Deionised water should be used as a calibrator 0.

Calibration stability depends on the type of analyser used for analysis. The calibration curve should be prepared every 12 weeks (Biolis 30i).

Calibration is recommended in the following cases:

- after each change of lot,
- after instrument service,

- if controls lie outside the expected range,
- each time a new reagent kit is used.

If the quality control results do not fall within the expected values or within the range determined in the laboratory, despite a successful calibration procedure, do not report results. In this case, please take the following actions:

- verify reagents are not out of expiration date.
- verify that the required maintenance has been carried out.
- verify that the procedure has been performed in accordance with the instructions for use.
- contact the Service Department or distributor for assistance.

REFERENCE VALUES ⁹

serum	37°C	
female	up to 31 U/l	up to 0.53 µkat/l
male	up to 35 U/l	up to 0.60 µkat/l

It is recommended for each laboratory to establish its own reference ranges for local population and method principle.

The assay results should be used in conjunction with other data, such as symptoms, results of other tests and clinical history to make clinical decisions. It is not recommended to make clinical diagnosis based on a single result.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using Multi+ for manual assay and automatic analyser Biolis 30i and/or BS-400. Results may vary if a different instrument or a manual procedure is used.

- **LoB (Limit of Blank):**
1.8 U/l (0.03 µkat/l) – BS-400
- **LoD (Limit of Detection):**
3.0 U/l (0.05 µkat/l) – BS-400
- **LoQ (Limit of Quantitation):**
9.0 U/l (0.15 µkat/l) – Biolis 30i
7.0 U/l (0.12 µkat/l) – Multi+
- **Linearity:**
up to 770 U/l (12.8 µkat/l) – Biolis 30i
up to 650 U/l (10,8 µkat/l) – Multi +
- **Measurement range**
9,0 U/l (0,15 µkat/l) - 770 U/l (12,8 µkat/l) – Biolis 30i
7,0 U/l (0,12 µkat/l) – 650 U/l (10,8 µkat/l) – Multi +

▪ Specificity / Interferences

Ascorbic acid up to 62 mg/l, triglycerides up to 1000 mg/dl, bilirubin up to 20 mg/dl, haemoglobin up to 0.31 g/dl (in samples containing low activity of ASAT) and up to 2.5 g/dl (in samples containing high activity of ASAT) do not interfere with the test.

▪ Precision (Biolis 30i)

Repeatability (run to run) n = 20	Mean [U/l]	SD [U/l]	CV [%]
level 1	44.6	1.56	3.5
level 2	200.4	1.33	0.7

Reproducibility (day to day) n = 80	Mean [U/l]	SD [U/l]	CV [%]
level 1	40.6	1.36	3.4
level 2	210.4	1.99	0.9

▪ Method comparison

A comparison between ASAT values determined at **Biolis 30i** (y) and at **ADVIA SIEMENS 1800** (x) using 59 serum samples gave following results:

$$y = 1.0272x - 1.1659 \text{ U/l};$$

$$R = 0.996 \quad (R - \text{correlation coefficient})$$

A comparison between ASAT values determined at **Multi+** (y) and at **BECKMAN COULTER AU680** (x) using 22 serum samples gave following results:

Sample Start method

$$y = 0.8982 + 2.524 \text{ U/l};$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

Reagent Start method

$$y = 0.9081 + 0.8803 \text{ U/l};$$

$$R = 0.999 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT ¹⁰

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 07
- Empty packages: 15 01 02
- Wastewater from the analyzer: 18 01 03*

INCIDENTS ¹¹

Any serious incident that has occurred in relations to the device shall be reported to the manufacturer (website address: incidents@cornay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

Serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:

- the death of a patient, user or other person,
- the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
- a serious public health threat.

LITERATURE

1. Pagana K. Pagana T.J., Mosby's Diagnostic and Laboratory Test Reference, 10th ed., Elsevier, St. Louis, USA, 27 (2013).
2. Lippi G., Plebani M., Laboratory abnormalities in patients with COVID-2019 infection, Clin Chem Lab Med, 58(7):1131-1134 (2020).
3. Bergmeyer H.U., IFCC methods for the measurement of catalytic concentrations of enzymes: Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2), Clin. Chem. Acta 105, 147-172 (1980).
4. Bergmeyer HU, Horder M, Rej R. International Federation of Clinical Chemistry (IFCC) Scientific Committee, Analytical Section: approved recommendation (1985) on IFCC methods

for the measurement of catalytic concentration of enzymes. Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2). J Clin Chem Clin Biochem. 1986 Jul;24(7):481-95. PMID: 3734711.

5. Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 911 (1989).
6. NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition. NCCLS document H18-A3 (ISBN 1-56238-555-0). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA (2004).
7. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 413, 1800 (2018).
8. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64
9. Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2257, (2006).
10. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
11. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

LIST OF CHANGES

Previous version: 04	Current version: 05
Sections updated: <i>INTENDED USE; SUMMARY; CONCENTRATION OF THE ACTIVE INGREDIENTS; WARNING AND NOTES; EQUIPMENT REQUIRED, BUT NOT PROVIDED; SPECIMEN; ASSAY PROCEDURE; QUALITY CONTROL; REFERENCE VALUES; PERFORMANCE CHARACTERISTICS; WASTE MANAGEMENT; LITERATURE</i>	
Sections added: <i>REAGENT STABILITY, INCIDENTS</i>	

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- Pagana K. Pagana T.J., Mosby's Diagnostic and Laboratory Test Reference, 10th ed., Elsevier, St. Louis, USA, 27 (2013).
- Lippi G., Plebani M., Laboratory abnormalities in patients with COVID-2019 infection, Clin Chem Lab Med, 58(7):1131-1134 (2020).
- Bergmeyer H.U., IFCC methods for the measurement of catalytic concentrations of enzymes: Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2), Clin. Chem. Acta 105, 147-172 (1980).
- Bergmeyer HU, Horder M, Rej R. International Federation of Clinical Chemistry (IFCC) Scientific Committee, Analytical Section: approved recommendation (1985) on IFCC methods for the measurement of catalytic concentration of enzymes. Part 3. IFCC method for alanine aminotransferase (L-alanine: 2-oxoglutarate aminotransferase, EC 2.6.1.2). J Clin Chem Clin Biochem. 1986 Jul;24(7):481-95. PMID: 3734711.
- Kaplan L.A., Pesce A.J., Clinical Chemistry: Theory, analysis, and correlation, 2nd ed., The C.V. Mosby Company, St. Louis, USA, 911 (1989).

6. NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline - Third Edition. NCCLS document H18-A3 (ISBN 1-56238-555-0). NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA (2004).

7. Rifai N., Horvath A.R., Wittwer C., ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics 6th ed., Elsevier, St. Louis, USA, 413, 1800 (2018).

8. World Health Organization (2002). Use of Anticoagulants in Diagnostic Laboratory Investigations (Stability of Blood, Plasma and serum samples); WHO/DIL/LAB/99.1 Rev.2: 1-64

9. Burtis C.A., Ashwood E.R., Bruns D.E. ed. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Elsevier Saunders, St. Louis, USA, 2257, (2006).

10. European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.

11. Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices.

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