ASO

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INTENDED USE

Diagnostic kit for determination of anti-streptolysin O levels used both for manual assay and in several automatic analysers. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions

INTRODUCTION

Most people infected with hemolytic streptococcus produce anti-streptolysin O (ASO), antibodies against streptolysin O (SLO), an exotoxin of Streptococcus. Measuring the level of ASO is effective for diagnosing, judging the progress of medical treatment, and assessing recovery from diseases caused by hemolytic streptococcus such as rheumatic fever, acute glomerulonephritis, scarlatina and tonsillitis.

METHOD PRINCIPLE

When an antigen-antibody reaction occurs between ASO in a sample and SLO which has been sensitized to latex particles, agglutination results. This agglutination is detected as an absorbance change (572 nm), with the magnitude of the change being proportional to the quantity of ASO in the sample. The actual concentration is then determined by interpolation from a calibration curve prepared from calibrators of known concentration.

REAGENTS

Package	CORMAY ASO	CORMAY ASO 60
1-REAGENT	1 x 28,5 ml	2 x 53ml
2-REAGENT	1 x 46 ml	3 x 60 ml
	OS-ASO	B50-ASO
1-REAGENT	3 x 19 ml	2 x 11.5 ml
2-REAGENT	3 x 31 ml	2 x 18.5 ml

The reagent is stable up to the kit expiry date printed on the package when stored at 2-8°C. On board stability of the reagents depends on type of analyser used for analysis.

Concentrations in the test

suspension of latex particles sensitized with SLO (pH 8.2) 0.17 w/v%
 glycine buffer solution (pH 8.3)
 preservative

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Reagent bottles should be shaken before use by gently inverting several times.
- After measurements are taken, reagent bottles should be capped and kept at 2-10°C. Care should be taken not to interchange the caps of reagent bottles.

- Reagents with different lot numbers should not be interchanged or mixed.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product

SPECIMEN

Serum or plasma (Na-EDTA, K-EDTA, Na-Heparin, Li-Heparin, citric acid). After blood has clotted thoroughly, the sample is centrifuged and the serum is separated from blood cells and fibrins.

If the test cannot be done immediately, the sample should be placed in a tightly sealable container and stored at -20°C. Repeated freezing and thawing should be avoided. Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

The reagents are ready to use.
 Applications for analyzers are available on request.

REFERENCE VALUES

serum, plasma	< 160 IU/ml
---------------	-------------

It is recommended for each laboratory to establish its own reference ranges for local population. Diagnosis should only be made after taking clinical symptoms and the results of other tests into consideration.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL I (Cat. No 4-288) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY ASO CALIBRATOR kit (Cat. No 4-278) is recommended. Calibrator and 0.9% NaCl should be used for calibration.

Calibration stability depends on type of analyser used for analysis. The calibration curve should be prepared with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using the automatic analysers Hitachi 912 and Hitachi 917. Results may vary if a different instrument is used.

- Analytical range:** 20 – 1200 IU/ml.
For higher concentrations dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.
- Specificity / Interferences**
Haemoglobin up to 0.5 g/dl, bilirubin up to 20 mg/dl and triglycerides up to 500 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	55.62	2.79	5.02
level 2	288.91	2.34	0.81

Reproducibility (day to day) n = 12	Mean [IU/ml]	SD [IU/ml]	CV [%]
level 1	48.8	2.83	5.81
level 2	78.8	2.60	3.30
level 3	219.8	5.24	2.38

Method comparison

A comparison between ASO values determined at Hitachi 912 (y) and at Advia 1650 (x) using 37 samples gave following results:

$$y = 0.9926x - 7.5251 \text{ IU/ml};$$

$$R = 0.9953 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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