

DIAGNOSTIC KIT FOR DETERMINATION OF HAPTOGLOBIN CONCENTRATION



OS – HAPTOGLOBIN

INTRODUCTION

Haptoglobin is an acute phase protein whose primary function consists in binding free haemoglobin in serum. The complex is removed within minutes by the reticulo-endothelial system where its components are metabolised to free aminoacids and iron. Haptoglobin consequently plays a major role in preventing the loss of haemoglobin in urine and the consequent iron loss from the iron pool. The haptoglobin levels are increased during the acute phase and in such conditions as burns or nephrotic syndrome. Haptoglobin levels are abnormally high in intravascular hemolysis and when haemoglobin turnover is increased such as during haemolytic anaemia, transfusion reactions and malaria.

METHOD PRINCIPLE

The haptoglobin presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum measured at $\lambda=340$ nm is proportional to haptoglobin concentration in the sample.

REAGENTS

Package

- 1-Reagent 1 x 35.5 ml
- 2-Reagent 1 x 9 ml

Buffer (1-Reagent) stored at 2-25°C and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package. Store closed and avoid contamination.

Concentrations in the test

Imidazole buffer (pH 7.0); PEG; sodium chloride; anti human haptoglobin antiserum; HEPES buffer (pH 7.4); sodium azide (< 1 g/l); stabilizers.

Warnings and notes

- Products for in vitro diagnostic use only.
- The reagents must be used only for the intended purpose, by suitably qualified laboratory personnel, under appropriate laboratory conditions.
- Products from human source have been tested for the HIV antibody, HbsAg and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Products contain sodium azide (< 1 g/l) as a preservative. Avoid contact with skin and mucous membranes.

SPECIMEN

Serum.
Sample may be stored several days at 2-8°C. Samples frozen at -20°C can be stored longer.
Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

These reagents may be used in automatic analysers Olympus AU400/AU640.

1-Reagent and 2-Reagent are ready to use.
For reagent blank 0.9% NaCl is recommended.

APPLICATION

Reagent ID: 065

Specific Test Parameters										
General		LIH	ISE	Range						
Test name:		HAPTO				Type:	Serum		Operation: Yes	
Sample: Volume	7	μL	Dilution	0	μL	Pre-Dilution Rate:	1			
Reagents: R1 Volume	160	μL	Dilution	0	μL	Min OD	Max OD			
R2 Volume	32	μL	Dilution	0	μL	L	-2.0000	H	2.5000	
Wavelength: Pri.						340	Sec.	700		
Method:						END				
Reaction Slope:						+				
Measuring Point 1: First			0	Last		27				
Measuring Point 2: First			0	Last		10				
Linearity:								Correlation Factor:		
No-Lag-Time:								A		1.000
								B		0.000
								On-board Stability Period:		

Specific Test Parameters											
General		LIH	ISE	Range							
Test name:		HAPTO				Type:	Serum				
Value/Flag:		#		Level L:		#		Level H:		#	
Normal Ranges:											
		Sex	Age L	Age H			L	H			
			Year	Month	Year	Month					
1.	#	#	#	#	#	#	#	#			
2.	#	#	#	#	#	#	#	#			
3.	#	#	#	#	#	#	#	#			
4.	#	#	#	#	#	#	#	#			
5.	#	#	#	#	#	#	#	#			
6.	#	#	#	#	#	#	#	#			
7. None Selected							#	#			
8. Out of Range							#	#			
							L	H			
Panic Value:		#		#		Unit:	g/l		Decimal Places: 2		

Calibration Specific											
General		ISE									
Test name:		HAPTO				Type:	Serum				
Calibration Type:		6AB		Formula:		Spline		Counts: 1		Process: CONC	
	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H						
Point 1:	#		**	-2.0000	2.5000						
Point 2:	#		*	-2.0000	2.5000						
Point 3:	#		*	-2.0000	2.5000						
Point 4:	#		*	-2.0000	2.5000						
Point 5:	#		*	-2.0000	2.5000						
Point 6:	#		*	-2.0000	2.5000						
Point 7:	#		*	-2.0000	2.5000						
1-Point Cal.Point:				with CONC=0		Slope Check:		None		Advanced Calibration: #	
MB Type Factor:				Calibration Stability Period:							

- # User defined
- * Calibrator value
- ** Saline should be used as calibrator 1

REFERENCE VALUES ⁴

adults	0.26 – 1.85 g/l
newborns	0.05 – 0.48 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL III (Cat. No 4-291) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended. The calibration curve should be prepared every 4 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using an automatic analyser Cobas Mira. Results may vary if a different instrument or manual procedure is used.

▪ **Analytical range:** 0.05 g/l do 4 g/l.

▪ **Interferences:**

Hemoglobin up to 0.32 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 2000 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

▪ **Precision**

Repeatability (run to run) n = 10	Mean [mg/dl]	SD	CV [%]
level 1	55.9	0.7	1.3
level 2	110.8	1.2	1.1
level 3	137.7	1.5	1.1

Reproducibility (day to day) n = 10	Mean [mg/dl]	SD	CV [%]
level 1	60.0	2.1	3.6
level 2	113.0	4.6	4.1
level 3	141.1	4.9	3.5

▪ **Method comparison**

A comparison between CORMAY reagent (y) and commercially available assay (x) using 68 samples gave following results:

$$y = 0.86 x + 11.1 \text{ mg/dl};$$

$$R = 0.946 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Kaplan L.A., Pesce A.J.: Clinical Chemistry, Third Edition, Mosby, 731 (1996).
2. Jacobs, D. S. et al., Laboratory test Handbook, Mosby, St Louis, (1984).
3. Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, (1994).
4. Alan H.B. Wu, ed.: Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B. Saunders Company., 512, (2006).

Date of issue: 05. 2016.

MANUFACTURER

PZ CORMAY S.A.
22 Wiosenna Street,
05-092 Łomianki, POLAND
tel.: +48 (0) 22 751 79 10
fax: +48 (0) 22 751 79 14
<http://www.cormay.pl>

05/16/05/16