



**DIAGNOSTIC KIT  
FOR DETERMINATION OF  
α 1-ANTITRYPSIN CONCENTRATION**

**OS – ALPHA 1-ANTITRYPSIN**

**INTRODUCTION**

α 1-antitrypsin is an acute phase protein showing anti-protease activity. The main function of this protein consists in neutralising lysosomal elastase released upon phagocytosis by polymorphonuclear leukocytes. Inherited deficiency of the protein is associated with lung and liver diseases. Low levels are encountered in neonatal respiratory distress syndrome and in severe protein losing disorders. Increased levels are more common, particularly during the acute phase.

**METHOD PRINCIPLE**

The α 1-antitrypsin presents in a sample form with the specific antibody an immunological complex. The increase of turbidity after the addition of antiserum is proportional to α 1-antitrypsin concentration in the sample.

**REAGENTS**

**Package**

- 1-Reagent 1 x 45.5 ml
- 2-Reagent 1 x 9 ml

Buffer (1-Reagent) and antiserum (2-Reagent) stored at 2-8°C are stable until expiry date printed on the package. Protect from light and avoid contamination!

**Concentrations in the test**

TRIS buffer (pH 8); PEG; sodium chloride; anti human α 1-antitrypsin antiserum; HEPES buffer (pH 7.4); stabilizers; detergents.

**Warnings and notes**

- Products for in vitro diagnostic use only.
- The reagents must be used only for the intended purpose, by suitably qualified laboratory personnel, under appropriate laboratory conditions.
- Products from human source have been tested for the HBsAg and antibodies to HIV and HCV and found to be non-reactive. However this material should be handled as though capable of transmitting infectious disease.
- Products contain sodium azide (< 1 g/l) as a preservative. Avoid contact with skin and mucous membranes.

**SPECIMEN**

Serum or plasma (EDTA).  
Sample may be stored up to 24 hours at 20-25°C, up to 7 days at 4-8°C or up to 3 months at -70°C.  
Nevertheless it is recommended to perform the assay with freshly collected samples!

**PROCEDURE**

These reagents may be used in automatic analysers Olympus AU400/AU640.  
1-Reagent and 2-Reagent are ready to use.  
For reagent blank 0.9% NaCl is recommended.

**APPLICATION**

Reagent ID: 063

Specific Test Parameters									
General		LIH	ISE	Range					
Test name:	ATT					Type:	Serum	Operation:	Yes
Sample: Volume	2	μL	Dilution	0	μL	Pre-Dilution Rate:	1		
Reagents: R1 Volume	210	μL	Dilution	0	μL	Min OD	Max OD		
R2 Volume	30	μL	Dilution	0	μL	L	-2.0000	H	2.5000
						Reagent OD Limit:			
Wavelength: Pri.	340		Sec.	700		First L	-2.0000	First H	2.5000
Method:	END					Last L	-2.0000	Last H	2.5000
Reaction Slope:	+					Dynamic Range:			
Measuring Point 1: First	0		Last	27		L		H	
Measuring Point 2: First	0		Last	10		Correlation Factor:			
Linearity:		%				A	1.000	B	0.000
No-Lag-Time:						On-board Stability Period:			

Specific Test Parameters									
General		LIH	ISE	Range					
Test name:	ATT					Type:	Serum		
Value/Flag:	#		Level L:	#		Level H:	#		
Normal Ranges:									
	Sex	Age L		Age H		L		H	
		Year	Month	Year	Month				
1.	#	#	#	#	#	#	#	#	#
2.	#	#	#	#	#	#	#	#	#
3.	#	#	#	#	#	#	#	#	#
4.	#	#	#	#	#	#	#	#	#
5.	#	#	#	#	#	#	#	#	#
6.	#	#	#	#	#	#	#	#	#
7.	None Selected					#	#	#	#
8.	Out of Range					#	#	#	#
						L		H	
Panic Value:	#		#			Unit:	g/l	Decimal Places:	2

Calibration Specific									
General		ISE							
Test name:	ATT					Type:	Serum		
Calibration Type:	6AB		Formula:	Spline		Counts:	1	Process:	CONC
	Cal. No.	OD	CONC	Factor/OD-L	Factor/OD-H				
Point 1:	#		**	-2.0000	2.5000				
Point 2:	#		*	-2.0000	2.5000				
Point 3:	#		*	-2.0000	2.5000				
Point 4:	#		*	-2.0000	2.5000				
Point 5:	#		*	-2.0000	2.5000				
Point 6:	#		*	-2.0000	2.5000				
Point 7:	#		*	-2.0000	2.5000				
1-Point Cal.Point:	<input type="checkbox"/>	<input type="checkbox"/>	with CONC=0	Slope Check:	None	Advanced Calibration:	<input type="checkbox"/>		
MB Type Factor:			Calibration Stability Period:						

- # User defined
- \* Calibrator value
- \*\* Saline should be used as calibrator 1

**REFERENCE VALUES<sup>3</sup>**

adults	0.78 – 2.00 g/l
newborns	1.45 – 2.70 g/l

It is recommended for each laboratory to establish its own reference ranges for local population.

**QUALITY CONTROL**

For internal quality control it is recommended to use the CORMAY IMMUNO-CONTROL III (Cat. No 4-291) with each batch of samples.

For the calibration of automatic analysers systems the CORMAY IMMUNO-MULTICAL (Cat. No 4-287) is recommended.

The calibration curve should be prepared every 4 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

## PERFORMANCE CHARACTERISTICS

These metrological characteristics have been obtained using the automatic analyser Hitachi 917. Results may vary if a different instrument or manual procedure is used.

▪ **Measurement range:** 0.007 g/l to 3.50 g/l.

▪ **Interferences:**

Hemoglobin up to 0.322 g/dl, bilirubin up to 29.5 mg/dl, triglycerides up to 1682 mg/dl, heparin up to 0.5 g/l, sodium fluoride up to 4 g/l, EDTA up to 5 g/l, sodium citrate up to 5 g/l do not interfere with the test.

▪ **Precision**

Repeatability (run to run) n = 30	Mean [mg/dl]	SD	CV [%]
level 1	58.3	0.7	1.1
level 2	117.7	1.0	0.9
level 3	176.7	2.0	1.1

▪ **Method comparison**

A comparison between CORMAY reagent (y) and commercially available assay (x) using 40 samples gave following results:

$y = 0.9668x + 3.6562$  mg/dl;

$R = 0.9786$  (R – correlation coefficient)

## WASTE MANAGEMENT

Please refer to local legal requirements.

## LITERATURE

1. Pesce A.J., Kaplan L.A.: Methods in Clinical Chemistry., 962-965, (1996).
2. Bergstrom, K. & Lefvert, A.K.: Scand.J.clin.Lab.Invest. 40 (1980) 637.
3. Alan H.B. Wu, ed.: Tietz Clinical Guide to Laboratory Tests, 4th ed. W.B. Saunders Company., 138, (2006).
4. Tietz Textbook of Clinical Chemistry, W.B. Saunders, Philadelphia, (1994).
5. Roitt, I.: Essential Immunology, Blackwell, Oxford, (1991).
6. Use of Anticoagulants in Diagnostic Laboratory Investigations & Stability of blood, plasma and serum samples. Publication WHO/DIL/LAB/99.1 Rev. 2. Jan. 2002.

**Date of issue:** 05. 2016.

## MANUFACTURER

**PZ CORMAY S.A.**  
22 Wiosenna Street,  
05-092 Łomianki, POLAND  
tel.: +48 (0) 22 751 79 10  
fax: +48 (0) 22 751 79 14  
<http://www.cormay.pl>

05/16/05/16