CALCIUM ARSENAZO

Nazwa zestawu	(PL) Nr kat.
Liquick Cor-CALCIUM ARSENAZO 30	3-254
Liquick Cor-CALCIUM ARSENAZO 60	3-255
Liquick Cor-CALCIUM ARSENAZO 120	3-256
HC- CALCIUM ARSENAZO	4-547
OS- CALCIUM ARSENAZO	9-456
B50- CALCIUM ARSENAZO	5-537

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INTENDED USE

Diagnostic kit for determination of calcium concentration used for manual assay and in several automatic analysers.

The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Calcium plays an essential role in many cell functions: intracellularly in muscle contraction and glycogen metabolism, extracellularly, in bone mineralization, in blood coagulation and in transmission of nerve impulses. Calcium is present in plasma in three forms: free, bound to proteins or complexed with anions as phosphate, citrate and bicarbonate. Decreased total calcium levels can be associated with diseases of the bone apparatus (especially osteoporosis), kidney diseases (especially under dialysis), defective intestinal absorption and hypoparathyroidism. Increased total calcium can be measured in hyperparathyroidism, malignant diseases with metastases and sarcoidosis. Calcium measurements also help in monitoring of calcium supplementation mainly in the prevention of osteoporosis.

METHOD PRINCIPLE

Photometric test using arsenazo III.

Calcium with arsenazo III at neutral pH yields a blue colored complex, whose intensity is proportional to the calcium concentration. Interference by magnesium is eliminated by addition of 8-hydroxyquinoline-5-sulfonic acid.

REAGENTS

Package	Liquick Cor-CALCIUM ARSENAZO 30	Liquick Cor-CALCIUM ARSENAZO 60	Liquick Cor-CALCIUM ARSENAZO 120
1-CALCIUM ARSENAZO	6 x 30 ml	6 x 60 ml	6 x 120 ml
	HC-CALCIUM ARSENAZO	OS-CALCIUM ARSENAZO	B50-CALCIUM ARSENAZO
1-CALCIUM ARSENAZO	2 x 97.5 ml	6 x 53 ml	4 x 58.5 ml

The reagent, stored at 2-8°C is stable up to expiry date printed on the package.

Concentrations in the test

phosphate buffer (pH 7.5)	≤ 80 mmol/l
8-hydroxyquinoline-5-sulfonic acid	≤ 6 mmol/l
arsenazo III	≤ 120 μmol/l
detergents, preservative	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Do not freeze the reagent.
- Contaminated glassware is the greatest source of error. The use of disposable plastic ware is recommended. Glassware should be soaked for a few hours in 2M HCl solution and then thoroughly rinsed with distilled water.

ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 650 nm (630-670 nm);
- thermostat at 37°C;
- general laboratory equipment;

SPECIMEN

Serum. Random or 24-hours urine.

Serum can be stored up to 7 days at 20-25°C or up to 3 weeks at 4-8°C. Samples frozen at -20°C can be stored up to 8 months. Discard contaminated specimens.

24-hours urine preparation: To prevent calcium salt precipitation specimens should be collected in 10 ml of 6M HCl. In case of presence of precipitants they can be solved by lowering pH of the urine to below 2.0. Prior to determination dilute the sample with 0.9% NaCl in the ratio of 1 to 1. Multiply the result by the dilution factor.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

1-Reagent is ready to use.

Applications for analyzers are available on request.

Manual procedure

wavelength	650 nm (630 – 670 nm)
temperature	37°C
cuvette	1 cm

Sample Start method

Pipette into the cuvettes:

	reagent blank (RB)	test (T)	calibrator (C)
1-CALCIUM ARSENAZO	1000 μl	1000 μl	1000 μl
Bring up to the temperature of determination. Then add:			
calibrator	-	-	10 μl
sample	-	10 μl	-

Mix well, incubate for 2 min. Read the absorbance of calibrator samples A(C) and test sample A(T) against reagent blank (RB).

Calculation

$$\text{calcium concentration [mg/dl]} = \frac{A(T)}{A(C)} \times \text{calibrator concentration [mg/dl]}$$

REFERENCE VALUES ⁴

serum	mg/dl	mmol/l
adult	8.6 – 10.3	2.15 – 2.57
random urine	mg/dl	mmol/l
male	0.9 – 37.9	0.225 – 9.47
female	0.5 – 35.7	0.125 – 8.92
24-hours urine	mg/24h	mmol/24h
adult	100 – 300	2.5 – 7.5

It is recommended for each laboratory to establish its own reference ranges for local population.

Calcium concentration 24-hours urine – calculation

$$\text{calcium concentration in 24-hours urine [mg/24h]} = \frac{\text{calcium concentration in sample of 24-hours urine [mg/dl]} \times \text{urine volume of 24-hours urine [dl/24h]}}{\text{urine volume of 24-hours urine [dl/24h]}}$$

QUALITY CONTROL

For internal quality control it is recommended to use, with each batch of samples, the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) and LEVEL 2 (Cat. No 5-162) for determination in urine.

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) and LEVEL 2 (Cat. No 5-175; 5-177) are recommended.

For calibration when using the manual methods the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176) is recommended.

The calibration curve should be prepared with every change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using the automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- LoB (Limit of Blank):** 0.06 mg/dl (0.015 mmol/l)

- LoD (Limit of Detection):** 0.09 mg/dl (0.023 mmol/l)

- LoQ (Limit of Quantitation):** 0.30 mg/dl (0.075 mmol/l)

- Linearity:** up to 23 mg/dl (5.75 mmol/l)

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbate up to 62 mg/l, bilirubin up to 40 mg/dl, triglycerides up to 2000 mg/dl and magnesium up to 20 mg/dl do not interfere with the test.

Precision

Repeatability (run to run) n = 20	Mean [mg/dl]	SD [mg/dl]	CV [%]
level 1	9.67	0.09	0.91
level 2	12.8	0.10	0.79
Reproducibility (day to day) n = 80	Średnia [mg/dl]	SD [mg/dl]	CV [%]
level 1	9.5	0.14	1.4
level 2	12.7	0.13	1.0

Method comparison

A comparison between calcium values determined at **Biolis 24i Premium** (y) and at **BECKMAN COULTER AU680 0** (x) using 30 serum samples gave following results:

$$y = 0.9595 x + 0.3577 \text{ mg/dl;}$$

$$R = 0.991 \quad (R - \text{correlation coefficient})$$

A comparison between calcium values determined at **Biolis 24i Premium** (y) and at **BECKMAN COULTER AU680** (x) using 39 urine samples gave following results:

$$y = 1.0338 x - 0.2229 \text{ mg/dl;}$$

$$R = 1.000 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

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