CREA ENZYMATYC

Nazwa zestawu	(PL) Nr kat.
Liquick Cor-CREA ENZYMATYC 30	2-257
Liquick Cor-CREA ENZYMATYC 60	2-267
Liquick Cor-CREA ENZYMATYC 120	3-333
HC-CREA ENZYMATYC	4-537
OS-CREA ENZYMATYC	9-470
B50-CREA ENZYMATYC	5-514

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LITERATURA

- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed, St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006, p.316.
- Mazzachi BC, Peake M, Erhardt V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

Data wydania: 06. 2021.

CREA ENZYMATIC

Kit name	(EN) Cat. No	sarcosine oxidase	≤ 0.01%
Liquick Cor-CREA ENZYMATIC 30	2-257	ascorbate oxidase	≤ 1%
Liquick Cor-CREA ENZYMATIC 60	2-267	detergents, stabilizers and preservatives	
Liquick Cor-CREA ENZYMATIC 120	3-333	2-CREA ENZYMATIC	
HC-CREA ENZYMATIC	4-537	Good's buffer	≤ 5%
OS-CREA ENZYMATIC	9-470	creatininase	≤ 1%
B50-CREA ENZYMATIC	5-514	peroxidase	≤ 5%
		4-ammoantipyrene (4-AA)	≤ 0.01%
		stabilizers and preservatives	

INTENDED USE





Diagnostic kit for determination of creatinine concentration intended to use for manual assay and in several automatic analyzers. The reagents must be used only for *in vitro* diagnostic, by suitably qualified laboratory personnel, only for the intended purpose, under appropriate laboratory conditions.

INTRODUCTION

Creatinine is a product of creatine nonenzymatic dehydration in skeletal muscle. The amount of creatinine generated and excreted by kidney is proportional to muscle mass and usually is higher in men than women. Daily creatinine generation remains fairly constant, with the exception of crushing injury or degenerative diseases that cause massive damage to muscle. Creatinine blood and urine level depends on glomerular filtration so creatinine clearance is excellent index of renal function.

METHOD PRINCIPLE

Enzymatic, colorimetric method.

The colour intensity measured at 546 nm is proportional to the creatinine concentration.

REAGENTS

Package	Liquick Cor - CREA ENZYMATIC 30	Liquick Cor - CREA ENZYMATIC 60	Liquick Cor - CREA ENZYMATIC 120
1-CREA ENZYMATIC	3 x 30 ml	3 x 60 ml	4 x 112.5 ml
2-CREA ENZYMATIC	1 x 30 ml	1 x 60 ml	2 x 75 ml
	HC - CREA ENZYMATIC	OS - CREA ENZYMATIC	B50 - CREA ENZYMATIC
1-REAGENT	3 x 48 ml	2 x 48.5 ml	2 x 48 ml
2-REAGENT	3 x 15.8 ml	2 x 18.5 ml	2 x 18.2 ml

The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 8 weeks on board the analyser at 2-10°C.

Components and concentrations

1-CREA ENZYMATIC	
Good's buffer	≤ 5%
creatinase	≤ 5%
N-ethyl-N-(3-sulfopropyl)-3-methylaniline (ESPMT)	≤ 5%

CREA ENZYMATIC

51_03_08_028_03

sarcosine oxidase	≤ 0.01%
ascorbate oxidase	≤ 1%
detergents, stabilizers and preservatives	
2-CREA ENZYMATIC	
Good's buffer	≤ 5%
creatininase	≤ 1%
peroxidase	≤ 5%
4-ammoantipyrene (4-AA)	≤ 0.01%
stabilizers and preservatives	

Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.

ADDITIONAL EQUIPMENT

- automatic analyzer or photometer able to read at 546 nm (550 nm);
- thermostat at 37°C;
- general laboratory equipment;

SPECIMEN

Serum and urine.

Urine preparation: before analysis urine sample should be diluted with 0.9% NaCl 2-10 times. Multiply the result by dilution factor. Serum can be stored up to 1 day at 2-8°C. For longer storage samples should be frozen at -20°C.

Urine can be stored up to 1 day at 20-25°C, 4 days at 2-8°C. For longer storage samples should be frozen at -20°C.

Nevertheless it is recommended to perform the assay with freshly collected samples!

PROCEDURE

Applications for analysers are available on request.

NOTE:

1-CREA ENZYMATIC volume, and respectively the sample and 2-CREA ENZYMATIC volumes should be adjusted to measuring photometer capacity.

Manual procedure

wavelength	546 nm (550 nm)
temperature	37°C
cuvette	1 cm

Pipette into the cuvettes:

	standard (S)	test (T)	blank (B)
1-CREA ENZYMATIC	900 µl	900 µl	900 µl

Bring up to the temperature of determination. Then add:

standard	30 µl	-	-
sample	-	30 µl	-
distilled water	-	-	30 µl

Mix well and incubate for 5 minutes at 37 °C. Read the absorbance A₁ of standard samples A(S) and test sample A(T) against reagent blank (B). Then add:

2-CREA ENZYMATIC	300 µl	300 µl	300 µl
------------------	--------	--------	--------

Mix well. Incubate for 5 minutes at 37 °C. Read the absorbance A₂ of standard samples A(S) and test sample A(T) against reagent blank (B). Calculate ΔA (A₂-A₁) for the test and standard.

Calculation

$$\Delta A(T) = (A_2 - A_1)T \times K$$

$$\Delta A(S) = (A_2 - A_1)S \times K$$

$$\text{Creatinine concentration [mg/dl]} = \frac{\Delta A(T)}{\Delta A(S)} \times \frac{\text{standard / calibrator concentration}}{\text{concentration}}$$

$$K = \frac{(\text{sample volume} + R1 \text{ volume})}{(\text{sample volume} + R1 \text{ volume} + R2 \text{ volume})}$$

$$K = 0.756$$

REFERENCE VALUES 2,3

serum / plasma	mg/dl	µmol/l
newborns	0.3 – 1.0	26.5 – 88.4
infants	0.2 – 0.4	17.7 – 35.4
children	0.2 – 0.8	17.7 – 70.7
female	0.5 – 1.0	44.2 – 88.4
male	0.7 – 1.2	61.9 – 106.1
urine (morning)	mg/dl	mmol/l
female	29 – 226	2.56 – 20.0
male	40 – 278	3.54 – 24.6

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY SERUM HN (Cat. No 5-172) and CORMAY SERUM HP (Cat. No 5-173) for determination in serum or CORMAY URINE CONTROL LEVEL 1 (Cat. No 5-161) or LEVEL 2 (Cat. No 5-162) for determination in urine with each batch of samples.

For the calibration of manual assay the CREATININE STANDARD 2 (Cat. No 5-123).

For the calibration of automatic analysers systems the CORMAY MULTICALIBRATOR LEVEL 1 (Cat. No 5-174; 5-176), LEVEL 2 (Cat. No 5-175; 5-177) is recommended.

The calibration curve should be prepared every 4 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analyser Biolis 24i Premium. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 0.08 mg/dl (7.07 µmol/l).

- Linearity:** up to 24 mg/dl (2122 µmol/l).

Specificity / Interferences

Haemoglobin up to 5 g/dl, ascorbic acid up to 62 mg/l, bilirubin up to 20 mg/dl, triglycerides up to 1000 mg/dl and creatine up to 20 mg/dl both at serum and urine do not interfere with the test.

N-acetyl-p-benzoquinone imine (NAPQI), a metabolite of paracetamol (acetaminophen), may cause falsely low results for patients with a toxic level of paracetamol.

Precision

Repeatability (run to run) n = 20	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	1.23	0.02	1.23
level 2	5.63	0.04	0.67
Reproducibility (day to day) n = 80	Mean	SD	CV
	[mg/dl]	[mg/dl]	[%]
level 1	1.17	0.04	3.63
level 2	5.51	0.30	5.42

CREA ENZYMATIC

Method comparison

A comparison between creatinine values determined at **Biolis 24i Premium** (y) and at **Prestige 24i** (x) using 31 samples gave following results:

$$y = 0.9661x + 0.0226 \text{ mg/dl}$$

$$R = 0.9903 \quad (R - \text{correlation coefficient})$$

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed. St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006, p.316.
- Mazzachi BC, Peake M, Erhardt V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

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- Newman DJ, Pnce CP, Renal function and nitrogen metabolites. In: Burtis CA, Ashwood ER, Bruns DE, editors. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 4th ed. St. Louis: W.B Saunders Company; 2006. p. 797-801.
- Alan H.B. Wu. editor. Tietz Clinical Guide to Laboratory Tests, 4th ed. St. Louis: W.B Saunders Company; 2006, p.316.
- Mazzachi BC, Peake M, Erhard V, Reference range and method comparison for enzymatic and Jaffe Creatinine assays in plasma and serum and early morning urine. Clin Lab 2000; 46: 53-5.
- Susumu Osawa, Medical Technology 1982, Vol. 10, No 7, 575-579.
- Minoru Konno, Medical Technology 1984, Vol. 12, No 3, 270-276.
- Schlebusch H, Liappis N, Klein G. Ultrasensitive CRP and Creatinine: Reference intervals from infancy to childhood. Clin Chem Lab Med. 2001; 39 Special supplement pp S1-S448; May 2001. PO-T042.

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