## BILE ACIDS

|                            | (PL)           | HC- BILE ACIDS | OS- BILE ACIDS |
|----------------------------|----------------|----------------|----------------|
| <b>Nazwa zestawu</b>       | <b>Nr kat.</b> |                |                |
| Liquick Cor-BILE ACIDS 30  | 2-338          | 1 x 36,62 ml   | 2 x 14,5 ml    |
| Liquick Cor-BILE ACIDS 60  | 2-339          | 1 x 11,76 ml   | 2 x 5,5 ml     |
| Liquick Cor-BILE ACIDS 120 | 2-340          |                |                |
| HC- BILE ACIDS             | 4-597          |                |                |
| OS- BILE ACIDS             | 9-475          |                |                |

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The reagents when stored at 2-8°C are stable up to expiry date printed on the package. The reagents are stable for 7 weeks on board the analyser at 2-10°C.

### Concentrations in the test

|                     |            |
|---------------------|------------|
| <b>1-BILE ACIDS</b> |            |
| Thio-NAD            | > 0.1 mmol |
| Buffer              |            |
| <b>2-BILE ACIDS</b> |            |
| 3- $\alpha$ -HSD    | > 2 kU/l   |
| NADH                | > 0.1 mmol |
| Buffer              |            |

### Warnings and notes

- Protect from direct sunlight and avoid contamination!
- Avoid contact with skin and mucous membranes.
- Yellow or yellow-brown color of the reagent does not affect the reagents performance.
- Reagents from different lots must not be interchanged.
- Samples from patients treated with ursodeoxycholic acid (UDCA) are not suitable for the determination of total bile acid concentrations.

### SPECIMEN

Serum.  
Total bile acids concentration is increased after meals, therefore samples should be collected under fasting conditions. Serum and plasma samples are stable for a 7 days at 4 °C or for 3 month at -20 °C.  
Nevertheless it is recommended to perform the assay with freshly collected samples!

### ADDITIONAL EQUIPMENT

- automatic analyser or photometer able to read wavelength at 405 nm;
- thermostat at 37°C;
- general laboratory equipment;

### PROCEDURE

Applications for analyzers are available on request.

### Manual procedure

|             |        |
|-------------|--------|
| wavelength  | 405 nm |
| temperature | 37°C   |
| cuvette     | 1 cm   |

Pipette into the cuvettes:

|              | test<br>(T) | calibrator<br>(C) |
|--------------|-------------|-------------------|
| 1-BILE ACIDS | 900 $\mu$ l | 900 $\mu$ l       |
| 2-BILE ACIDS | 300 $\mu$ l | 300 $\mu$ l       |

Bring up to the temperature of determination. Then add:

|            |            |            |
|------------|------------|------------|
| calibrator | -          | 20 $\mu$ l |
| sample     | 20 $\mu$ l | -          |

Mix well and after 2 min. of incubation read the absorbance of calibrator (C) and test (T) against water or air. After next 1, 2, and 3 minutes repeat absorbance reading and calculate the mean absorbance change ( $\Delta A$ ) for calibrator and sample.

### Calculation

$$\text{bile acids concentration} = \frac{\Delta A(T)}{\Delta A(C)} \times \text{calibrator concentration}$$

### REFERENCE VALUES<sup>3</sup>

|       |   |
|-------|---|
| serum | 2.5 – 6.8 $\mu$ mol/l (1.25 – 3.4 $\mu$ g/ml) |
|-------|---|

It is recommended for each laboratory to establish its own reference ranges for local population.

### QUALITY CONTROL

For internal quality control it is recommended to use the CORMAY BILE ACIDS CONTROLS (Cat. No 5-149) with each batch of samples.

For calibration CORMAY BILE ACIDS CALIBRATOR (Cat. No 3-125) is recommended.

The calibration curve should be prepared every 7 weeks, with change of reagent lot number or as required e.g. quality control findings outside the specified range.

### PERFORMANCE CHARACTERISTICS

The following results have been obtained using automatic analysers Biolis 24i Premium and Hitachi 717. Results may vary if a different instrument or a manual procedure is used.

- Sensitivity:** 2.9  $\mu$ mol/l (1.45  $\mu$ g/ml).
- Linearity:** up to 180  $\mu$ mol/l (90  $\mu$ g/ml).

For higher concentration, dilute the sample with 0.9% NaCl and repeat the assay. Multiply the result by dilution factor.

### Specificity / Interferences

Haemoglobin up to 0.5 g/dl, bilirubin up to 50 mg/dl, ascorbic acid up to 50 mg/dl and triglycerides up to 750 mg/dl do not interfere with the test.

### Precision

| Repeatability<br>(run to run) n = 20   | Mean<br>[ $\mu$ mol/l] | SD<br>[ $\mu$ mol/l] | CV<br>[%] |
|--|------------------------|----------------------|-----------|
| level 1                                | 30.72                  | 0.34                 | 1.11      |
| level 2                                | 47.96                  | 0.64                 | 1.34      |
| Reproducibility<br>(day to day) n = 80 | Mean<br>[ $\mu$ mol/l] | SD<br>[ $\mu$ mol/l] | CV<br>[%] |
| level 1                                | 8.12                   | 0.24                 | 2.9       |
| level 2                                | 23.0                   | 0.61                 | 2.6       |

### Method comparison

A comparison between bile acids values determined at Biolis 24i Premium (y) and at OLYMPUS AU400 (x) using 45 samples gave following results:

$$y = 1.0813 x - 0.0198 \mu\text{mol/l};$$

$$R = 0.9997 \quad (R - \text{correlation coefficient})$$


### WASTE MANAGEMENT

Please refer to local legal requirements.

### LITERATURE

- LaRusso, N.F. et al., Dynamics of Enterohepatic Circulation of Bile Acids, New Engl J M, 291, 689-692, (1974).
- Skrede S. et al: Bile acids measured in serum during fasting as a test for liver disease, Clin Chem 24: 1095-1099, 1978
- Wu, Alan H.B. Tietz Clinical Guide to Laboratory Tests. 4th ed. St. Louis, MO: Saunders/Elsevier, 2006. 170-171.
- Dembińska-Kieć A., Naskalski J.W.: Diagnostyka laboratoryjna z elementami biochemii klinicznej, Volumed, 261-262, (1998).

Date of issue: 06. 2021.

The rate of thio-NADH formation can be monitored at 405 nm and is proportional to the bile acids activity.

### REAGENTS

#### Package

|              | Liquick<br>Cor-BILE<br>ACIDS<br>30 | Liquick<br>Cor-BILE<br>ACIDS<br>60 | Liquick<br>Cor- BILE<br>ACIDS<br>120 |
|--------------|------------------------------------|------------------------------------|--------------------------------------|
| 1-BILE ACIDS | 3 x 30 ml                          | 3 x 50 ml                          | 3 x 100 ml                           |
| 2-BILE ACIDS | 1 x 30 ml                          | 1 x 50 ml                          | 1 x 100 ml                           |

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