

Dia-MULTI P**ABNORMAL KONTROLL**

Cat. No.: 90212 12x1ml

Cat. No.: 90201 1x1ml

PRODUCT NAME

Dia-MULTI P

RANGE OF APPLICATION

For quality assurance or control of the coagulationanalyses: antithrombin (AT) For In Vitro Diagnostic use.

PRODUCT DESCRIPTION

Dia-MULTI P (Cat.No.:90212,90201) is a freeze-dried, citrate-anticoagulated humanplasma, with contents of coagulationanalyses within or near the reference values for normal individuals. This control is valid for AT only.

PRECAUTIONS

This product contains substance of human origin and is to be considered as infectious. The product must be treated as plasmasample from a patient. For the purpose of minimizing the danger of infection, every single unit of plasma that has been used during the production has been checked and found free from antibodies against HIV, hepatitis B virus and hepatitis C virus.

RECONSTITUTION

- 1) Before the bottle is opened, knock the bottle against a underlayer to let the powdered material cumulate in the bottom.
- 2) Add 1,00 mL water of good laboratory quality. The water temperature should be 15 to 25 °C.

- 3) After addition of water, reseal the bottle and rotate it ten or more times between your palms in an upright position, backwards and forwards, To dissolve powder residues on the interior surface of the bottle.
- 4) After approx. 15 minutes the control plasma is ready for use. Control visually that all material is completely dissolved, if not use a mixer.

STABILITY AND STORAGE

Unopened reagent is stable until the expiration date shown on the label when stored at 2-8°C.

Stability after reconstitution: 12 hours at 15-25°C in the original vial.

12 hours at 2-8°C in the original vial. (Make sure the control reaches room temperature before use).

MANUFACTURER 

DIAGON Ltd.

Baross u. 52. 1047 Budapest, Hungary

Tel.: +36 1 369 6500,

Fax: +36 1 369 6301

E-mail: diagon@diagon.comWeb: www.diagon.com