Fibronectin

Diagnostic reagent for the quantitative in vitro determination of Fibronectin in human plasma by turbidimetric assay

REF	Content
A00518	1x 5 mL Fibronectin Antibody Reagent 2x 25 mL PEG6 Buffer

Additionally offered:
 A02712 1x 1 mL Fibronectin Calibrator

GENERAL INFORMATION

Method	Immunoturbidimetric
Reaction	Nonlinear, endpoint
Wavelength	340 nm
Assay Temperature	18 – 37 °C
Sample	Plasma
Measuring Range	approx. 0 – 70 mg/dL
Sensitivity	5 mg/dL (Cobas Mira)
Hook Effect	> 1,150 mg/dL (Cobas Mira)

Manual Test Procedure Tests/Kit*

without sample dilution	50
with sample dilution	71

Automated Test Procedure

Instrument dependent – please ask for applications

* calculated on amount of antibody reagent; additional buffer on request

REAGENT COMPOSITION

COMPONENTS FINAL CONCENTRATION

Fibronectin Antibody Reagent	
Turbidimetric grade antibody raised in the goat, monospecific for Fibronectin	variable
Sodium azide	0.095 %

PEG6 Buffer

Phosphate buffered saline	
PEG	6 %
Sodium azide	0.095 %

REAGENT PREPARATION

The reagents are liquid and ready to use.

REAGENT STABILITY AND STORAGE

Conditions: Protect from light. Close immediately after use.

Stability:	at 2 – 8 °C	up to the expiration date
	at 18 – 25 °C	1 month

Do not freeze!

SAMPLE STABILITY AND STORAGE

Stability:	at 2 – 8 °C	48 hours
	at – 20 °C	3 months

Freeze only once!

MANUAL TEST PROCEDURE

Test Procedure without Sample Dilution:

Samples/Controls: ready to use

Calibration curve: Use Fibronectin Calibrator to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent. Use 0.9% saline as zero point.

Pipette into test tubes:	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrs/Samples	15 µL	15 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	100 µL	100 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

Test Procedure with Sample Dilution:

Sample/ Control: dilute 1:10 in saline 0.9%

Calibration curve: Use Fibronectin Calibrator to generate a calibration curve by making 1:10, 1:20, 1:40, 1:80 and 1:160 dilutions with 0.9% saline as diluent. Use 0.9% saline as zero point.

Pipette into test tubes	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrs/Samples	60 µL	60 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	70 µL	70 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

CALCULATION

Calculate and plot $\Delta A = (A2 - A1)$ of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate ΔA optical densities of samples and control(s) and read values in mg/dL on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

REFERENCE RANGE

25 – 40 mg/dL

It is recommended that each laboratory establishes its own normal range.

TEST PRINCIPLE

The assay of Fibronectin is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes. The formation of the complexes is accelerated and enhanced by PEG

DIAGNOSTIC IMPLICATIONS

Fibronectins are large adhesive glycoproteins present in many tissues and in most body fluids. The insoluble form, called cellular fibronectin, is produced by fibroblasts, endothelial cells, macrophages and blood platelets and is a widespread component of connective tissue and basement membranes. The soluble form is referred to as plasma fibronectin. A variety of functions have been ascribed to fibronectin, including roles in opsonization, cell adhesion, cell motility, tissue repair and coagulation. Decreased plasma concentrations of fibronectin occur e.g. in shock, severe infections, hepatic cirrhosis, malnutrition, burns, disseminated intravascular coagulation (DIC), acute pancreatitis as well as after trauma and major surgery. The relative drop in plasma fibronectin concentrations serves as a severity index, with a rise in fibronectin concentration as a prognostically favourable sign. Increased levels of plasma fibronectin have been observed in established gestational proteinuric hypertension.

PERFORMANCE CHARACTERISTICS

SENSITIVITY

5 mg/dL (Cobas Mira)

ACCURACY

A control from Behring was assayed on the Cobas Mira.

Control	Assigned Value (mg/dL)	Measured Value (mg/dL)
Behring	26.5 (22.5 – 30.5)	28.1

PRECISION

Intra-Assay Precision

Two sera (low, medium) were consecutively measured 20 times on the Cobas Mira.

Expected Value	n	Mean	S.D.	C.V.
Low	20	18.3	0.44	2.40
Medium	20	32.4	0.73	2.24

Inter-Assay Precision

After calibration of the SPACE autoanalyser, Fibronectin was measured in two samples at regularly time intervals during 3 weeks. Samples were stored at 4 °C.

Sample	n	Mean	S.D.	C.V.
1	23	45.9	2.90	6.33
2	23	29.9	1.41	4.71

METHOD COMPARISON

A comparison with Nephelometry gave the following results:

$y = 1.1071x + 2.3834$; $r = 0.9251$

INTERFERING SUBSTANCES

No interference up to:

Bilirubin	20 mg/dL	Hemoglobin	1000 mg/dL
Triglycerides	2500 mg/dL	Heparin	50 mg/dL
Sodium Citrate	1000 mg/dL		

QUALITY CONTROL

All commercially available Control sera with Fibronectin values measured by this method may be used.

CALIBRATION

The assay requires the use of Fibronectin Calibrators. We recommend the Dialab Fibronectin Calibrator.

AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

WARNINGS AND PRECATIONS

- The Fibronectin reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

WASTE MANAGEMENT

Please refer to local requirements.

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