# CRP (C-Reactive Protein)

Diagnostic reagent for the quantitative in vitro determination of CRP (C-Reactive Protein) in human serum by turbidimetric assay

| REF                   | Content   |
|-----------------------|---|
| <b>A00519</b>         | 1x 10 mL CRP Antibody Reagent<br>5x 25 mL PEG4 Buffer |
| Additionally offered: |   |
| A00713                | 5x 1 mL CRP Calibrator 5 level series                 |
| A00538                | 1x 1 mL CRP Calibrator High                           |
| A00710                | 1x 5 mL CRP Calibrator High                           |
| A00711                | 1x 1 mL CRP Calibrator Super High                     |
| A00712                | 1x 5 mL CRP Calibrator Super High                     |
| A00822                | 1x 1 mL CRP Calibrator Low                            |
| A00709                | 1x 5 mL CRP Calibrator Low                            |
| A00520                | 1x 1 mL CRP Control Low                               |
| A00803                | 1x 5 mL CRP Control Low                               |
| A00804                | 1x 1 mL CRP Control High                              |
| A00805                | 1x 5 mL CRP Control High                              |
| A00590                | 1x 1 mL Protein Control                               |
| A00800                | 1x 5 mL Protein Control                               |
| A08591                | 1x 1 mL Protein Control Low                           |
| A08823                | 1x 5 mL Protein Control Low                           |
| A04823                | 1x 1 mL Triple Control (ASO, CRP, RF)                 |
| A04824                | 1x 5 mL Triple Control (ASO, CRP, RF)                 |

## GENERAL INFORMATION

|                          |  |
|--------------------------|--|
| <b>Method</b>            | Immunoturbidimetric                              |
| <b>Reaction</b>          | Nonlinear, endpoint                              |
| <b>Wavelength</b>        | 340 nm   |
| <b>Assay Temperature</b> | 18 – 37 °C                                       |
| <b>Sample</b>            | Serum  |
| <b>Measuring Range</b>   | approx. 0 – 22 mg/dL                             |
| <b>Sensitivity</b>       | 0.3 mg/dL (Cobas Mira)                           |
| <b>Hook Effect</b>       | without sample dilution: > 43 mg/dL (Cobas Mira) |

**Manual Test Procedure Tests/Kit\***  
without sample dilution 133

**Automated Test Procedure**  
Instrument dependent – please ask for applications

\* calculated on amount of antibody reagent; additional buffer on request

## REAGENT COMPOSITION

| COMPONENTS                          | FINAL CONCENTRATION |
|-------------------------------------|---------------------|
| <b>CRP Antibody Reagent</b>         |                     |
| Polyclonal goat anti-h CRP antibody | variable            |
| Sodium azide                        | 0.095 %             |
| <b>PEG4 Buffer</b>                  |                     |
| Phosphate buffered saline, pH 7.4   |                     |
| PEG                                 | 4 %                 |
| Sodium azide                        | 0.095 %             |

## REAGENT PREPARATION

The reagents are liquid and ready to use.

## REAGENT STABILITY AND STORAGE

|                    |  |
|--------------------|--|
| <b>Conditions:</b> | Protect from light. Close immediately after use.               |
| <b>Stability:</b>  | at 2 – 8 °C up to the expiration date<br>at 18 – 25 °C 1 month |

Do not freeze!  
Incompetent handling will release DIALAB from any responsibility.

## SAMPLE STABILITY AND STORAGE

|                   |   |
|-------------------|---|
| <b>Stability:</b> | at 2 – 8 °C 48 hours<br>at – 20 °C 3 months |
|-------------------|---|

Freeze only once!

## MANUAL TEST PROCEDURE

### Test Procedure without Sample Dilution:

Samples/Controls: ready to use  
Calibration curve: Use CRP Calibrator High to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent or use the 5 level calibrator series. Use 0.9% saline as zero point.

|  |             |                  |
|--|-------------|------------------|
| Pipette into test tubes:   | Calibrators | Samples/Controls |
| Buffer   | 900 µL      | 900 µL           |
| Cal./Ctrls/Samples   | 60 µL       | 60 µL            |
| Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:  |             |                  |
| Antibody Reagent   | 75 µL       | 75 µL            |
| Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$ |             |                  |

## CALCULATION

Calculate and plot  $\Delta A = (A2 - A1)$  of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate  $\Delta A$  optical densities of samples and control(s) and read values in mg/dL on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

## REFERENCE RANGE

0 – 1 mg/dL  
It is recommended that each laboratory establishes its own normal range.

## TEST PRINCIPLE

The assay of CRP is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes. The formation of the complexes is accelerated and enhanced by PEG.

## DIAGNOSTIC IMPLICATIONS

C-Reactive Protein (CRP) is an acute marker of inflammatory processes. In case of an acute inflammation the concentration of CRP increases and decreases more rapidly than the red cell sedimentation rate. The increase of CRP occurs in a non-specific way in different kinds of tissular aggression, as for example in infectious states, rheumatoid arthritis, myocardial infarction, malignant tumor, etc.

## PERFORMANCE CHARACTERISTICS

**SENSITIVITY**  
0.3 mg/dL (Cobas Mira)

## ACCURACY

| Analyzer    | Control    | Assigned Value (mg/dL) | Measured (mg/dL) |
|-------------|------------|------------------------|------------------|
| Cobas Mira  | DIALAB     | 1.80 (1.53 - 2.07)     | 1.82             |
| Cobas Mira  | DIALAB     | 3.50 (2.98 - 4.03)     | 3.34             |
| Cobas Mira  | Seronorm   | 3.10 (2.63 - 3.57)     | 2.78             |
| Cobas Mira  | CRM/470    | 3.92 (3.83 - 4.51)     | 3.90             |
| Express 550 | Behring    | 8.3 (7.05 - 9.55)      | 7.85             |
| Express 550 | Biokit Lot | 4.0 (3.4 - 4.6)        | 3.94             |

## PRECISION

### Intra-Assay Precision

Patient sera with positive CRP values were measured.

| Analyzer   | n  | Mean  | S.D. | C.V. |
|------------|----|-------|------|------|
| Cobas Mira | 20 | 1.23  | 0.5  | 4.06 |
|            | 20 | 3.49  | 0.9  | 2.57 |
|            | 20 | 10.15 | 3.5  | 3.44 |

### Inter-Assay Precision

DIALAB CRP Control serum was measured daily on the MEGA during one month after calibration. Furthermore a serum with a high CRP level was divided in aliquots and stored at -20 °C. After calibration it was measured on 6 occasional days on the Hitachi 717.

| Analyzer    | n  | Mean | S.D. | C.V. |
|-------------|----|------|------|------|
| MEGA        | 26 | 1.53 | 0.07 | 4.4  |
| Hitachi 717 | 6  | 8.4  | 0.3  | 3.5  |

## METHOD COMPARISON

A comparison with Nephelometry gave the following results:  
 $y = 0.9981x + 0.014254$ ;  $r = 0.9917$

## INTERFERING SUBSTANCES

No interference up to:

|               |            |                |            |
|---------------|------------|----------------|------------|
| Triglycerides | 2500 mg/dL | Hemoglobin     | 1000 mg/dL |
| Bilirubin     | 20 mg/dL   | Sodium Citrate | 1000 mg/dL |
| Heparin       | 50 mg/dL   |                |            |

## QUALITY CONTROL

All Control sera with CRP values measured by this method may be used. We recommend the Dialab CRP Control Low and CRP Control High, the Protein Control and the Protein Control Low and the Triple Control.

## CALIBRATION

The assay requires the use of CRP serum Calibrators. We recommend the Dialab CRP Calibrator 5 Level Series, the CRP Calibrator High, the CRP Calibrator Super High or the CRP Calibrator Low.

The use of other commercially available CRP calibrators is not recommended. It is sufficient to perform calibration once a month.

## AUTOMATION

Applications for automated systems (with and without sample dilution) are available upon request.

## WARNINGS AND PRECAUTIONS

- The CRP reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA
- Avoid eyes and skin contact. If contact, flush with a large amount of water. If irritation persists, consult a physician

## WASTE MANAGEMENT

Please refer to local requirements.

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