Microalbumin

Diagnostic reagent for the quantitative in vitro determination of Microalbumin in human urine by turbidimetric assay

REF

Content

A00528 1x 10 mL Microalbumin Antibody Reagent
 5x 25 mL Microalbumin Buffer

Additionally offered:

A00955 5x 1 mL Microalbumin Calibrator 5 Level Series
 A00707 1x 1 mL Microalbumin Calibrator
 A00708 1x 5 mL Microalbumin Calibrator
 A00808 1x 1 mL Microalbumin Control
 A00809 1x 5 mL Microalbumin Control

GENERAL INFORMATION

Method Immunoturbidimetric
Reaction Nonlinear, endpoint
Wavelength 340 nm
Assay Temperature 18 – 37 °C
Sample Urine
Measuring Range approx. 0 – 400 mg/L
Sensitivity 12.5 mg/L (Cobas Mira)
 High sensitive application: 3 mg/L (Hitachi 911)
 > 6,000 mg/L (Hitachi 911)
Hook Effect
Manual Test Procedure Tests/Kit*
 without sample dilution 100
 without sample dilution 66 (high sensitivity application)
Automated Test Procedure
 Instrument dependent – please ask for applications
 * calculated on amount of antibody reagent; additional buffer on request

REAGENT COMPOSITION

COMPONENTS	FINAL CONCENTRATION
Microalbumin Antibody Reagent	
Polyclonal goat anti-human Albumin antibody	variable
Sodium azide	0.095 %
Microalbumin Buffer	
Saline	9 g/L
Enhancer	
Sodium azide	0.095 %

REAGENT PREPARATION

The reagents are liquid and ready to use.

REAGENT STABILITY AND STORAGE

Conditions: Protect from light. Close immediately after use.
Stability: at 2 – 8 °C up to the expiration date
 at 18 – 25 °C 1 month
 Do not freeze!

SAMPLE STABILITY AND STORAGE

Collect urine during 24 hours or as a random midstream sample.
Stability: at 2 – 8 °C 48 hours
 at – 20 °C 3 months
 Centrifuge the urine before assay. Freeze only once!

MANUAL TEST PROCEDURE

Test Procedure without Sample Dilution:

Samples/Controls: ready to use
 Calibration curve: Use Microalbumin Calibrator High to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent or use the 5 level calibrator series. Use 0.9% saline as zero point.

Pipette into test tubes:	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	20 µL	20 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	100 µL	100 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

High Sensitive Test Procedure without Sample Dilution:

Samples/Controls: ready to use
 Calibration curve: Use Microalbumin Calibrator High to generate a calibration curve by making 1:2 serial dilutions of the calibrator with 0.9% saline as diluent or use the 5 level calibrator series. Use 0.9% saline as zero point.

Pipette into test tubes	Calibrators	Samples/Controls
Buffer	900 µL	900 µL
Cal./Ctrls/Samples	60 µL	60 µL
Mix. Read A1 of calibrators and samples/controls at 340 nm. Then add:		
Antibody Reagent	150 µL	150 µL
Mix. Incubate 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 340 nm. Calculate: $\Delta A = (A2 - A1)$		

CALCULATION

Calculate and plot $\Delta A = (A2 - A1)$ of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate ΔA optical densities of samples and control(s) and read values in mg/L on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

REFERENCE RANGE

0 – 25 mg/L (IFCC)
 It is recommended that each laboratory establishes its own normal range.

TEST PRINCIPLE

The assay of Microalbumin is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes.

DIAGNOSTIC IMPLICATIONS

Diabetic nephropathy, which is accompanied by irreversible kidney damage and persistent proteinuria, is a major cause of death in persons with insulin-dependent diabetes mellitus. An early sign of diabetic nephropathy are small Albumin secretions in urine, i.e. Microalbuminuria. Therefore, detection of kidney (glomerular) damage that is minimal and reversible is important.

PERFORMANCE CHARACTERISTICS

SENSITIVITY

12.5 mg/L (Cobas Mira)
 High sensitive application: 3 mg/L (Hitachi 911)

ACCURACY

A Control from DIALAB was assayed in duplicate on the Cobas Mira and Hitachi 911 (high sensitive assay).

Control	Assigned Value (mg/L)	Measured Value (mg/L)
DIALAB	200 (170 – 230)	201

High sensitive assay:

Control	Assigned Value (mg/L)	Measured Value (mg/L)
DIALAB	196 (166 – 225)	196.7

PRECISION

Intra-Assay Precision

3 control solutions (low-medium-high) were consecutively measured 20x on the Cobas Mira and Hitachi 911 (high sensitive assay) and the variation coefficient was calculated.

Expected Value	n	Mean	S.D.	C.V
Low	20	16.94	1.43	5.90
Medium	20	55.56	2.22	3.99
High	20	200.94	5.36	2.67

High sensitive assay:

Expected Value	n	Mean	S.D.	C.V
Low	20	32.29	0.29	0.91
Medium	20	92.46	0.82	0.88
High	20	196.01	2.98	1.52

Inter-Assay Precision

A Control was measured daily on the Cobas Mira and Hitachi 911 (high sensitive assay) after calibration. The variation coefficient was calculated.

Expected Value	n	Mean	S.D.	C.V
DIALAB	32	204.9	6.97	3.35

High sensitive assay:

Expected Value	n	Mean	S.D.	C.V
Low	20	22.29	0.65	2.93
Medium	20	90.69	0.60	0.66
High	20	192.27	1.03	0.53

METHOD COMPARISON

A comparison with Nephelometry gave the following results:
 $y = 1.1702x + 1.4811$; $r = 0.9879$

High sensitive assay: $y = 1.0096x - 0.2344$; $r = 0.9978$

INTERFERING SUBSTANCES

No interference up to:

Triglycerides	2500 mg/dL	Hemoglobin	500 mg/dL
Bilirubin	15 mg/dL	Sodium Citrate	1000 mg/dL
Heparin	100 mg/dL	Turbidity	2.5 %

QUALITY CONTROL

All commercially available Control sera with Microalbumin values measured by this method may be used. We recommend the Dialab Microalbumin control.

CALIBRATION

The assay requires the use of Microalbumin Calibrators. We recommend the Dialab Microalbumin Calibrator 5 Level Series or the Microalbumin Calibrator.

AUTOMATION

Applications for automated systems are available upon request.

WARNINGS AND PRECATIONS

- The Microalbumin reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

WASTE MANAGEMENT

Please refer to local requirements.

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