## Protein Calibrator Low

Liquid human based calibrator for use as a calibrator in immunoturbidimetric assays

REF	Cont.	
A00701	1x 1 mL	Protein Calibrator Low
A00702	1x 5 mL	Protein Calibrator Low

### COMPOSITION

The Protein Calibrator Low is a defibrinated human plasma, liquid stabilized and filtered through a 0.2 µ pore-size filter.  
 Sodium azide 0.095 %

### PREPARATION

The Calibrator is liquid, ready to use.

### STABILITY AND STORAGE

Storage: at 2 – 8 °C  
 Stability: up to the expiration date  
 Stability after opening: 6 weeks at 2 – 8 °C  
 DO NOT FREEZE!

### DIRECTIONS FOR USE

- Use this Calibrator in the same manner as specified with the reagent and/or instrument being used.
- Allow the Calibrator to equilibrate to room temperature (22 – 28 °C) prior to use.
- The Calibrator should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use

### WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

### WASTE MANAGEMENT


Please refer to local legal requirements.

### LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: MPL007	EXP: 2017/02
Parameter	Assigned Value *
Albumin	5,963.0 mg/dL
α-1 Antitrypsin	201.0 mg/dL
α-2 Macroglobulin	380.0 mg/dL
α-1 Acid Glycoprotein	113.0 mg/dL
Ceruloplasmin	51.0 mg/dL
Haptoglobin	176.0 mg/dL
Transferrin	401.0 mg/dL
IgA	286.0 mg/dL
IgG	1,448.0 mg/dL
IgM	180.0 mg/dL
Complement C3	218.0 mg/dL
Complement C4	37.00 mg/dL
Antithrombin III	41.00 mg/dL
C1 Esterase Inhibitor **	41.00 mg/dL
Kappa Light Chain **	357.0 mg/dL
Lambda Light Chain **	194.0 mg/dL
Prealbumin	38.0 mg/dL

\* Values traceable to ERM-DA470k from IFCC.  
 \*\* Values based on Siemens standard material.

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