# CRP Calibrator 5 Level Series

*Calibrator set for preparation of reference curves for quantitative immunoturbidimetric determination of CRP*

**REF**

**Content**

**A00713** - 5x 1 mL CRP Calibrator 5 Level Series

**For professional in vitro diagnostic use only.**

**GENERAL INFORMATION**

**Method** Immunoturbidimetry  
**Shelf life** 36 months from date of production  
**Storage** at 2 – 8 °C  
**Other Info** Liquid, ready to use

**INTENDED USE**

DIALAB CRP Calibrator 5 Level Series is used for the preparation of reference curves for quantitative immunoturbidimetric determination of CRP.

**DIAGNOSTIC SIGNIFICANCE**

Please refer to the corresponding section of the used CRP reagent.

**TEST PRINCIPLE**

Please refer to instructions for use of used CRP reagent.

**COMPOSITION**

Dialab CRP Calibrators of the CRP Calibrator 5 Level Series are a dilution of defibrinated human plasma and recombinant human C-reactive protein with phosphate buffered saline. The product is stabilised with 0.095 % sodium azide as preservative.

**MATERIAL REQUIRED BUT NOT PROVIDED**

- DIALAB CRP reagent
- Clinical chemistry analyser or photometer
- DIALAB CRP Controls

**REAGENT PREPARATION**

Calibrator 5 Level Series is ready to use.

**STORAGE AND STABILITY**

**Storage:** at 2 - 8°C  
**Stability:** up to expiration date  
**Stability after opening:** 6 weeks at 2 - 8°C if stored tightly closed

The calibrators should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature

**DO NOT FREEZE!**

**WARNINGS AND PRECAUTIONS**

- For in Vitro Diagnostic use.
- Keep vials tightly capped at all time when not in use to avoid any microbial contamination and evaporation.
- Keep vials refrigerated when not in use.
- Each individual donation intended for use in manufacture of Calibrator 5 Level Series was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests and only donations with negative findings were used for its manufacture. Nevertheless, every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.
- Reagents containing sodium azide must be handled with caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

**TEST PROCEDURE**

1. Allow the calibrators to equilibrate to room temperature (22 – 28 °C) prior to use.
2. Use and dilute the calibrators as described in the documentation of the used CRP reagent and according to the specific requirements (instrument manual) of your instrument.

**INTERPRETATION OF RESULTS**

The results of the different dilutions should be within their expectation. Use controls to verify your calibration.

**QUALITY CONTROL AND CALIBRATION**

For calibration please refer to the instructions for use of the used CRP reagent, calibrator and instrument. It is recommended to use quality controls after each calibration, at new measuring sessions, after each change of reagents and every time where a potential increased risk of malfunction can be assumed or when results indicate a potential malfunction. Each laboratory should establish its own quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

**PERFORMANCE CHARACTERISTICS**

Please refer to the section "PERFORMANCE CHARACTERISTICS" of the used CRP reagent

**TRACEABILITY**

The values of Calibrator 5 Level Series have been assigned to ERM-DA474/IFCC.

**EXPECTED VALUES**

For the biological reference range please refer to the instruction for use of the used CRP reagent.

**LIMITATIONS**

For information about biological limitations please refer to the description for use of the used CRP reagent.

**WASTE MANAGEMENT**

Please refer to local legal requirements.

**LITERATURE**

Please refer to the corresponding section of the used CRP reagent

**USED SYMBOLS**

**Symbol**

**Description**


Dispose product, residuals and packaging appropriately

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REF	Inhalt
A00713	- 5x 1 mL CRP Calibrator 5 Level Series

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