α-1 Microglobulin Calibrator

Liquid human based calibrator for use as a calibrator in immunoturbidimetric assays for α-1 Microglobulin

REF	Cont.	
A00737	1x 1 mL	α-1 Microglobulin Calibrator
A00738	1x 5 mL	α-1 Microglobulin Calibrator

COMPOSITION

Dialab α-1 Microglobulin Calibrator is based in human urine.
 Sodium azide 0.095 %

PREPARATION

The calibrator is liquid and ready to use.

STABILITY AND STORAGE

Storage: at 2 – 8 °C
 Stability: up to the expiration date
 Stability after opening: 6 weeks at 2 – 8 °C
 DO NOT FREEZE!

DIRECTIONS FOR USE

1. Use this Calibrator in the same manner as specified with the reagent and/or instrument being used.
2. Allow the Calibrator to equilibrate to room temperature (22 – 28 °C) prior to use.
3. The calibrator should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use.

WARNINGS AND PRECATIONS

1. For in Vitro Diagnostic use.
2. Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
3. Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

WASTE MANAGEMENT

Please refer to local legal requirements.

LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: S0114		EXP: 2016/05	
Parameter	Assigned Value		
α-1 Microglobulin	49.80 mg/L		
Values are traceable to Dade Behring.			

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