APO A1/B Control

(en) English

REF

Content

A00806 - 1x 1 mL lyophilized human based control

For professional in vitro diagnostic use only.

INTENDED USE

APO A1/B Control is a human based control for use as a control in immunoturbidimetric assays for APO A1 and APO B.

DIAGNOSTIC SIGNIFICANCE

Please refer to the corresponding section of the used APO A1 / APO B reagent

TEST PRINCIPLE

Please refer to the corresponding section of the used APO A1 / APO B reagent

REAGENT COMPOSITION

Lyophilized preparation of stabilized, buffered defibrinated human plasma.
 Sodium azide 0.095 %

MATERIAL REQUIRED BUT NOT PROVIDED

- APO A1 or APO B reagent, e.g.: Dialab APO A1 5+1 or APO B 5+1
- APO A1 or APO B calibrator, e.g.: APO A1/B Calibrator High
- Clinical chemistry analyser or photometer

REAGENT PREPARATION

The control is supplied in dry form. Dissolve the vial contents in exactly 1 mL distilled water and let stand at 15 - 25°C for 30 minutes. Invert gently to mix. Avoid foam formation and vigorous shaking.

STORAGE AND STABILITY

Storage: at 2 - 8 °C
 Stability: up to the expiration date

After reconstitution, the APO A1/B Control is stable for 7 days at 2 - 8 °C.

The APO A1/B Control should not be allowed to stand for repeated long periods of time at room temperature. Keep vials tightly capped all time when not in use to avoid evaporation and microbial contamination. Keep vials refrigerated when not in use DO NOT FREEZE!

WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Each individual donation intended for use in manufacture of A1/B Control was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2, Hepatitis B, Hepatitis C virus or other pathogens, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
- Reagents containing sodium azide must be handled with caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.
- Please refer to the safety data sheets (SDS) and take the necessary precautions for the use of laboratory reagents.
- In case of product malfunction or altered appearance that could affect the performance, contact the manufacturer.
- In the event of an incident related to the device, report it to the manufacturer and your competent authority as required.

SPECIMEN COLLECTION AND STORAGE

Please refer to the corresponding section of the used APO A1 / APO B reagent.

TEST PROCEDURE

- Use the APO A1/B Control in the same manner as specified with the reagent

and/or instrument being used.

- Allow the APO A1/B Control to equilibrate to room temperature (22 - 28 °C) prior to use.

INTERPRETATION OF RESULTS

The assay values for APO A1/B Control have been determined by calibration with IFCC APO A1 Reference (BCR-393) / IFCC APO B Reference (SP3-07) material in accordance with established protocols. The assay results with the controls should be within the specific range as shown in the table below and are specific for this lot number of each control only. Assay values may vary slightly with different reagents and methodologies used.

QUALITY CONTROL AND CALIBRATION

Please refer to the instructions for use of the used APO A1/B reagent, calibrator and instrument. It is recommended to use quality controls after each calibration, at new measuring sessions, after each change of reagents and every time where a potential increased risk of malfunction can be assumed or when results indicate a potential malfunction and to repeat the calibration where required. Each laboratory should establish its own quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

PERFORMANCE CHARACTERISTICS

Please refer to the corresponding section of the used APO A1 / APO B reagent

TRACEABILITY

Values of this control for APO A1 were titred against IFCC APO A1 Reference (BCR-393), values of APO B were titred against IFCC APO B Reference (SP3-07).

EXPECTED VALUES

For the biological reference range of APO A1 and APO B please refer to the instruction for use of the used APO A1/B reagent.

LIMITATIONS

Do not use the product after the expiry date.

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE


Please refer to the corresponding section of the used APO A1/B reagent

USED SYMBOLS

Symbol	Description
	Content
	Dispose product, residual and packaging appropriately

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