# CRP Calibrator Low

Calibrator for preparation of reference curves for quantitative immunoturbidimetric determination of CRP

REF	Content
A00822	- 1x 1 mL CRP Calibrator Low
A00709	- 1x 5 mL CRP Calibrator Low

For professional in vitro diagnostic use only.

## GENERAL INFORMATION

Method	Immunoturbidimetry
Shelf life	36 months from date of production
Storage	at 2 – 8 °C
Other Info	Liquid, ready to use

## INTENDED USE

DIALAB CRP Calibrator Low is used for the preparation of reference curves for quantitative immunoturbidimetric determination of CRP.

## DIAGNOSTIC SIGNIFICANCE

Please refer to the corresponding section of the used CRP reagent.

## TEST PRINCIPLE

Please refer to instructions for use of used CRP reagent.

## COMPOSITION

Dialab CRP Calibrator Low is a dilution of defibrinated human plasma and recombinant human C-reactive protein with phosphate buffered saline. It is stabilised with 0.095 % sodium azide as preservative

## MATERIAL REQUIRED BUT NOT PROVIDED

- DIALAB CRP reagent
- Clinical chemistry analyser or photometer
- DIALAB CRP Controls

## REAGENT PREPARATION

CRP Calibrator Low is ready to use.

## STORAGE AND STABILITY

Storage:	at 2 - 8°C
Stability:	up to expiration date
Stability after opening:	6 weeks at 2 - 8°C if stored tightly closed

The calibrator should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature

DO NOT FREEZE!

## WARNINGS AND PRECAUTIONS

- For in Vitro Diagnostic use.
- Keep vials tightly capped at all time when not in use to avoid any microbial contamination and evaporation.
- Keep vials refrigerated when not in use.
- Each individual donation intended for use in manufacture of CRP Calibrator Low was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C virus (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests and only donations with negative findings were used for its manufacture. Nevertheless, every product obtained from human body fluids should be handled with appropriate care in accordance with recommended procedures for biohazardous materials since absence of infectious agents can never be proven.

- Reagents containing sodium azide must be handled with caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

## TEST PROCEDURE

1. Allow the calibrator to equilibrate to room temperature (22 – 28 °C) prior to use.
2. Use and dilute this calibrator as described in the documentation of the used CRP reagent and according to the specific requirements (instrument manual) of your instrument.

## INTERPRETATION OF RESULTS

The results of the different dilutions should be within their expectation. Use controls to verify your calibration.

## QUALITY CONTROL AND CALIBRATION

For calibration please refer to the instructions for use of the used CRP reagent, calibrator and instrument. It is recommended to use quality controls after each calibration, at new measuring sessions, after each change of reagents and every time where a potential increased risk of malfunction can be assumed or when results indicate a potential malfunction. Each laboratory should establish its own quality control scheme and corrective actions if controls do not meet the acceptable tolerances.

## PERFORMANCE CHARACTERISTICS

Please refer to the section "PERFORMANCE CHARACTERISTICS" of the used CRP reagent

## TRACEABILITY

The value of CRP Calibrator Low has been assigned to BCR-ERM-DA474/IFCC.

## EXPECTED VALUES

For the biological reference range please refer to the instruction for use of the used CRP reagent.

## LIMITATIONS

For information about biological limitations please refer to the description for use of the used CRP reagent.

## WASTE MANAGEMENT

Please refer to local legal requirements.

## LITERATURE

Please refer to the corresponding section of the used CRP reagent

## USED SYMBOLS

Symbol	Description
	Dispose product, residuals and packaging appropriately

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