Fibrinogen Calibrator

Lyophilised human based calibrator for use as a calibrator in immunoturbidimetric assays for Fibrinogen

REF

A01520

Cont.

1x 0.5 mL Fibrinogen Calibrator (lyo.)

COMPOSITION

Dialab Fibrinogen Calibrator is a lyophilised preparation of stabilised, buffered human plasma containing less than 5% bovine plasma.

Sodium azide 0.095 %

PREPARATION

1. Reconstitute the Fibrinogen Calibrator with 0.5 mL of distilled or deionized water.
2. Close vial and allow to stand until fully reconstituted (a minimum of 5 minutes).
3. Gently mix by rotation or inversion until dissolution is complete. DO NOT SHAKE.
4. Allow to equilibrate for a total of 30 minutes before use.

STABILITY AND STORAGE

Storage: at 2 – 8 °C
 Stability: up to the expiration date
 Stability after reconstitution: 48 hours at 2 – 8 °C
 DO NOT FREEZE!

DIRECTIONS FOR USE

1. Use this Calibrator in the same manner as specified with the reagent and/or instrument being used.
2. Allow the Calibrator to equilibrate to room temperature (22 – 28 °C) prior to use.
3. The Calibrator should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use.

WARNINGS AND PRECATIONS

1. For in Vitro Diagnostic use.
2. Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
3. Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

WASTE MANAGEMENT

Please refer to local legal requirements.

LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.

LOT: 07022	EXP: 2016/07
Parameter	Assigned Value
Fibrinogen	446.0 mg/dL
Values are traceable to Dade Behring standardisation.	

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