Ferritin

Diagnostic reagent for the quantitative in vitro determination of Ferritin in human serum by turbidimetric assay

REF	Content
A13550	2x 10 mL Ferritin Latex Reagent 2x 30 mL Ferritin Buffer

Additionally offered:

A06560	5x 1 mL Ferritin Calibrator 5 level series
A00610	1x 1 mL Ferritin Control High
A00821	1x 5 mL Ferritin Control High
A00570	1x 1 mL Ferritin Control Low
A00820	1x 5 mL Ferritin Control Low
A00590	1x 1 mL Protein Control
A00800	1x 5 mL Protein Control
A08591	1x 1 mL Protein Control Low
A08823	1x 5 mL Protein Control Low

GENERAL INFORMATION

Method	Immunoturbidimetric
Reaction	Nonlinear, endpoint
Wavelength	600 nm
Assay Temperature	18 – 37 °C
Sample	Serum
Measuring Range	approx. 0 – 500 ng/mL
Sensitivity	5 ng/mL (Hitachi 911)
Hook Effect	no risk

Manual Test Procedure Tests/Kit*
without sample dilution 61

Automated Test Procedure

Instrument dependent – please ask for applications

* calculated on amount of antibody reagent; additional buffer on request

REAGENT COMPOSITION

COMPONENTS	FINAL CONCENTRATION
Ferritin Latex Reagent	
Suspension of suspended latex microparticles sensitized with rabbit IgG anti-human Ferritin	variable
Ferritin Buffer	
Hepes Buffer	
Sodium azide	0.095 %

REAGENT PREPARATION

The reagents are liquid and ready to use.

REAGENT STABILITY AND STORAGE

Conditions:	Protect from light. Close immediately after use.
Stability:	at 2 – 8 °C up to the expiration date at 18 – 25 °C 1 month

Do not freeze!

SAMPLE STABILITY AND STORAGE

Stability:	at 2 – 8 °C 48 hours at – 20 °C 3 months
-------------------	---

Freeze only once!

MANUAL TEST PROCEDURE

Test Procedure without Sample Dilution:

Samples/Controls: ready to use

Calibration curve: Use Ferritin Calibrator 5 level series to generate a calibration curve. Use 0.9% saline as zero point.

Pipette into test tubes:	Calibrators	Samples/Controls
Buffer	975 µL	975 µL
Cal./Ctrls/Samples	130 µL	130 µL
Mix and incubate for 1 min. Then add:		
Latex Reagent	325 µL	325 µL
Mix and incubate for 30 seconds. Read A1 of calibrators and samples/controls at 600 nm.		
Incubate for 5 minutes at assay temperature. Read A2 of calibrators and samples/controls at 600 nm. Calculate: $\Delta A = (A2 - A1)$		

CALCULATION

Calculate and plot $\Delta A = (A2 - A1)$ of the calibrators versus assigned concentration values on a linear-linear graph paper. Calculate ΔA optical densities of samples and control(s) and read values in ng/mL on the reference curve. Samples yielding absorbances above highest calibrator should be retested after further dilution.

REFERENCE RANGE

Men:	20 – 300 ng/mL
Women:	15 – 200 ng/mL

It is recommended that each laboratory establishes its own normal range.

TEST PRINCIPLE

The assay of Ferritin is based on turbidimetric measurement. Turbidity is caused by the formation of antigen-antibody insoluble immuno complexes.

DIAGNOSTIC IMPLICATIONS

The plasma Ferritin concentration declines very early in the development of iron deficiency. On the other hand, a large number of chronic diseases result in increased serum Ferritin concentrations. These diseases include chronic infections, chronic inflammatory disorders such as rheumatoid arthritis or renal disease, Gaucher's disease, and numerous types of malignancies, especially lymphomas, leukaemia's, breast cancer and neuroblastoma. Increase in plasma Ferritin concentration also occurs in viral hepatitis or following toxic liver injury as a result of release of Ferritin from damaged liver cells. Plasma Ferritin concentration is also increased with increases of iron stores, as seen in patients with hemosiderosis or hemochromatosis. Besides the use of Ferritin as an iron metabolism parameter, Ferritin as also gained importance as a tumor marker for therapeutic drug monitoring and follow-up.

PERFORMANCE CHARACTERISTICS

The performance characteristics for the Ferritin reagents were measured on a clinical chemistry analyzer according to the procedure R1: 150 µL and R2: 50 µL.

SENSITIVITY

5 ng/mL

ACCURACY

Controls were assayed in duplicate for Ferritin to verify proper assay recovery.

Control	Assigned Value (ng/mL)	Measured Value (ng/mL)
DIALAB Control Low	85.2 (72.4-98.0)	92.1
DIALAB Control High	224 (190-258)	221.1
SIEMENS	87.7 (70.2-105.2)	82.5

PRECISION

	Low % CV	Medium % CV	High % CV
Intra-Run	0.76	0.63	0.61
Inter-Run	4.14	3.80	4.11

METHOD COMPARISON

A comparison with external ferritin reagents the following results:
 $y = 0.9261x - 15.596$; $r = 0.9868$

INTERFERING SUBSTANCES

No interference up to:

Bilirubin	30 mg/dL	Hemoglobin	250 mg/dL
Triglycerides	2500 mg/dL	Sodium Citrate	1000 mg/dL
EDTA	5 mg/mL	Heparin	50 mg/dL

QUALITY CONTROL

All commercially available Control sera with Ferritin values measured by this method may be used. We recommend the Dialab Ferritin Control High, the Ferritin Control Low, the Protein Control and the Protein Control Low.

CALIBRATION

The assay requires the use of Ferritin Calibrators. We recommend the Dialab Ferritin Calibrator 5 level series..

AUTOMATION

Applications for automated systems are available upon request.

WARNINGS AND PRECATIONS

- The Ferritin reagents are intended for in vitro diagnostic use only.
- Sodium azide has been reported to form lead or copper azide in laboratory plumbing which may explode on percussion.
- Each donor unit used in the preparation of the standards and controls was found to be negative for the presence of HIV antibodies, as well as for Hepatitis B surface antigen, using a method approved by the FDA

WASTE MANAGEMENT

Please refer to local requirements.

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