CRP Uhs Calibrator Set for the high sensitivity range

For the use as calibrators in the determination of CRP in the high sensitivity range by immunoturbidimetric assay

REF

D04571

Cont.

5x 1 mL

CRP Uhs Calibrator Set
(for the high sensitivity range)

COMPOSITION

The CRP Uhs Calibrator Set for the high sensitivity range is a set of 5 liquid calibrators based on human serum with 5 different concentrations.

PREPARATION

The Calibrators are liquid and ready to use.

STABILITY AND STORAGE

Storage: at 2 – 8 °C
 Stability: up to the expiration date
 Stability after opening: 3 months at 2 – 8 °C
DO NOT FREEZE!

DIRECTIONS FOR USE

- Use these Calibrators in the same manner as specified with the reagent and/or instrument being used.
- Allow the Calibrators to equilibrate to room temperature (22 – 28 °C) prior to use.
- The Calibrators should not be allowed to stand for repeated long periods of time (up to 6 hours) at room temperature. Keep vials tightly capped at all times when not in use to avoid microbial contamination. Keep vials refrigerated when not in use.

WARNINGS AND PRECATIONS

- For in Vitro Diagnostic use.
- Each individual donation intended for use in manufacture of protein calibrator was tested for hepatitis B surface antigen (HBsAg), anti-hepatitis C (anti-HCV) and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.
- Reagents containing sodium azide must be handled with due caution: Do not ingest or allow to contact skin or mucous membranes! Sodium azide can form explosive azides when contacting heavy metals such as copper or lead.

WASTE MANAGEMENT

Please refer to local legal requirements.

LOT SPECIFIC ASSAY DATA

Values and expiry date are lot specific.
 The values of the CRP Calibrator Set are based on IFCC/BCR/CAP reference material for 15 plasma proteins CRM470.

LOT: 16733		EXP: 2014/04	
Level	LOT	Assigned Value	
		mg/dL	mg/L
Level 1	16729	0.050 mg/dL	0.50 mg/L
Level 2	16730	0.245 mg/dL	2.45 mg/L
Level 3	16731	0.845 mg/dL	8.45 mg/L
Level 4	16732	1.245 mg/dL	12.45 mg/L
Level 5	16733	2.045 mg/dL	20.45 mg/L

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