DIACON LIPIDS (Lipid Control Serum Normal)

Lyophilized control serum for the control of accuracy and precision of quantitative in vitro determination of lipids on photometric systems.

REF

Cont.

D99486 **3 x 3 ml** Control Serum
D99486SV **1 x 3 ml** Control Serum

COMPOSITION

Lyophilized human-based serum
 Additives of purified material of human origin

PREPARATION

1. Open the vial very carefully, avoiding any loss of the lyophilized material.
2. Add exactly 3 ml of dist. water (inaccurate reconstitution of the control can cause erroneous results).
3. Close the vial carefully and allow the control to stand for 30 min.
4. Dissolve contents completely by occasionally swirling gently, avoiding the formation of foam. Do not shake!
5. Transfer the quantity needed for determination into a clean sample vial and handle like a patient sample.

Frozen aliquots:

Leave frozen aliquots of the reconstituted control in the dark at room temperature (18 – 25 °C) until they are completely unfrozen. To homogenize, slightly swivel aliquots and immediately afterwards use them for determination.

STABILITY AND STORAGE

Storage: at 2 – 8 °C
 Stability: until the end of the indicated month of expiry.
 Stability after reconstitution (tightly closed):
 15 –25°C 2 – 8 °C - 20°C
 NEFA: 8 hours 7 days not possible
 Other analytes: 8 hours 7 days 30 days
CLOSE IMMEDIATELY AFTER USE! FREEZE ONLY ONCE!

Avoid bacterial contamination. Do not use the product if there is visible evidence of microbial growth in the vial.

Improper handling and/or storage of the control can affect results.

LOT SPECIFIC CONTROL VALUES AND RANGES

The analyte concentrations contained in Diacon Lipids are specific and only valid for the corresponding lot. Please refer to the enclosed table. All values have been established within standardized conditions with the method stated in the value sheet using the corresponding Dialab reagents.

Ranges of acceptance were calculated as assigned value ± the maximum tolerable deviation of a single value according to the Guidelines of the German Federal Medical Council (Rilibäk) from 2003 [3]. For analytes not listed in the Guidelines of the German Federal Medical Council (Rilibäk) ranges are indicated with a deviation of ± 20% from the given mean.

Each laboratory should check if the reference range is transferable to its own conditions and regulatory requirements. Where appropriate, the values stated in this document are to be regarded as a guideline and own mean values and acceptance limits have to be determined.

Each laboratory should establish corrective action in case of deviations in control recovery. These deviations might be attributed to varying laboratory procedures, to instruments, reagents or modified test methods.

TEST PROCEDURE

Please refer to the reagent package insert for instructions for use. Each laboratory has to establish procedures for quality assurance which fulfil all national requirements.

LIMITATION

Compatibility of Diacon Lipids is only guaranteed if those methods stated in the enclosed table are used. When using Diacon Lipids with other methods, proceed with extreme caution. In such case, values found might vary considerably from those stated.

All changes made in method standardization, in applications, reagent compositions or other influences may generate value deviations.

WARNINGS AND PRECAUTIONS

Diacon Lipids is intended for professional in vitro diagnostics.


1. Only blood donations of European origin which were found to be non-reactive when tested with approved methods for HBsAg, anti-HIV 1+2 and anti-HCV were used for the production of Diacon Lipids. Moreover, HCV and HIV were additionally tested by PCR. As there is no possibility to exclude definitely that products derived from human blood transmit infectious agents, it is recommended to handle the control with the same precautions used for patient specimens.
2. Please refer to the safety data sheets and take the necessary precautions for the use of calibrators and controls.

WASTE MANAGEMENT

Please refer to local legal requirements.

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