DIAQUICK DOA Dipsticks

for human urine samples

	REF	Content
DIAQUICK AMP Dipstick	Z02504CE	- 50 Tests (50x REF Z02504B)
DIAQUICK BAR Dipstick	Z02506CE	- 50 Tests (50x REF Z02506B)
DIAQUICK BZO Dipstick	Z02501CE	- 50 Tests (50x REF Z02501B)
DIAQUICK COC Dipstick	Z02510CE	- 50 Tests (50x REF Z02510B)
DIAQUICK ETG Dipstick	Z15101CE	- 50 Tests (50x REF Z15101B)
DIAQUICK LSD Dipstick	Z15105CE	- 50 Tests (50x REF Z15105B)
DIAQUICK MDMA Dipstick	Z05610CE	- 50 Tests (50x REF Z05610B)
DIAQUICK MET Dipstick	Z02500CE	- 50 Tests (50x REF Z02500B)
DIAQUICK MOP Dipstick	Z02505CE	- 50 Tests (50x REF Z02505B)
DIAQUICK MTD Dipstick	Z02550CE	- 50 Tests (50x REF Z02550B)
DIAQUICK PCP Dipstick	Z02560CE	- 50 Tests (50x REF Z02560B)
DIAQUICK TCA Dipstick	Z05045CE	- 50 Tests (50x REF Z05045B)
DIAQUICK THC Dipstick	Z02502CE	- 50 Tests (50x REF Z02502B)
DIAQUICK TRA Dipstick	Z10514CE	- 50 Tests (50x REF Z10514B)
DIAQUICK ZOL Dipstick	Z15103CE	- 50 Tests (50x REF Z15103B)

All tests are individually packed. All products contain a package insert.
 For in vitro diagnostic use only. For use by medical professionals only.
 For diagnosis and therapeutic monitoring only.

INTENDED USE

The DIAQUICK DOA Dipsticks (urine) are rapid, lateral flow chromatographic immunoassays for the qualitative detection of the following drugs and their metabolites:

Parameter	Code	Calibrator Substance	Cut-off
Amphetamine	AMP	d-Amphetamine	1 000 ng/mL
Barbiturates	BAR	Secobarbital	300 ng/mL
Benzodiazepines	BZO	Oxazepam	300 ng/mL
Cocaine	COC	Benzoylcegonine	300 ng/mL
Ethylglucuronide	ETG	Ethyl-β-D-Glucuronide	500 ng/mL
LSD	LSD	Lysergic acid diethylamide	20 ng/mL
Ecstasy	MDMA	(±)-3,4-Methylenedioxyamphetamine HCl	500 ng/mL
Methamphetamine	MET	d-Methamphetamine	1 000 ng/mL
Opiate, Morphine, Heroine	MOP	Morphine	300 ng/mL
Methadone	MTD	Methadone	300 ng/mL
Opiate, Morphine, Heroine	OPI	Morphine	2 000 ng/mL
Phencyclidine	PCP	Phencyclidine	25 ng/mL
Tricyclic Antidepressants	TCA	Nortriptyline	1 000 ng/mL
Marijuana/Cannabis	THC	11-nor-Δ9-THC-9-COOH	50 ng/mL
Tramadol	TRA	cis-Tramadol	100 ng/mL
Zolpidem	ZOL	Zolpidem	50 ng/mL

This test will detect other related compounds; please refer to the Analytical Specificity table in this insert. This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be applied to any drug of abuse test result, particularly when preliminary positive results are obtained.

TEST PRINCIPLE

The DIAQUICK DOA Dipsticks (urine) are immunoassays based on the principle of competitive binding. Drugs which may be present in the urine specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a urine specimen migrates upward by capillary action. A drug, if present in the urine specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody coated on the particles. The antibody coated particles will then be captured by the immobilized drug conjugate and a visible coloured line will show up in the test line region of the specific drug strip. The coloured line will not form in the test line region if the drug level is above its cut-off concentration because it will saturate all the binding sites of the antibody coated on the particles. A drug-positive urine specimen will not generate a coloured line in the specific test line region of the strip because of drug competition, while a drug-negative urine specimen or a specimen containing a drug concentration less than the cut-off will generate a line in the test line region. To serve as a procedural control, a coloured line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

Each test line in the test contains mouse monoclonal antibody-coupled particles and corresponding drug-protein conjugates. A goat antibody is employed in each control line.

WARNINGS AND PRECAUTIONS

- For medical and other in vitro diagnostic use only. Do not use after the expiration date.
- The test strips should remain in the sealed pouch until use.
- All specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.
- The used test strips should be discarded according to federal, state and local regulations.

STORAGE

The DIAQUICK DOA Dipsticks can be stored refrigerated or at room temperature (2 – 30 °C). The tests are stable through the expiration date printed on the sealed pouch. The test strips must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SAMPLE COLLECTION AND PREPARATION

The urine must be collected in a clean and dry container. Urine collected at any time of the day may be used. Urine specimens exhibiting visible precipitations should be centrifuged, filtered or allowed to settle to obtain a clear specimen for testing. Urine specimens may be stored at 2 – 8 °C for up to 48 h prior to testing. For prolonged storage, specimens may be frozen and stored below –20 °C. Frozen specimens should be thawed and mixed before

testing.


MATERIALS REQUIRED BUT NOT PROVIDED

- Specimen collection container
- Timer

ASSAY PROCEDURE

Allow the test strip, urine specimen, and/or controls to equilibrate to room temperature (15 – 30 °C) prior to testing.

- Bring the pouch to room temperature before opening it. Remove the test strip from the sealed pouch and use it as soon as possible.
- With arrows pointing toward the urine specimen, immerse the test strip vertically in the urine specimen for at least 10-15 seconds. Do not pass the maximum line (MAX) on the test strip when immersing the strip. See the illustration below.
- Place the test strip on a non-absorbent flat surface, start the timer and wait for the colored line(s) to appear. Read results at 5 minutes. Do not interpret the result after 10 minutes.

INTERPRETATION OF RESULTS

NEGATIVE: Two lines appear. One coloured line should be in the control line region (C), and another visibly coloured line should be in the test line region (T). This negative result indicates that the drug concentration is below the detectable level.

***NOTE:** The shade of colour in the test region (T) may vary, but it should be considered negative whenever there is even a faint coloured line.

POSITIVE: One coloured line appears in the control line region (C). No line appears in the test line region (T). This positive result indicates that the drug concentration exceeds the detectable level.

INVALID: Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test. If the problem persists, discontinue using the lot immediately and contact your local distributor.

QUALITY CONTROL

A procedural control is included in the test. A coloured line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

LIMITATIONS

- The DIAQUICK DOA Dipsticks provide only a preliminary analytical result. A more specific chemical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.^{1,2}
- It is possible that technical or procedural errors, as well as other interfering substances in the urine specimen may cause erroneous results.
- Adulterants, such as bleaching agents in urine specimens may produce erroneous results regardless of the analytical method used. If adulteration is suspected, the test should be repeated with another urine specimen.
- A positive result indicates presence of the drug or its metabolites but does not indicate the level of intoxication, administration route or concentration in urine.
- A negative result may not necessarily indicate drug-free urine. Negative results can be obtained if a drug is present but below the cut-off level of the test.
- The DIAQUICK DOA Dipsticks do not distinguish between drugs of abuse and certain medications.
- A positive result might be obtained from certain foods or food supplements.

PERFORMANCE CHARACTERISTICS

ACCURACY

A side-by-side comparison of the DIAQUICK DOA Dipsticks and a commercially available rapid drug test was conducted. Testing was performed on approx. 100 specimens previously collected from subjects present for drug screen testing. The agreement was > 99.9 % for all tests.

A side-by-side comparison of the DIAQUICK DOA Dipsticks and GC/MS at the cut-off level of the tests was conducted. Testing was performed on 250 specimens previously collected from subjects present for drug screen testing. The following results were tabulated:

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ANALYTICAL SPECIFICITY

The following table lists the concentration of compounds (ng/mL) that are detected positive in urine by the DIAQUICK DOA Dipsticks at 5 minutes.

AMPHETAMINE	AMP	BARBITURATES	BAR
D,L-Amphetamine sulfate	300	Amobarbital	5 000
L-Amphetamine	25 000	5,5-Diphenylhydantoin	8 000
(±) 3,4-Methylenedioxyamphetamine	500	Allobarbitol	600
Phentermine	800	Barbital	8 000
Maprotiline	50 000	Talbutal	200
Methoxyphenamine	6 000	Butalbitol	8 000
D-Amphetamine	1 000	Phenobarbital	300
BENZODIAZEPINES	BZO	Cyclopentobarbital	30 000
Alprazolam	100	Pentobarbital	8 000
a-hydroxyalprazolam	1 500	Alphenol	600
Bromazepam	900	Aprobarbital	500
Chlordiazepoxide	900	Butabarbitol	200
Clobazam	200	Butethal	500
Clonazepam	500	Secobarbital	300
Clorazepate dipotassium	500	COCAINE	COC
Delorazepam	900	Benzoyllecgonine	300
Desalkylflurazepam	200	Cocaine HCl	200
Diazepam	300	Cocacethylene	20 000
Estazolam	6 000	Ecgonine HCl	30 000
Flunitrazepam	200	ETHYLGLUCURONIDE	ETG
(±) Lorazepam	3 000	Ethyl-β-D-Glucuronide	500
RS-Lorazepam glucuronide	200	Propyl-β-D-Glucuronide	50 000
Midazolam	6 000	Morphine-3-β-Glucuronide	100 000
Nitrazepam	200	Morphine-6-β-Glucuronide	100 000
Norchlordiazepoxide	100	Glucuronic Acid	100 000
Nordiazepam	900	Ethanol	100 000
Oxazepam	300	Methanol	100 000
Temazepam	100	LSD	LSD
Triazolam	3 000	Lysergic Acid Diethylamide	20
MORPHINE	MOP	Fentanyl	30
Codeine	200	ECSTASY	MDMA
Ethylmorphine	6 000	(±) 3,4-Methylenedioxyamphetamine HCl	500
Hydrocodone	50 000	(±) 3,4-Methylenedioxyamphetamine HCl (MDA)	3 000
Hydromorphone	3 000	3,4-Methylenedioxyethylamphetamine (MDE)	300
Levorphanol	1 500	METHAMPHETAMINE	MET
6-Monoacetylmorphine	300	p-Hydroxymethamphetamine	25 000
Morphine 3-β-D-glucuronide	800	D-Methamphetamine	1 000
Morphine	300	L-Methamphetamine	20 000
Norcodeine	6 000	(±)-3,4-Methylenedioxyamphetamine	12 500
Normorphine	50 000	Mephentermine	50 000
Oxycodone	30 000	METHADONE	MTD
Oxymorphone	50 000	Methadone	300
Procaine	15 000	Doxylamine	100 000
Thebaine	6 000	Cis-tramadol	300 000
OPIATES	OPI	TRICYCLIC ANTIDEPRESSANTS	TCA
Codeine	2 000	Nortriptyline	1 000
Ethylmorphine	3 000	Nordoxepine	500
Hydrocodone	50 000	Trimipramine	3 000
Hydromorphone	15 000	Amitriptyline	1 500
Levorphanol	25 000	Promazine	3 000
6-Monoacetylmorphine	3 000	Desipramine	200
Morphine 3-β-D-glucuronide	2 000	Cyclobenzaprine	2 000
Morphine	2 000	Imipramine	400
Norcodeine	25 000	Clomipramine	50 000
Normorphine	50 000	Doxepine	2 000
Oxycodone	25 000	Maprotiline	2 000
Oxymorphone	25 000	Promethazine	50 000
Procaine	50 000	Perphenazine	50 000
Thebaine	25 000	Dithiaden	10 000
CANNABIS	THC	TRAMADOL	TRA
Cannabinol	35 000	n-Desmethyl-cis-tramadol	200
11-nor-Δ ⁹ -THC-9 COOH	30	Cis-tramadol	100
11-nor-Δ⁹-THC-9 COOH	50	Procyclidine	100 000
Δ ⁹ -THC	17 000	o-Desmethyl-cis-tramadol	10 000
Δ ⁸ -THC	17 000	Phencyclidine	100 000
ZOLPIDEM	ZOL	d,l-O-Desmethyl venlafaxine	50 000
Zolpidem	50	PHENCYCLIDINE	PCP
		4-Hydroxyphencyclidine	12 500
		Phencyclidine	25

CROSS-REACTIVITY


A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free or drug positive urine. The following compounds did not show a cross-reactivity when tested with the DIAQUICK DOA Dipsticks at a concentration of 100 µg/mL.

Non Cross-Reacting Compounds:

Acetophenetidin	Cortisone	Zomepirac	d-Pseudoephedrine
N-Acetylprocainamide	Creatinine	Ketoprofen	Quinidine
Acetylsalicylic acid	Deoxycorticosterone	Labelalol	Quinine
Aminopyrine	Dextromethorphan	Loperamide	Salicylic acid
Amoxicillin	Diclofenac	Meprobamate	Serotonin
Ampicillin	Diflunisal	Methoxyphenamine	Sulfamethazine
l-Ascorbic acid	Digoxin	Methylphenidate	Sulindac
Apomorphine	Diphenhydramine	Nalidixic acid	Tetracycline
Aspartame	Ethyl-p-aminobenzoate	Naproxen	Tetrahydrocortisone,
Atropine	β-Estradiol	Niacinamide	3-acetate
Benzilic acid	Estrone-3-sulfate	Nifedipine	Tetrahydrocortisone
Benzoic acid	Erythromycin	Norethidrone	Tetrahydrozoline
Bilirubin	Fenoprofen	Noscapine	Thiamine
d,l-Brompheniramine	Furosemide	d,l-Octopamine	Thioridazine
Caffeine	Gentisic acid	Oxalic acid	d,l-Tyrosine
Cannabidiol	Hemoglobin	Oxolinic acid	Tolbutamide
Chloral hydrate	Hydralazine	Oxymetazoline	Triamterene
Chloramphenicol	Hydrochlorothiazide	Papaverine	Trifluoperazine
Chlorothiazide	Hydrocortisone	Penicillin-G	Trimethoprim
d,l-Chlorpheniramine	o-Hydroxyhippuric acid	Perphenazine	d,l-Tryptophan
Chlorpromazine	3-Hydroxytyramine	Phenelzine	Uric acid
Cholesterol	d,l-Isoproterenol	Prednisone	Verapamil
Clonidine	Isoxsuprine	d,l-Propranolol	

REFERENCES

- Baselt, RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company, 1986; 1735
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

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	REF	Inhalt
DIAQUICK AMP Dipstick	Z02504CE	- 50 Tests (50x REF Z02504B)
DIAQUICK BAR Dipstick	Z02506CE	- 50 Tests (50x REF Z02506B)
DIAQUICK BZO Dipstick	Z02501CE	- 50 Tests (50x REF Z02501B)
DIAQUICK COC Dipstick	Z02510CE	- 50 Tests (50x REF Z02510B)
DIAQUICK ETG Dipstick	Z15101CE	- 50 Tests (50x REF Z15101B)
DIAQUICK LSD Dipstick	Z15105CE	- 50 Tests (50x REF Z15105B)
DIAQUICK MDMA Dipstick	Z05610CE	- 50 Tests (50x REF Z05610B)
DIAQUICK MET Dipstick	Z02500CE	- 50 Tests (50x REF Z02500B)
DIAQUICK MOP Dipstick	Z02505CE	- 50 Tests (50x REF Z02505B)
DIAQUICK MTD Dipstick	Z02550CE	- 50 Tests (50x REF Z02550B)
DIAQUICK PCP Dipstick	Z02560CE	- 50 Tests (50x REF Z02560B)
DIAQUICK TCA Dipstick	Z05045CE	- 50 Tests (50x REF Z05045B)
DIAQUICK THC Dipstick	Z02502CE	- 50 Tests (50x REF Z02502B)
DIAQUICK TRA Dipstick	Z10514CE	- 50 Tests (50x REF Z10514B)
DIAQUICK ZOL Dipstick	Z15103CE	- 50 Tests (50x REF Z15103B)

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