# ELECTROLYTE CONTROL

## CONTROLS FOR DETERMINATION OF ELECTROLYTES CONCENTRATION

Kit name	Kit size	Cat. No
ELECTROLYTE CONTROL LEVEL 1	30 x 1.8 ml	E-717
ELECTROLYTE CONTROL LEVEL 2	30 x 1.8 ml	E-718
ELECTROLYTE CONTROL LEVEL 3	30 x 1.8 ml	E-719
ELECTROLYTE CONTROL TRILEVEL	30 x 1.8 ml	E-705

### INTENDED USE

ELECTROLYTE CONTROL sets are intended for quality control of electrolytes measurements in ISE analysers.

ELECTROLYTE CONTROL are a buffered solution of electrolytes ( $\text{Na}^+$ ,  $\text{K}^+$ ,  $\text{Cl}^-$ ,  $\text{Ca}^{2+}$ ,  $\text{Li}^+$ ) and contains no human-based materials.

### STABILITY AND STORAGE

- Unopened controls are stable at 18-25°C until expiry date printed on the product label.
- Avoid freezing and exposure to temperatures greater than 30°C.
- Controls can be also store at 4-25°C. Then allow the controls to equilibrate to room temperature prior to use.
- Do not store once opened ampoules.

### PROCEDURE

Controls are ready for use.

- Mix controls gently before use.
- After opening immediately introduce the liquid from the ampoule to the analyser, following the instrument instruction for sampling a control material.

### NOTES

- Product for in vitro diagnostic use only.
- Controls should be used by suitably qualified laboratory personnel only in accordance with intended purpose.
- These controls are not intended for use as a calibrator.
- Do not use controls after expiry date.
- Controls are not a blood-based material, it may not detect certain malfunctions, which would affect the testing of blood.

### ASSIGNED VALUES

Assigned Values are calculated from results of multiple determinations performed on randomly selected samples from each lot. It is recommended that each laboratory establish its own means and acceptable ranges. The mean value established should fall within the Assigned Values shown on the Values Sheet.

### WASTE MANAGEMENT

Please refer to local legal requirements.

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**PZ CORMAY S.A.**