Instructions for Use

17-OH Progesterone ELISA

IVD

REF EIA-1292

96

DRG Instruments GmbH, Germany
Frauenbergstraße 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the Instructions for Use provided with the kit.
 Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
 Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
 Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.
 Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos / Table des matières

1	INTENDED USE.....	2	1	FINALIDAD PREVISTA.....	26
2	PRINCIPLE OF THE TEST	2	2	FUNDAMENTO DEL ENSAYO	26
3	WARNINGS AND PRECAUTIONS	3	3	PRECAUCIONES	26
4	REAGENTS.....	4	4	COMPONENTES DEL KIT.....	27
5	SPECIMEN COLLECTION AND PREPARATION.....	5	5	MUESTRAS	28
6	ASSAY PROCEDURE.....	5	6	PROCEDIMIENTO DE ENSAYO	28
7	EXPECTED NORMAL VALUES.....	7	7	VALORES ESPERADOS.....	30
8	QUALITY CONTROL.....	7	8	CONTROL DE CALIDAD	30
9	PERFORMANCE CHARACTERISTICS.....	8	9	CARACTERÍSTICAS DEL ENSAYO.....	31
10	LIMITATIONS OF USE.....	9	10	LIMITACIONES DE USO	31
11	LEGAL ASPECTS	10	11	ASPECTOS LEGALES	32
1	ZWECKBESTIMMUNG	11	1	UTILISATION PRÉVUE	33
2	TESTPRINZIP	11	2	PRINCIPE DU TEST	33
3	VORSICHTSMAßNAHMEN.....	12	3	PRECAUTIONS D'UTILISATION.....	33
4	BESTANDTEILE DES KITS	13	4	COMPOSITION DU KIT	34
5	PROBENVORBEREITUNG.....	14	5	ECHANTILLON	35
6	TESTDURCHFÜHRUNG	14	6	RÉALISATION DU TEST	35
7	ERWARTETE WERTE	16	7	VALEURS ATTENDUES.....	37
8	QUALITÄTSKONTROLLE	16	8	CONTROLE DE QUALITE	37
9	ASSAY-CHARAKTERISTIKA.....	17	9	CARACTERISTIQUES DU TEST	38
10	GRENZEN DES TESTS	17	10	LIMITES D'UTILISATION.....	38
11	RECHTLICHE GRUNDLAGEN	18	11	ASPECTS LEGAUX.....	39
1	DESTINAZIONE D'USO	19	12	REFERENCES / LITERATURE	40
2	PRINCIPIO DEL TEST	19			
3	PRECAUZIONI.....	19			
4	COMPONENTI DEL KIT.....	20			
5	CAMPIONI.....	21			
6	ATTUAZIONE DEL TEST.....	21			
7	VALORI NORMALI.....	23			
8	CONTROLLO QUALITÀ.....	23			
9	CARATTERISTICHE DEL TEST	24			
10	LIMITAZIONE DEL TEST	24			
11	ASPETTI LEGALI	25			
				SYMBOLS USED.....	41

1 INTENDED USE

The **DRG 17-OH Progesterone ELISA** is an enzyme immunoassay for the quantitative *in vitro diagnostic* measurement of 17- α -OH Progesterone (17- α -OHP) in serum or plasma (EDTA, lithium heparin or citrate plasma).

1.1 Summary and Explanation

The steroid 17- α -Hydroxyprogesterone (17- α -OHP) is produced by both the adrenal cortex and gonads. Even though 17- α -OHP has relatively low progestational activity, it is of intense clinical interest because it is the immediate precursor to 11-desoxycortisol (Cpd-S). Because Cpd-S is produced by 21-hydroxylation of 17- α -OHP, measurement of 17- α -OHP is a useful indirect indicator of 21-hydroxylase activity. In congenital 21-hydroxylase deficiency, the most common variety of congenital adrenal hyperplasia (CAH), 17- α -OHP is secreted in abundant excess. It is moderately elevated in the 11- β -hydroxylase deficiency as well. Measurement of 17- α -OHP is therefore valuable in the initial diagnosis of CAH.

1.2 Clinical Physiology

Adult non-pregnant women:

In adult non-pregnant women in the childbearing age group, 17- α -OHP concentrations vary over the menstrual cycle with luteal phase concentrations being higher than follicular phase concentrations. This is because 17- α -OHP is secreted parallel with progesterone from maturing follicles or from the corpus luteum. There is also a diurnal variation of 17- α -OHP concentrations.

This rhythm is parallel with adrenal cortisol secretion such that maximum 17- α -OHP concentrations are measured in samples obtained between midnight and 8:00 am.

Adult males:

There is little information available on the systematic variability of 17- α -OHP concentration in adult males.

Pregnant women and newborn children:

The steroid 17- α -OHP is produced in large amounts by the fetus and the adrenals. It is secreted in abundance into both the fetal and maternal circulation. The maternal concentrations of 17- α -OHP increase very sharply after 32 weeks gestational age to about 4-fold above basal concentrations at term.

1.3 Clinical Applications

Congenital adrenal hyperplasia:

The principal application of the 17- α -OHP ELISA is in the diagnosis of CAH in newborns with ambiguous genitalia and in virilized adolescent girls. Since 17- α -OHP is the immediate precursor to 11-desoxycortisol, basal 17- α -OHP concentrations are sharply elevated in patients with 21-hydroxylase deficiency and to a lesser degree in patients with 11-hydroxylase deficiency.

Because 17- α -OHP concentrations are so markedly elevated in newborns and adolescent girls afflicted with CAH, a single basal measurement is all that is normally required to make the diagnosis.

Late onset adrenal hyperplasia:

More recently, 17- α -OHP concentrations have been utilized in the evaluation of androgenized women where late onset of 21-hydroxylase deficiency is suspected. This condition is clinically very subtle and since the presentation is the same as classical polycystic ovarian disease, basal plasma 17- α -OHP concentrations, unlike classical congenital adrenal hyperplasia, are normal. The diagnosis is made by administration of an ACTH stimulation test.

Other applications:

Measurement of 17- α -OHP concentrations is also utilized in evaluation of both men and women with acne vulgaris, male pattern baldness and in some subtle forms of infertility. Experiences with these applications are very limited.

2 PRINCIPLE OF THE TEST

The DRG 17-OH Progesterone ELISA is a solid phase enzyme-linked immunosorbent assay (ELISA), based on the **principle of competitive binding**.

The microtiter wells are coated with a polyclonal [rabbit] antibody directed towards antigenic sites of the 17- α -OHP molecule.

Samples are pre-incubated in the coated wells.

During the second incubation, 17- α -OHP in the added sample competes with the added enzyme conjugate, which is 17- α -OHP conjugated to horseradish peroxidase, for binding to the coated antibody.

After a washing step to remove all unbound substances, the solid phase is incubated with the substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured. The intensity of colour is inversely proportional to the concentration of the analyte in the sample.

A standard curve is constructed by plotting OD values against concentrations of standards, and concentrations of unknown samples are determined using this standard curve.

3 WARNINGS AND PRECAUTIONS

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of the package insert provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (21 °C to 26 °C) before starting the test. Temperature will affect the absorbance readings of the assay.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with *Stop Solution* containing 0.5 M H₂SO₄. It may cause skin irritation and burns.
18. Some reagents contain Proclin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

4 REAGENTS

4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;
Wells coated with anti-17- α -OHP antibody (polyclonal).
2. **Standard (Standard 0 - 6)**, 7 vials, 1 mL each, ready to use;
Concentrations: 0; 0.15; 0.5; 1.5; 3; 7.5; 20 ng/mL
0; 0.45; 1.5; 4.5; 9.1; 22.7; 60.6 nmol/L
Conversion: ng/mL \times 3.03 = nmol/L.
The standards are calibrated against the following reference material: *Certified Reference Material Cerilliant H-085*
Contain non-mercury preservative.
3. **Control Low & High**, 2 vials, 1 mL each, ready to use;
For control values and ranges please refer to vial label or Certificate of Analysis.
Contain non-mercury preservative.
4. **Enzyme Conjugate**, 1 vial, 25 mL, ready to use,
17- α -OHP conjugated to horseradish peroxidase;
Contains non-mercury preservative.
5. **Substrate Solution**, 1 vial, 25 mL, ready to use,
Tetramethylbenzidine (TMB).
6. **Stop Solution**, 1 vial, 14 mL, ready to use,
contains 0.5 M H₂SO₄,
Avoid contact with the stop solution. It may cause skin irritations and burns.
7. **Wash Solution**, 1 vial, 30 mL (40X concentrated);
See "Reagent Preparation".

Note: Additional *Standard 0* for sample dilution is available upon request.

4.2 Materials required but not provided

- A calibrated microtiter plate reader (450 nm, with reference wavelength at 620 nm to 630 nm)
- Calibrated variable precision micropipettes
- Manual or automatic equipment for rinsing microtiter plates
- Absorbent paper
- Distilled water
- Timer
- Semi logarithmic graph paper or software for data reduction

4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date.

Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

Opened kits retain activity for 8 weeks if stored as described above.

4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (20 °C to 26 °C) prior to use.

Wash Solution

Add distilled water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL distilled water to a final volume of 1200 mL.

The diluted Wash Solution is stable for 1 week at room temperature.

4.5 Disposal of the Kit

The disposal of the kit must be made according to the national regulations. Special information for this product is given in the Safety Data Sheet.

4.6 Damaged Test Kits

In case of any severe damage to the test kit or components, DRG has to be informed in writing, at the latest, one week after receiving the kit. Severely damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed according to the official regulations.

5 SPECIMEN COLLECTION AND PREPARATION

Serum or plasma (EDTA, lithium heparin or citrate plasma) can be used in this assay.

Please note: Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use haemolytic, icteric, or lipaemic specimens. For further information refer to chapter "Interfering Substances".

5.1 Specimen Collection

Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

Plasma:

Whole blood should be collected into centrifuge tubes containing anti-coagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 7 days at 2 °C to 8 °C prior to assaying.

Specimens stored for a longer time (up to 12 months) should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

5.3 Specimen Dilution

If in an initial assay, a specimen is found to contain more than the highest standard, the specimens can be diluted with *Standard 0* and re-assayed as described in Assay Procedure.

For the calculation of the concentrations this dilution factor has to be taken into account.

Example:

- a) dilution 1:10: 10 µL sample + 90 µL *Standard 0* (mix thoroughly)
b) dilution 1:100: 10 µL dilution a) 1:10 + 90 µL *Standard 0* (mix thoroughly).

6 ASSAY PROCEDURE

6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Optical density is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.

6.2 Test Procedure

Each run must include a standard curve.

1. Secure the desired number of Microtiter wells in the frame holder.
2. Dispense **25 µL** of each **Standard, Control** and **sample** with new disposable tips into appropriate wells.
3. Incubate for **5 minutes** at room temperature
4. Dispense **200 µL Enzyme Conjugate** into each well.
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
5. Incubate for **60 minutes** at room temperature.
6. Briskly shake out the contents of the wells.
Rinse the wells **3 times** with diluted Wash Solution (400 µL per well). Strike the wells sharply on absorbent paper to remove residual droplets.
Important note:
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
7. Add **200 µL of Substrate Solution** to each well.
8. Incubate for **30 minutes** at room temperature.
9. Stop the enzymatic reaction by adding **100 µL of Stop Solution** to each well.
10. Measure the optical density of the solution in each well at **450 nm (reading) and at 620 nm to 630 nm (background subtraction, recommended)** with a microtiter plate reader.
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

6.3 Calculation of Results

1. Calculate the average optical density (OD) values for each set of standards, controls and patient samples.
2. Using semi-logarithmic graph paper, construct a standard curve by plotting the mean OD obtained from each standard against its concentration with OD value on the vertical (Y) axis and concentration on the horizontal (X) axis.
3. Using the mean OD value for each sample determine the corresponding concentration from the standard curve.
4. Automated method: The results in the Instructions for Use have been calculated automatically using a 4-Parameter curve fit. (4 Parameter Rodbard or 4 Parameter Marquardt are the preferred methods.) Other data reduction functions may give slightly different results.
5. The concentration of the samples can be read directly from this standard curve. Samples with concentrations higher than that of the highest standard have to be further diluted or reported as > 20 ng/mL. For the calculation of the concentrations this dilution factor has to be taken into account.

6.3.1 Example of Typical Standard Curve

The following data is for demonstration only and **cannot** be used in place of data generations at the time of assay.

Standard	Optical Density (450 nm)
Standard 0 (0 ng/mL)	2.15
Standard 1 (0.15 ng/mL)	1.77
Standard 2 (0.5 ng/mL)	1.28
Standard 3 (1.5 ng/mL)	0.77
Standard 4 (3.0 ng/mL)	0.49
Standard 5 (7.5 ng/mL)	0.25
Standard 6 (20 ng/mL)	0.12

7 EXPECTED NORMAL VALUES

It is strongly recommended that each laboratory should determine its own normal and abnormal values.

In a study conducted with newborns and children, using the DRG 17-OH Progesterone ELISA, the following values are observed:

Newborns (boys and girls)	n		Range (min - max) (ng/mL)	Mean (ng/mL)	Median (ng/mL)	2.5 th - 97.5 th Percentile (ng/mL)
	26	1 st month after birth	0 - 17.3	7.2	6.7	1.0 - 17.0
	43	2 nd month after birth	0.32 - 13.7	4.9	4.6	1.6 - 9.8
	21	3 rd month after birth	0.06 - 4.2	2.3	2.3	0.5 - 4.1
	12	4 th month after birth	0.2 - 4.6	2.1	2.3	0.2 - 4.3

	n	Age (years)	Range (min - max) (ng/mL)	Mean (ng/mL)	Median (ng/mL)	2.5 th - 97.5 th Percentile (ng/mL)
Children	75	1 - 10	0.03 - 2.85	1.04	0.88	0.08 - 2.58
Adolescent	3	11 - 14	0.06 - 1.38	0.65	0.50	0.07 - 1.34
	10	15 - 18	0.41 - 2.35	1.24	1.26	0.42 - 2.26

In a study conducted with apparently normal healthy donors, using the DRG 17-OH Progesterone ELISA, the following values are observed:

Adult females	Follicular phase	0.1 - 0.8 ng/mL
	Luteal phase	0.6 - 2.3 ng/mL
	Ovulation	0.3 - 1.4 ng/mL
	Post ACTH	< 3.2 ng/mL
	Third trimester	2.0 - 12 ng/mL
	Postmenopausal	0.13 - 0.51 ng/mL
Adult males		0.5 - 2.1 ng/mL

The results alone should not be the only reason for any therapeutic consequences. The results should be correlated to other clinical observations and diagnostic tests.

8 QUALITY CONTROL

Good laboratory practice requires that controls be run with each calibration curve. A statistically significant number of controls should be assayed to establish mean values and acceptable ranges to assure proper performance.

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

The controls and the corresponding results of the QC-Laboratory are stated in the QC certificate added to the kit. The values and ranges stated on the QC sheet always refer to the current kit lot and should be used for direct comparison of the results.

It is also recommended to make use of national or international Quality Assessment programs in order to ensure the accuracy of the results.

Employ appropriate statistical methods for analysing control values and trends. If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

9 PERFORMANCE CHARACTERISTICS

9.1 Assay Dynamic Range

The range of the assay is between 0.042 – 20 ng/mL.

9.2 Specificity of Antibodies (Cross Reactivity)

The following substances were tested for cross-reactivity of the assay:

Substance	Conc. Range (ng/mL)	Mean Bias (ng/mL)	Mean Bias %
17-Benzoate Estradiol	2 - 2000	0.05	< 0.01
17-Cypionate Estradiol	2 - 2000	< 0.01	< 0.01
17-Valerate Estradiol	2 - 2000	< 0.01	< 0.01
Aldosterone	2 - 2000	0.06	< 0.01
Androstenedione	2 - 2000	< 0.01	< 0.01
Corticosterone	2 - 200	< 0.01	< 0.01
Cortisol	2 - 200	< 0.01	< 0.01
Cortisone	2 - 200	< 0.01	< 0.01
DHEA	2 - 2000	< 0.01	< 0.01
DHEA-S(a)	2 - 2000	< 0.01	< 0.01
Estradiol	2 - 2000	< 0.01	< 0.01
Estriol	2 - 2000	< 0.01	< 0.01
Estrone	2 - 2000	< 0.01	< 0.01
Progesterone	2 - 20	0.05	1.31
Testosterone	2 - 2000	< 0.01	< 0.01

9.3 Sensitivity

The analytical sensitivity of the DRG ELISA was calculated by subtracting 2 standard deviations from the mean OD of 20 replicate analyses of *Standard 0* and was found to be 0.013 ng/mL.

The Limit of Blank (LoB) is 0.006 ng/mL.

The Limit of Detection (LoD) is 0.042 ng/mL.

The Limit of Quantification (LoQ) is 0.156 ng/mL.

9.4 Reproducibility

9.4.1 Intra Assay

The within assay variability is shown below:

Sample	n	Mean (ng/mL)	CV (%)
1 (Li-Heparin plasma)	10	0.23	4.6
2 (Serum)	10	2.10	4.0
3 (Serum)	10	7.74	3.0
4 (Citrate plasma)	10	12.11	4.2

9.4.2 Inter Assay

The between assay variability is shown below:

Sample	n	Mean (ng/mL)	CV (%)
1 (Li-Heparin plasma)	30	0.21	8.7
2 (Serum)	30	1.99	6.3
3 (Serum)	30	7.41	6.3
4 (Citrate plasma)	30	11.53	6.2

9.4.3 Inter-Lot

The inter-assay (between-lots) variation was determined by measuring each sample 6 times with 3 different kit lots:

Sample	n	Mean (ng/mL)	CV (%)
1	18	0.303	9.4
2	18	0.864	2.5
3	18	1.883	5.7
4	18	7.221	6.7

9.5 Recovery

Samples have been spiked by adding 17- α -OHP solutions with known concentrations.

The recovery (%) was calculated by multiplying the ratio of measured and expected values with 100.

	Li-Heparin plasma	Citrate plasma	EDTA plasma	Serum
Concentration (ng/mL)	0.17	0.46	0.72	1.92
Average Recovery (%)	97.9	96.0	95.7	93.8
Range of Recovery (%)	from	92.6	89.2	93.3
	to	103.8	103.2	98.6

9.6 Linearity

Samples were measured undiluted and in serial dilutions with standard 0. The recovery (%) was calculated by multiplying the ratio of expected and measured values with 100.

	Citrate plasma	Serum	EDTA plasma	Li-Heparin plasma
Concentration (ng/mL)	1.31	6.99	10.90	15.20
Average Recovery (%)	104.1	104.0	108.0	105.7
Range of Recovery (%)	from	93.9	92.2	102.4
	to	111.4	108.8	113.8

10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

10.2 Drug Interferences

Substance	Concentration	Mean Bias
	ng/mL	ng/mL
Coumestrol	5 – 500	-0.34
Daidzein	5 – 500	0.59
Ethisterone	23.4 – 2340	-0.48
Fulvestrant	5 – 500	-0.11
Genistein	5 – 500	0.28
Levonorgestrel	0.2 – 20	0.27
Mifepristone	23.4 – 2340	-0.53
Prednisolone	0.2 – 20	-0.12
Prednisone	0.2 – 20	-0.59
Secoisolariciresinol	5 – 500	0.23

Daidzein will increase the measured 17-OH Progesterone concentration in a sample on average by more than 0.5 ng/mL. Mifepristone and Prednisolone will decrease the measured 17-OH Progesterone concentration in a sample on average by more than 0.5 ng/mL.

All other tested substances will not change the 17-OH Progesterone concentration by more than ± 0.5 ng/mL.

10.3 High-Dose-Hook Effect

"High Dose Hook Effect" is not detected in the range between 0 - 640 ng/mL.

11 LEGAL ASPECTS**11.1 Reliability of Results**

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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


12 REFERENCES / LITERATURE

1. Chrousos GP et al. Late onset 21-hydroxylase deficiency mimicking idiopathic hirsutism or polycystic ovarian disease. *Annals Intern Med.* 1982; 96,143-148.
2. Choi J-H, and Yoo H-W. Management issues of congenital adrenal hyperplasia during transition from pediatric to adult care. *Korean J. Pediatr.* 2017; 60(2), 31-37.
3. Eisenhoefer G et al. Reference intervals for plasma concentrations of adrenal steroids measured by LC-MS/MS: Impact of gender, age, oral contraceptives, body mass index and blood pressure status. *Clin. Chim Acta.* 2017; 470, 115-124.
4. Honour JW. 17-Hydroxyprogesterone in children, adolescents and adults. *Annals of Clinical Biochemistry,* 2014; 51(4), 424-440.
5. Hughes IA, Riad-Fahmy D, and Griffith K. Plasma 17OH-progesterone concentrations in newborn infants. *Archives of Disease in Childhood.* 1979; 54, 347-349.
6. Ishii T et al. Guidelines for diagnosis and treatment of 21-hydroxylase deficiency (2014 revision). *Clin. Pediatr. Endocrinol.* 2015; 24(3), 77-105.
7. Kang MJ et al. Relationships of basal level of serum 17-hydroxyprogesterone with that of serum androstenedione and their stimulated responses to a low dose of ACTH in young adult patients with congenital adrenal hyperplasia due to 21-hydroxylase deficiency. *J. Korean Med. Sci.* 2011; 26(11), 1454-60.
8. Maas KH et al. Relationship between 17-Hydroxyprogesterone responses to human chorionic gonadotropin and markers of ovarian follicle morphology in women with polycystic ovary syndrome. *J. Clin. Endocrinol. Metab.* 2015; 100(1), 293-300.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitätskennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico in vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	No. di Cat.	No de catálogo	Référence
	Batch code *	Fertigungslosnummer, Charge *	Lotto no	Número de lote	No. de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques
	Caution *	Achtung *	Attenzione	Precaución	Attention
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement
<i>Volume/No.</i>	Volume / No.	Volumen / Anzahl	Volume / Quantità	Volumen / Número	Volume / Quantité
<i>Microtiterwells</i>	Microtiterwells	Mikrotiterwells	Micropozzetti	Placas multipocillo	Microplaques
<i>Antiserum</i>	Antiserum	Antiserum	Antisiero	Antisuero	Antisérum
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Tracciante enzimatico	Conjugado enzimático	Conjugué enzymatique
<i>Enzyme Complex</i>	Enzyme Complex	Enzymkomplex	Complesso enzimatico	Complejo enzimático	Complexe enzymatique
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat
<i>Stop Solution</i>	Stop Solution	Stopplösung	Soluzione d'arresto	Solución de parada	Solution d'arrêt
<i>Zero Standard</i>	Zero Standard	Nullstandard	Standard zero	Estándar cero	Zero Standard
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle
<i>Assay Buffer</i>	Assay Buffer	Assaypuffer	Tampone del test	Tampón de ensayo	Tampon d'essai
<i>Wash Solution</i>	Wash Solution	Waschlösung	Soluzione di lavaggio	Solución de lavado	Solution de lavage
<i>1N NaOH</i>	1N NaOH	1N NaOH	1N NaOH (idrossido di sodio 1N)	1N NaOH	1N NaOH
<i>1 N HCl</i>	1 N HCl	1 N HCl		1 N HCl	1N HCl
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungsmedium	Diluyente dei campioni	Solución para dilución de la muestra	Solution pour dilution de l'échantillon
<i>Conjugate Diluent</i>	Conjugate Diluent	Konjugatverdünnungsmedium	Diluyente del tracciante	Solución para dilución del conjugado	Solution pour dilution du conjugué