Instructions for Use

TPS (Cyk18) ELISA

IVD

REF EIA-2354

 96

DRG 

DRG Instruments GmbH, Germany
Frauenbergstraße. 18, 35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG 

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the Instructions for Use provided with the kit.
 Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
 Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
 Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.
 Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

Table of Contents / Inhaltsverzeichnis / Tabella die Contenuti / Tabla de Contenidos

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1 INTENDED USE

TPS (Cyk18) ELISA is an *in vitro* immunoassay intended for quantitative determination of soluble cytokeratin 18 fragments in human serum.

The result may be valuable as an aid in diagnosis, prognosis and follow up of cancer patients. Assays for TPS are frequently used to manage patients with malignancies such as epithelial breast, prostate and ovarian cancers.

2 SUMMARY AND EXPLANATION OF THE ASSAY

Cytokeratins (CK) are intracellular proteins expressed in most epithelial cells. The cytokeratin pattern is usually preserved during transformation of normal cells into malignant cells. Soluble fragments of cytokeratins are released into serum and other body fluids due to apoptosis or necrosis. CK8, CK18 and CK19 are among the most abundant cytokeratins in simple epithelial cells. The quantity of CK fragments is correlated to cell growth and death. An elevated level of soluble cytokeratin fragments in patient samples is an indication of epithelial tumor cell activity.

TPS (Cyk18) ELISA uses specific monoclonal antibodies that bind to defined epitopes on soluble fragments of CK18 in patient samples. The values generated with TPS (Cyk18) ELISA shall be interpreted in conjunction with all available clinical information and diagnostic data for the patient.

3 PRINCIPLE OF THE ASSAY

TPS (Cyk18) ELISA is a one-step enzyme linked sandwich immunoassay. Standards, controls and samples react during incubation simultaneously with a solid phase monoclonal catcher antibody and the HRP-conjugated monoclonal detector antibody (M3). After washing, the TMB substrate is added and after an incubation time the reaction is stopped and the absorbance at 450 nm is measured. The developed color is directly proportional to the concentration of the analyte.

4 SPECIMEN COLLECTION AND STORAGE

Serum samples or heparinized plasma samples are recommended.

The required sample volume for duplicate measurement is 2 x 50 µL.

Specimens may be stored refrigerated (2 °C - 8 °C) up to 24 h.

For longer storage the specimen should be kept frozen (≤-18 °C). There is no indication of loss of activity if samples are freeze-thawed but it is recommended to avoid repeated thawing and freezing.

Do not use samples that are highly haemolysed.

5 PRECAUTIONS

- TPS (Cyk18) ELISA is for *in vitro* diagnostic use only.
- Do not use the kit after expiry date.
- Do not mix reagents from different lots.
- The accuracy of the test is related to adherence to the assay procedure and accurate volume pipetting.
- Standards, controls and samples in duplicates are recommended.
- All patient specimens should be regarded as contagious and must be handled and disposed of according to appropriate regulations.
- Wear protective gloves and eyewear.
- Avoid microbiological contamination of reagents.
- A Safety Data Sheet (SDS) is available on request.

6 MATERIALS REQUIRED BUT NOT PROVIDED

- Routine laboratory equipment, e.g. precision pipettes and vortex
- Shaker for incubation with a recommended oscillation ~450 rpm
- Deionized or distilled water
- Wash equipment for microplates
- Microplate reader (wavelength 450 nm)

7 KIT CONTENT

Materials supplied for 96 determinations.

TPS (Cyk18) ELISA Microstrips

1 plate, 96 dry wells (12 strips), coated with monoclonal antibodies.
Packed in aluminium bag with desiccating device. Ready for use.

TPS (Cyk18) ELISA HRP Conjugate

1 vial, 11 mL, monoclonal antibodies conjugated with HRP in protein stabilized buffer, pH 7.5.
Colour blue. Preservative added. Ready for use.

TPS (Cyk18) ELISA Diluent (Standard 0 U/L)

1 vial, 5 mL, sample diluent and standard 0 U/L, protein stabilized buffer, pH 7.5.
Colour yellow. Preservative added. Ready for use.

TPS (Cyk18) ELISA Standard (30, 150, 500, 1200 U/L)

4 vials, 1 mL/vial, TPS (Cyk18) ELISA standard material in protein stabilized buffer, pH 7.5.
Concentrations as stated on vials.
Colour yellow. Preservative added. Ready for use.

TPS (Cyk18) ELISA Control (Low, High)

2 vials, 1 mL/vial, TPS (Cyk18) ELISA standard material in protein stabilized buffer, pH 7.5.
Colour yellow. Preservative added. Ready for use.

Wash tablet

One blisterpacked tablet, the tablet should be dissolved in 500 mL of fresh deionized water.

TMB Substrate

1 vial, 22 mL. Protect from light and keep lid tightly closed. Ready for use.
Do not use more than required for analysis.

Stop Solution

1 vial, 12 mL, 0.5 M sulfuric acid. Ready for use.

Sealing Tape

1 sheet Sealing Tape for Microstrips.

TPS (Cyk18) ELISA Certificate

Certificate of lot content and release data.

TPS (Cyk18) ELISA IFU

Instruction for use

8 KIT AND REAGENT STORAGE

The kit should be stored at 2 °C - 8 °C. Do not freeze.

Reagents should always be stored in their original containers. If not all strips are used at once, be sure to store remaining TPS (Cyk18) ELISA Coated Microstrips in the aluminium bag with the desiccant device.

The wash solution is stable for 4 weeks when stored at 2 °C - 8 °C.

9 PREPARATION OF REAGENTS

Preparation of the Wash Solution:

Dissolve one wash tablet in 500 mL of deionised or distilled water.

10 ASSAY PROCEDURE


The assay should be performed at room temperature, 18 °C - 26 °C

Allow all reagents and samples to adjust to room temperature. Mix all reagents on a vortex prior to use.

Step 1 and 2 should be performed sequentially without interruption.

1. Pipette 50 µL standards, controls or samples per defined well. Leave two empty wells for background absorbance measurement (optional). It is recommended to run duplicates of standards, controls and samples.
2. Add 100 µL TPS (Cyk18) ELISA HRP Conjugate per well except the blank wells.
3. Incubate for 2 h ± 10 min on a shaker at ~450 rpm. Setting of the shaker at 450 rpm is important for correct result.
4. Aspirate and wash the wells 3 times with 0.3 mL wash solution. Make sure there is no residual wash solution in the wells after washing.
5. Add 200 µL TMB substrate per well, including the blank wells.
6. Incubate in darkness for 20 ± 1 min.
7. Add 50 µL Stop Solution per well. Agitate on a shaker for 1 min.
8. Read the absorbance at 450 nm within 30 min after addition of the Stop Solution
9. Calculate the cytokeratin 18 concentration (U/L) of the samples.
Samples showing concentrations > 1200 U/L should be suitably diluted with TPS (Cyk18) ELISA Diluent (Standard 0 U/L) before repeated analysis.

Schematic assay procedure

11 CALCULATION OF RESULTS

Use computer software for process the raw data. Recommended curve fitting algorithm is spline smoothed. For quality control and for generation of valid data, ensure that included controls are within range.

NOTE: Different curve fitting algorithms may give slightly different patient sample concentrations; therefore it is important to always use the same algorithm when comparing results

12 PERFORMANCE CHARACTERISTICS

Measuring range:

The measuring range is 10 - 1200 U/L.

Normal range:

The 95th percentile for apparently healthy Swedish blood donors has been determined to 80 U/L. It is recommended that each laboratory establish its own normal range.

Detection limit:

The minimal detectable concentration in TPS (Cyk18) ELISA is < 6 U/L, defined as the concentration of TPS antigen that corresponds to the OD-value being two standard deviations from the OD-value of standard 0 U/L.

Interference:

Interference from added hemoglobin caused up to three times higher values than expected, 300% recovery. There was no detectable interference from bilirubin and lipids, recovery ranges from 97 - 105% and 97 - 107% respectively.

Precision:

The intra- and inter-assay precision of the assay, defined according to CLSI guidelines, ranges from 1 - 8% CV. The average within and between assays CV was 5% and 4% respectively.

Recovery:

Determined recovery was 95 - 102% after adding specified quantities of TPS antigen to human serum specimens.

Linearity:

Determined recovery was 95 - 113% after diluting high concentration samples with TPS (Cyk18) ELISA Diluent.

High-dose hook effect:

The assay does not show any high-dose hook effect up to 20 000 U/L.

13 LIMITATIONS OF THE ASSAY

TPS (Cyk18) ELISA is not intended to be used as a screening test.

TPS (Cyk18) ELISA values should be interpreted in conjunction with all available clinical information.

Elevated results may be found e.g. in cases of pregnancy, liver disease, renal failure and general infections. If a temporary infection is suspected, it may be necessary to repeat the test at a later occasion.

14 WARRANTY

Any change or modification in the procedure, not recommended by DRG, may affect the results. In such event, DRG disclaims all warranties expressed, implied or statutory, including the implied warranty of merchantability and the fitness for use.

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14 REFERENCES / LITERATUR / BIBLIOGRAFIA / BIBLIOGRAFÍA

1. Rydlander L *et al.* Molecular characterization of a Tissue Polypeptide Specific-antigen epitope and its relationship to human cytokeratin 18.
Eur J Biochem 1996; 241: 309-314.
2. Stigbrand T *et al.* Epitope specificity of 30 monoclonal antibodies against cytokeratin antigens: The ISOBM TD5-1 Workshop.
Tumor Biol 1998; 19:132-152.
3. van Dalen A *et al.* The prognostic significance of increasing marker levels in metastatic breast cancer patients with clinically complete remission, partial remission or stable disease.
Int J Biol Markers 1998; 13: 10-15.
4. Kil PJ *et al.* Tissue Polypeptide Specific antigen (TPS) determinations before and during intermittent Maximal Androgen Blockade in patients with metastatic prostatic carcinoma.
Eur Urol 2003; 43: 31-38.
5. Barak V *et al.* Clinical utility of cytokeratins as tumor markers.
Clin Biochem 2004; 37: 529-540.
6. van Dalen *et al.* Significance of serum CA125 and TPS antigen levels for determination of overall survival after three chemotherapy courses in ovarian cancer patients during long-term follow-up.
Eur. J. Gynaec. Oncol. 2009; 30(6): 609-15.
7. TPS; Xie S *et al.* Serum tissue polypeptide-specific antigen is an independent prediktor in breast cancer.
Acta Histochem 2014; 116: 372-376.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico In vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	Numero di Catalogo	Número de catálogo	Référence de catalogue
	Batch code *	Chargencode *	Codice del lotto	Codigo de lote	Numéro de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservacion	Température de conservation
	Use-by date *	Verwendbar bis *	Utilizzare prima del	Establa hasta	Utiliser jusque
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Caution *	Achtung *			
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité