## Instructions for Use

# $\alpha$ 1-Antitrypsin ELISA

IVD

REF EIA-5299

$\Sigma$  96

**DRG**

**DRG Instruments GmbH**, Germany  
Frauenbergstraße. 18, D-35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

Distributed by:

**DRG**

**DRG International, Inc.**, USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.  
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.  
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.  
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.**

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## 1 INTENDED USE

The assay is an enzyme immunoassay intended for the quantitative determination of α1-Antitrypsin in stool. It is for *in vitro* diagnostic use only.

## 2 INTRODUCTION

Intestinal protein loss is a serious consequence of various systemic or local gastrointestinal diseases (e. g. allergies, chronic inflammation, malignancies). These pathologies damage the mucosal integrity and/or cause lymphostasis, thereby leading to an increased transfer of plasma proteins into the bowel lumen. Subsequently, hypoproteinemia accompanied with edema may develop. This condition is diagnosed by exclusion of other sources of protein loss and by proof of an elevated α1-antitrypsin concentration in stool.

In serum, α1-antitrypsin represents the majority of serine protease inhibitors and protects tissues from protease damages during inflammation. The protein is synthesized primarily in the liver but also to a small extent in intestinal macrophages, monocytes, and intestinal epithelial cells. Since α1-antitrypsin is relatively resistant against enzymatic digestion, the secreted amount in stool reflects the internal concentration of the protein. An elevated α1-antitrypsin stool concentration is therefore a widely recognized marker for intestinal protein loss and for an increased mucosal permeability.

In clinical routine, the α1-antitrypsin clearance (ratio of the α1-antitrypsin ELISA values of stool and serum samples) has been established along with the sole determination of the 24h α1-antitrypsin secretion in stool. Thus the group of J. S. Fordtran reports that the sole determination of the α1-antitrypsin concentration in stool yielded false positive or false negative results in 21 % of the patients compared to the α1-antitrypsin clearance measurement (Strygler et al. 1990).

The analytical quality of DRG's α1-antitrypsin ELISA surpasses by far the conventional radial immunodiffusion (RID) technique in the determination of serum, stool and tissue culture supernatants. In direct comparison, the concentrations measured with the ELISA were approximately 30 % above the corresponding RID levels. Cell culture supernatants of an intestinal cell line as well as fecal samples of lymphostasis patients yielded negative results with RID. Our ELISA could detect α1-antitrypsin in all of these samples, in some of them even in very high concentrations.

These results clearly prove that the α1-antitrypsin ELISA is far more sensitive than the conventional method and that it recognizes not only hepatic but also enteral α1-antitrypsin. The discrepancy of both methods and hence the superiority of the ELISA to RID is especially striking in the analysis of extremely high enteral protein losses. The combination of two specific antibodies in our α1-antitrypsin ELISA widely excludes the possibility of false negative results thereby enabling a reliable diagnostics of enteral protein loss.

### Indications

- Suspected enteric protein loss
- Crohn's disease
- Necrotic enterocolitis
- Chronic mesenteric ischemia
- Viral, bacterial, allergic, or autoimmune-induced gastrointestinal inflammation

## 3 MATERIAL SUPPLIED

Label	Kit Components	Quantity
PLATE	Microtiter plate, pre-coated	12 x 8 wells
WASHBUF	Wash buffer concentrate, 10 x	2 x 100 mL
CONJ	Conjugate concentrate, (goat-anti-α1-antitrypsin, peroxidase-labelled)	1x 200 µL
STD	Standards*, lyophilized (0; 3.3; 10; 30; 90 µg/L)	2 x 5 vials
CTRL 1	Control, lyophilized (see specification for range)	2 x 1 vial
CTRL 2	Control, lyophilized (see specification for range)	2 x 1 vial
SUB	Substrate (tetramethylbenzidine), ready-to-use	1 x 15 mL
STOP	Stop solution, ready-to-use	1 x 15 mL
IDK Extract®	Extraction buffer concentrate IDK Extract® 2.5x	2 x 100 mL

\* The used standards have been calibrated on the WHO reference material CRM 470.

#### 4 MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultra pure water\*
- Stool sample application system such as cat. no.: EIA-5674
- Calibrated precision pipettors and 10 - 1000 µL single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeating pipets
- Vortex-Mixer
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

\* DRG recommends the use of Ultra Pure Water (Water Type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≤18.2 MΩ cm).

#### 5 STORAGE AND PREPARATION OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µL** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:**  
The **wash buffer concentrate (WASHBUF)** has to be diluted with ultra pure water **1:10** before use (100 mL WASHBUF + 900 mL ultra pure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at room temperature or in a water bath at 37 °C. The **WASHBUF** is stable at **2 °C - 8 °C** until the expiry date stated on the label.  
**Wash buffer** (1:10 diluted WASHBUF) can be stored in a closed flask at **2 °C - 8 °C for 1 month.**
- **Preparation of the extraction buffer:**  
The **extraction buffer concentrate IDK Extract®** has to be diluted with ultra pure water **1:2.5** before use (100 mL *IDK Extract®* + 150 mL ultra pure water), mix well. Crystals could occur due to high salt concentration in the concentrate. Before dilution, the crystals have to be redissolved at 37 °C in a water bath. The **IDK Extract®** is stable at **2 °C - 8 °C** until the expiry date stated on the label.  
Extraction buffer (1:2.5 diluted *IDK Extract®*) can be stored in a closed flask at **2 °C - 8 °C for 4 months.**
- The **lyophilized standards (STD)** and **controls (CTRL)** are stable at 2 °C - 8 °C until the expiry date stated on the label. Before use, STD and CTRL have to be reconstituted with **500 µL of ultra pure water and mixed by gentle inversion to ensure complete reconstitution.** Allow the vial content to dissolve for 10 minutes and then mix thoroughly. **Standards and controls** (reconstituted STD and CTRL) **can be stored at 2 °C - 8 °C for 4 weeks.**
- **Preparation of the conjugate:**  
Before use, the **conjugate concentrate (CON J)** has to be diluted **1:101** in wash buffer (100 µL CONJ + 10 mL wash buffer). The CONJ is stable at **2 °C - 8 °C** until the expiry date stated on the label.  
**Conjugate** (1:101 diluted CONJ) **is not stable and cannot be stored.**
- All other test reagents are ready to use. The test reagents are stable until the expiry date (see label of test package) when stored at **2 °C - 8 °C.**

## 6 STORAGE AND PREPARATION OF SAMPLES

### 6.1 Storage

The **sample stability** is as follows:

**Raw stool:** 3 days at room temperature (15 °C - 30 °C),  
3 days at 2 °C - 8 °C,  
or at least 4 weeks at -20 °C

**Stool extracts:** 9 days at room temperature, 2 °C - 8 °C or -20 °C, maximum 3 freeze-thaw cycles

### 6.2 Extraction of the stool sample

**Extraction buffer** (1:2.5 diluted IDK Extract<sup>®</sup>) is used as a sample extraction buffer.

We recommend the following sample preparation:

#### Stool Sample Application System (SAS) (Cat. No.: EIA-5674)

##### Stool sample tube – Instruction for use

Please note that the dilution factor of the final stool suspension depends on the amount of stool sample used and the volume of the buffer.

##### SAS with 1.5 mL sample extraction buffer:

Applied amount of stool: 15 mg  
Buffer Volume: 1.5 mL  
Dilution Factor: 1:100

Please follow the instructions for the preparation of stool samples using the SAS as follows:

- The raw stool sample has to be thawed. For particularly heterogeneous samples we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- Fill the empty stool sample tube with 1.5 mL sample extraction buffer** (1:2.5 diluted IDK Extract<sup>®</sup>) before using it with the sample.  
Important: Allow the sample extraction buffer to reach room temperature.
- Unscrew the tube (yellow part of cap) to open. Insert yellow dipstick into the sample. The lower part of the dipstick has notches which need to be covered completely with stool after inserting it into the sample. Place dipstick back into the tube. When putting the stick back into the tube, excess material will be stripped off, leaving 15 mg of sample to be diluted. Screw tightly to close the tube.
- Shake the tube well until no stool sample remains in the notches.  
Important: Please make sure that you have a maximally homogenous suspension after shaking. Especially with more solid samples, soaking the sample in the tube with sample extraction buffer for ~ 10 minutes improves the result
- Allow sample to stand for ~10 minutes until sediment has settled. Floating material like shells of grains can be neglected.
- Carefully unscrew the complete cap of the tube including the blue ring plus the dipstick. Discard cap and dipstick. Make sure that the sediment will not be dispersed again.

**Dilution I: 1:100**

### 6.3 Dilution of samples

The supernatant of the sample preparation procedure (dilution I) is further diluted **1:250** in wash buffer.

For example:

**20 µL** supernatant (dilution I) + **980 µL** wash buffer, mix well = **dilution II** (1:50)

**200 µL** dilution II + **800 µL** wash buffer, mix well = **dilution III** (1:5)

This results in a **final dilution of 1:25000**.

For analysis, pipette **100 µL** of **dilution III** per well.

## 7 ASSAY PROCEDURE

### 7.1 Principle of the test

The assay utilizes the sandwich technique with two selected polyclonal antibodies that bind to human α1-antitrypsin. Standards, controls and prediluted samples which are assayed for human α1-antitrypsin are added into the wells of a microplate coated with a high affine polyclonal anti-human α1-antitrypsin antibody. During the first incubation step, α1-antitrypsin is bound by the immobilized antibody. Then a peroxidase-conjugated polyclonal anti-human α1-antitrypsin antibody is added into each microtiter well and a sandwich of capture antibody – human α1-antitrypsin – peroxidase-conjugate is formed. Tetramethylbenzidine is used as peroxidase substrate. Finally, an acidic stop solution is added to terminate the reaction. The colour changes from blue to yellow.

The intensity of the yellow colour is directly proportional to the concentration of α1-antitrypsin. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. α1-Antitrypsin present in the samples, is determined directly from this curve.

### 7.2 Test procedure

Bring all **reagents and samples to room temperature** (15 °C - 30 °C) and mix well.

Mark the positions of standards / controls / samples on a protocol sheet.

Take as many microtiter strips as needed from kit. Store unused strips covered at 2 °C - 8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or DRG.

We recommend to carry out the tests in duplicate.

1. Before use, wash the wells **5 times** with **250 µL wash buffer** into each well. After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
2. Add each **100 µL standards / controls / samples** into the respective wells.
3. Cover the plate tightly and incubate for **1 hour at room temperature** (15 °C - 30 °C) on a **horizontal shaker\***.
4. Discard the contents of each well and wash **5 times** with **250 µL of wash buffer**. After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
5. Add **100 µL conjugate** (diluted CONJ) into each well.
6. Cover the strips and incubate for **1 hour at room temperature** (15 °C - 30 °C) on a **horizontal shaker\***.
7. Discard the contents of each well and wash **5 times** with **250 µL wash buffer**. After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
8. Add **100 µL substrate** (SUB) into each well.
9. Incubate for **10 - 20 minutes\*\* at room temperature** (15 °C - 30 °C) in the **dark**.
10. Add **100 µL stop solution** (STOP) into each well and mix well.
11. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at **405 nm** against 620 nm as a reference.

\* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

\*\* The intensity of the colour change is temperature sensitive. We recommend observing the colour change and stopping the reaction upon good differentiation.

## 8 RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-parameter-algorithm".

1. 4 Parameter algorithm  
It is recommended to use a linear ordinate for the optical density and a logarithmic abscissa for the concentration. When using a logarithmic abscissa, the zero standard must be specified with a value less than 1 (e. g. 0.001).
2. Point-to-point-calculation  
We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.
3. Spline algorithm  
We recommend a linear ordinate for the optical density and a linear abscissa for the concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the program used, the duplicate values should be evaluated manually.

### Stool samples

The obtained results have to be multiplied with the **dilution factor of 25 000** to get the actual concentration. In case **another dilution factor** has been used, multiply the obtained result by the dilution factor used.

## 9 LIMITATIONS

Samples with concentrations above the measurement range (see definition below) can be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

$$\text{highest concentration of the standard curve} \times \text{sample dilution factor to be used}$$

The lower limit of the measurement range can be calculated as:

$$\text{Analytical sensitivity} \times \text{sample dilution factor to be used}$$

Analytical sensitivity see chapter "Performance Characteristics".

## 10 QUALITY CONTROL

DRG recommends the use of external controls for internal quality control, if possible.

Control samples should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

### 10.1 Reference range

Based on studies of stool samples of apparently healthy persons (n = 76), a **cut off value of < 26.8 mg/dL** was estimated.

We recommend each laboratory to establish its own reference range.

## 11 PERFORMANCE CHARACTERISTICS

### 11.1 Analytical Sensitivity

The detection limit was set as  $B_0 + 2 \text{ SD}$ . The Zero-standard was measured 20 times.

The values were estimated in relation to the concentration of the calibration curve and resulted in a detection limit of 0.72 µg/L without consideration of the sample dilution factor.

### 11.2 Precision and reproducibility

#### Intra-Assay (n = 20)

The precision (intra-assay variation) of the α1-Antitrypsin ELISA was calculated from 20 replicate determinations on each one of two samples.

Sample	α1-antitrypsin [mg/dL]	CV [%]
1	15.2	4.5
2	42.2	13.1

#### Inter-Assay (n = 20)

The total precision (inter-assay variation) of the α1-Antitrypsin ELISA was calculated from data on 2 samples obtained in 20 different assays by three technicians on two different lots of reagents over a period of three months.

Sample	α1-antitrypsin [mg/dL]	CV [%]
1	16.15	9.8
2	54.46	14.8

### 11.3 Spiking Recovery

Two samples were spiked with different α1-antitrypsin standard amounts and measured with the assay.

Sample	Unspiked Sample [µg/L]	Spike [µg/L]	α1-antitrypsin expected [µg/L]	α1-antitrypsin measured [µg/L]
A	6.3	1.65	7.95	7.2
		5.0	11.3	11.4
		15.0	21.3	20.9
		45.0	51.3	46.7
B	5.6	1.65	7.3	7.0
		5.0	10.6	10.9
		15.0	20.6	17.5
		45.0	50.6	46.7

### 11.4 Specificity

No cross reactivity with other plasma proteins in stool was observed.

No cross reactivity with α1-antitrypsin in mouse serum was observed.

### 11.5 Dilution Recovery

Two patient serum samples were diluted with wash buffer. The results are shown below (n = 2):

Sample	Dilution	α1-antitrypsin expected [mg/dL]	α1-antitrypsin measured [mg/dL]
A	1:12 500	48.0	48.88
	1:25 000	24.5	23.25
	1:50 000	12.3	11.9
	1:100 000	6.1	6.0
B	1:12 500	158.4	158.4
	1:25 000	79.3	99.0
	1:50 000	39.6	33.0
	1:100 000	19.8	22.1

## 12 PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Proclin as bactericides. Sodium azide and Proclin are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

## 13 TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analyzed with each run.
- Reagents should not be used beyond the expiration date stated on kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

## 14 GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. DRG can therefore not be held responsible for any damage resulting from wrong use.

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## 15 REFERENCES/LITERATURE / LITERATUR

1. Amarri, S. et al., 2006. Changes of gut microbiota and immune markers during the complementary feeding period in healthy breast-fed infants. *Journal of pediatric gastroenterology and nutrition*, 42(5), pp.488–95.
2. Faust, D et al., 2001. Determination of alpha1-proteinase inhibitor by a new enzyme linked immunosorbant assay in feces, serum and an enterocyte-like cell line. *Zeitschrift für Gastroenterologie*, 39(9), pp.769–74.
3. Faust, D. et al., 2002. Regulation of alpha1-proteinase inhibitor release by proinflammatory cytokines in human intestinal epithelial cells. *Clinical and experimental immunology*, 128(2), pp.279–84.
4. Hsu, P.-I. et al., 2010. Diagnosis of gastric malignancy using gastric juice alpha1-antitrypsin. *Cancer epidemiology, biomarkers & prevention*, 19(2), pp.405–11.
5. Lamprecht, M. et al., 2012. Probiotic supplementation affects markers of intestinal barrier, oxidation, and inflammation in trained men; a randomized, double-blinded, placebo-controlled trial. *Journal of the International Society of Sports Nutrition*, 9(1), p.45.
6. Muss, C., Schütz, B. & Kirkamm, R., 2002. Alpha1-Antitrypsin - ein objektiver Verlaufsparemeter bei entzündlichen Darmerkrankungen. *Ärztezeitschrift für Naturheilverfahren*, 43(4).
7. Oswari, H. et al., 2013. Comparison of stool microbiota compositions, stool alpha1-antitrypsin and calprotectin concentrations, and diarrhoeal morbidity of Indonesian infants fed breast milk or probiotic/prebiotic-supplemented formula. *Journal of paediatrics and child health*, Epub ahead of print.
8. Quint, J.K. et al., 2011. SERPINA1 11478G>A variant, serum α1-antitrypsin, exacerbation frequency and FEV1 decline in COPD. *Thorax*, 66(5), pp.418–24.
9. Ragab, H.M. et al., 2009. Clinical utility of serum TNF alpha and alpha-1 anti-trypsin in predicting the stage and progression of lung cancer. *International Journal of Integrative Biology*, 7(1), pp.45–52.
10. Roedel, N. et al., 2009. High frequency of LMAN1 abnormalities in colorectal tumors with microsatellite instability. *Cancer research*, 69(1), pp.292–9.
11. Strygler, B. et al., 1990. α1-Antitrypsin Excretion in Stool in Normal Subjects and in Patients With Gastrointestinal Disorders. *Gastroenterology*, 99(5), pp.1380–1387.
12. Török, E. et al., 2011. Primary human hepatocytes on biodegradable poly(l-lactic acid) matrices: a promising model for improving transplantation efficiency with tissue engineering. *Liver transplantation*, 17(2), pp.104–14.

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico In vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	Numero di Catalogo	Número de catálogo	Référence de catalogue
	Batch code *	Chargencode *	Codice del lotto	Codigo de lote	Numéro de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis *	Utilizzare prima del	Establa hasta	Utiliser jusque
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Caution *	Achtung *			
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité