Instructions for Use

ox-LDL/MDA-Adducts ELISA

IVD

REF EIA-5656

 96

DRG 

DRG Instruments GmbH, Germany
Frauenbergstraße. 18, D-35039 Marburg
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50
Website: www.drg-diagnostics.de
E-mail: drg@drg-diagnostics.de

Distributed by:

DRG 

DRG International, Inc., USA
841 Mountain Ave., Springfield, NJ 07081
Phone: (973) 564-7555, Fax: (973) 564-7556
Website: www.drg-international.com
E-mail: corp@drg-international.com

**Please use only the valid version of the Instructions for Use provided with the kit.
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.**

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1 INTENDED USE

This assay is an enzyme immunoassay intended for the quantitative determination of ox-LDL in EDTA-plasma, serum and dried blood spots.

The test recognizes MDA-modified apolipoprotein B 100, even if it contains less than 60 MDA units per molecule. It is for in vitro diagnostic use only.

2 INTRODUCTION

Lipid peroxidation is a natural process essential for cell growth. However, when the oxidative stress overwhelms the antioxidative cell defense, the balance is disturbed and enhanced formation of lipid peroxidation products occurs. At present, lipid peroxidation is considered to be one of the basic mechanisms involved in the initiation and progression of many diseases. Various studies have provided evidence that oxidative stress resulting in lipid peroxidation and protein modification is involved in the pathogenesis of atherosclerosis and coronary heart disease.

Lipid peroxidation products are formed during normal cell metabolism via producing an excess of free radicals that can react with unsaturated fatty acids, in particular low-density lipoprotein (LDL), the major carrier of plasma cholesterol. LDL is eliminated by macrophages. Normally, receptor-mediated uptake of LDL is suppressed through down-regulation of LDL receptor expression in response to increasing cholesterol levels. Once LDL is oxidised, it is still internalised by macrophages but through scavenger receptors whose expression is not controlled by cholesterol loading. The binding of oxidised LDL (ox-LDL) is the step by which cholesterol accumulation in macrophages is induced transforming them into lipid-loaded foam cells. This process is accompanied by extensive cell proliferation and elaboration of extra cellular matrix components and contributes to the genesis and progression of atherosclerosis by promoting endothelial damage and amplifying the inflammatory response within the vessel wall. Cholesterol-loaded macrophage foam cells are present in the earliest detectable atherosclerotic lesions, the precursor of more complex atherosclerosis that cause stenosis and limited blood flow. These advanced lesions ultimately represent the sites of thrombosis leading to myocardial infarction.

3 MATERIAL SUPPLIED

Label	Kit components	Quantity
PLATE	Microtiter plate, pre-coated	12 x 8 wells
WASHBUF	Wash buffer concentrate, 10 x	2 x 100 mL
CONJ	Conjugate concentrate, goat-anti ox-LDL, peroxidase-labelled	1 x 150 µL
CONJBUF	Conjugate dilution buffer, ready-to-use	1 x 15 mL
STD	Standards, lyophilised (see specification for concentrations)	4 x 5 vials
CTRL1	Control, lyophilised (see specification for range)	4 x 1 vial
CTRL2	Control, lyophilised (see specification for range)	4 x 1 vial
SAMPLEBUF	Sample dilution buffer, ready-to-use	1 x 30 mL
SUB	Substrate (tetramethylbenzidine), ready-to-use	1 x 15 mL
STOP	Stop solution, ready-to-use	1 x 15 mL

4 MATERIAL REQUIRED BUT NOT SUPPLIED

- Ultrapure water*
- Dried blood spot carrier such as DrySpot cat. no.: DZ9020 or DZ902D
- Calibrated precision pipettors and 10 µL - 1000 µL single-use tips
- Foil to cover the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel pipets or repeater pipets
- Centrifuge 3000 g
- Vortex Mixer
- Standard single-use laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader (required filters see chapter 7)

* DRG recommends the use of ultrapure water (water type 1; ISO 3696), which is free of undissolved and colloidal ions and organic molecules (free of particles > 0.2 µm) with an electrical conductivity of 0.055 µS/cm at 25 °C (≥ 18.2 MΩ cm).

5 STORAGE AND PREPARATION OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at the conditions stated on the label. **Prepare only the appropriate amount necessary for each run.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µL** should be centrifuged before use to avoid loss of volume.
- **Preparation of the wash buffer:**
The **wash buffer concentrate (WASHBUF)** has to be diluted with ultrapure water **1:10** before use (100 mL WASHBUF + 900 mL ultrapure water), mix well.
Crystals could occur due to high salt concentration in the concentrate. The crystals must be redissolved at room temperature or in a water bath at 37 °C before dilution of the buffer solution.
The **WASHBUF** is stable at **2 °C - 8 °C** until the expiry date stated on the label.
Wash buffer (1:10 diluted WASHBUF) can be stored in a closed flask at **2 °C - 8 °C for 1 month**.
- The **lyophilised standards (STD)** and **controls (CTRL)** are stable at **2 °C - 8 °C** until the expiry date stated on the label.
Before use, STD and CTRL have to be reconstituted with **500 µL of ultrapure water** and mixed by gentle inversion to ensure complete reconstitution. Allow the vial content to dissolve for 10 minutes and then mix thoroughly.
Standards and controls (reconstituted STD and CTRL) **are not stable and cannot be stored**.
- **Preparation of the conjugate:**
Before use, the **conjugate concentrate (CONJ)** has to be diluted **1:101** in conjugate dilution buffer (100 µL CONJ + 10 mL CONJBUF).
The CONJ is stable at **2 °C - 8 °C** until the expiry date stated on the label.
Conjugate (1:101 diluted CONJ) **is not stable and cannot be stored**.
- All other test reagents are ready to use. Test reagents are stable until the expiry date (see label) when stored at **2 °C - 8 °C**.

6 STORAGE AND PREPARATION OF SAMPLES

6.1 EDTA plasma and Serum

6.1.1 Sample storage

Venous fasting blood is suited for this test system. Samples should be stored at -20 °C up to the measurement. The maximum storage time at -20 °C is two years.

Lipemic or hemolytic samples may give erroneous results and should not be used for analysis.

Samples with visible amounts of precipitates should be centrifuged (5 min at 3 000 g) prior to measurement and the resulting supernatant used in the test.

6.1.2 Sample preparation

EDTA plasma or serum samples must be diluted **1:10** before performing the assay, e.g. **30 µL** sample + **270 µL** sample dilution buffer (SAMPLEBUF), mix well.

100 µL of the dilution are used in the test per well.

6.2 Dried blood spots

6.2.1 Collection and storage of dried blood spots

50 µL whole blood dripped on a dried sample carrier cleared by DRG are suitable as sample material after complete drying. We recommend DrySpotD (catalogue no DZ9020 or DZ9021) as dried blood spot carrier.

The moistened cards are stable for 3 weeks at room temperature.

6.2.2 Preparation of the dried blood samples

1. Label 1.5-mL polypropylene tubes.
2. Remove filter from sampling device (always wear gloves since the sample is potentially infectious).
3. Put filter in a labelled tube.
4. Add **300 µL** sample dilution buffer (SAMPLEBUF) to each sample, allow sample to stand for **20 min** at room temperature (15 °C - 30 °C).
5. Vortex for **10 seconds**. The filter will decolourise.
6. Centrifuge the samples for **5 min** at **3000 g** to remove residual filter pieces.

For testing in duplicates, pipette 2 x 100 µL of each prepared sample per assay run.

7 ASSAY PROCEDURE

7.1 Principle of the test

This ELISA is designed for the quantitative determination of ox-LDL/MDA adducts.

This assay is a sandwich ELISA for the direct measurement of ox-LDL in human EDTA plasma, serum and dried blood spots.

Standards, controls and samples containing human ox-LDL are added to wells of microplate coated with high affinity antibodies. During the first incubation period, the antibodies immobilised on the wall of the microtiter wells capture the antigen in the patient samples. After washing away the unbound components from samples, a peroxidase-conjugated antibody is added to each microtiter well. Tetramethylbenzidine (TMB) is used as a substrate for peroxidase. Finally, an acidic stop solution is added to terminate the reaction. The intensity of the yellow colour is directly proportional to the ox-LDL concentration of sample. A dose response curve of the absorbance unit (optical density, OD at 450 nm) vs. concentration is generated, using the values obtained from the standard. Ox-LDL, present in the patient samples, is determined directly from this curve.

7.2 Test procedure

Bring all **reagents and samples to room temperature** (15 °C - 30 °C) and mix well.

Mark the positions of standards/controls/samples on a protocol sheet.

Take as many microtiter strips as needed from the kit. Store unused strips together with the desiccant bag in the closed aluminium packaging at 2 °C - 8 °C. Strips are stable until expiry date stated on the label.

For automated ELISA processors, the given protocol may need to be adjusted according to the specific features of the respective automated platform. For further details please contact your supplier or DRG.

We recommend to carry out the tests in duplicate.

1. **Before use**, wash the wells **5 times** with **250 µL wash buffer**.
After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
2. Add each **100 µL standards/controls/prepared samples** into the respective wells.
3. Cover the strips and incubate for **1 hour** at room temperature (15 °C - 30 °C) on a **horizontal shaker***.
4. Discard the content of each well and wash **5 times** with **250 µL wash buffer**.
After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
5. Add **100 µL conjugate** (diluted CONJ) into each well.
6. Cover the strips and incubate for **1 hour** at room temperature (15 °C - 30 °C) on a **horizontal shaker***.
7. Discard the content of each well and wash **5 times** with **250 µL wash buffer**.
After the final washing step, remove residual wash buffer by firmly tapping the plate on absorbent paper.
8. Add **100 µL substrate** (SUB) into each well.
9. Incubate for **10 - 20 minutes**** at room temperature (15 °C - 30 °C) in the dark.
10. Add **100 µL stop solution** (STOP) into each well and mix well.
11. Determine **absorption immediately** with an ELISA reader at **450 nm** against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at **405 nm** against 620 nm as a reference.

* We recommend shaking the strips at 550 rpm with an orbit of 2 mm.

** The intensity of the color change is temperature sensitive. We recommend observing the color change and stopping the reaction upon good differentiation.

8 RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend using the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

The plausibility of the duplicate values should be examined before the automatic evaluation of the results. If this option is not available with the programme used, the duplicate values should be evaluated manually.

8.1 EDTA-plasma, serum and dried blood spots

The obtained results have to be multiplied by the **dilution factor of 10** to get the actual concentrations.

In case **another dilution factor** has been used, multiply the obtained result by the dilution factor used.

9 LIMITATIONS

Samples with concentrations above the measurement range (see definition below) can be further diluted and re-assayed. Please consider this greater dilution when calculating the results.

Samples with concentrations lower than the measurement range (see definition below) cannot be clearly quantified.

The upper limit of the measurement range can be calculated as:

$$\text{highest concentration of the standard curve} \times \text{sample dilution factor to be used}$$

The lower limit of the measurement range can be calculated as:

$$\text{Analytical sensitivity} \times \text{sample dilution factor to be used}$$

Analytical sensitivity see chapter "Performance Characteristics".

10 QUALITY CONTROL

DRG recommends the use of external controls for internal quality control, if possible.

Control samples or serum pools should be analysed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid if within the same assay one or more values of the quality control sample are outside the acceptable limits.

10.1 Reference range

Based on DRG studies of evidently healthy persons (n=77; Germany) a mean value of 287 ng/mL was estimated.

Serum/Plasma (n = 77; Germany): 287 (41 - 2261) ng/mL

Median: 141 ng/mL

We recommend each laboratory to establish its own reference range. The values mentioned above are only for orientation and can deviate from other published data.

10.2 Additional reference ranges

Within a scientific study a mean value of 95.32 ± 37.85 ng ox-LDL/mL was estimated for control subjects (healthy, n=120; Tunisia) using EIA-5656.

Serum/Plasma (controls, n=120; Tunisia): 95.32 ± 37.85 ng/mL*

Furthermore, the obtained results demonstrate that a significantly elevated ox-LDL concentration (142.37 ± 49.84 ng ox-LDL/mL) was found in type 2 diabetes patients (n=86) compared with healthy controls.

In addition, higher ox-LDL values were detected in type 2 diabetes patients with hypertension, as compared with diabetic patients without hypertension.

The results of the study are summarized in the following table.

Sample	Ox-LDL [ng/mL]
Controls, healthy (n=120)	95.32 ± 37.85
type 2 diabetes patients (n=86)	142.37 ± 49.84
type 2 diabetes patients without hypertension	111.16 ± 33.42
type 2 diabetes patients with hypertension	157.4 ± 49.9

*Koubaa N et al. (2007) *Clin. Biochem.* 40, 1007-1014

11 PERFORMANCE CHARACTERISTICS

11.1 Analytical sensitivity

The following value has been estimated based on the concentrations of the standard without considering possibly used sample dilution factors.

Limit of blank, LoB 4.130 ng/mL

11.2 Accuracy – Precision

Repeatability (Intra-Assay); n = 42

The repeatability was assessed with 2 plasma samples under **constant** parameters (same operator, measurement system, day and kit lot). The results below were obtained without consideration of the sample dilution factor.

Sample	Mean value [ng/mL]	CV [%]
1	30.65	3.9
2	40.33	5.7

Reproducibility (Inter-Assay); n = 15

The reproducibility was assessed with 2 plasma samples under **varying** parameters (different operators, measurement systems, days and kit lots). The results below were obtained without consideration of the sample dilution factor.

Sample	Mean value [ng/mL]	CV [%]
1	45.35	11.8
2	34.07	9.9

11.3 Accuracy – Trueness

The trueness states the closeness of the agreement between the result of a measurement and the true value of the measurand. Therefore, ox-LDL spikes with known concentrations were added to 2 different plasma samples. The results below were obtained without consideration of the sample dilution factor.

Sample [ng/mL]	Spike [ng/mL]	Expected [ng/mL]	Obtained [ng/mL]	Recovery [%]
31.4	10.0	41.4	39.6	95.65
	25.0	56.4	56.0	99.29
	30.0	61.4	59.3	96.58
22.5	12.5	35.0	31.1	88.86
	25.0	47.5	43.6	91.79
	50.0	72.5	73.8	101.79

11.4 Linearity

The linearity states the ability of a method to provide results proportional to the concentration of analyte in the test sample within a given range. This was assessed according to CLSI guideline EP06-A with a serial dilution of 2 different plasma and 1 serum sample.

For ox-LDL in serum and plasma, the method has been demonstrated to be linear from 6.27 to 273.00 ng/mL based on the standard curve without considering possibly used sample dilution factors, showing a non-linear behaviour of less than $\pm 20\%$ in this interval.

Sample	Dilution	Expected [ng/mL]	Obtained [ng/mL]	Recovery [%]
Plasma 1	1:15	200.50	200.50	100.00
	1:30	100.25	96.87	96.63
	1:60	50.12	49.95	99.65
	1:120	25.06	26.79	106.88
	1:240	12.53	14.55	116.09
	1:480	6.27	5.94	94.81
Plasma 2	1:100	273.00	273.00	100.00
	1:200	136.50	143.50	105.13
	1:400	68.25	71.25	104.40
	1:800	34.13	32.25	94.51
	1:1 600	17.06	16.88	98.90
	1:3 200	8.53	9.91	116.12
Serum	1:20	190.59	190.59	100.00
	1:40	91.37	95.30	95.88
	1:80	39.82	47.65	83.57
	1:160	18.86	23.82	79.16
	1:320	9.10	11.91	76.39

11.5 Analytical specificity

The specificity of the antibody was tested by measuring the cross-reactivity against a range of compounds with structural similarity to ox-LDL. There was no cross-reactivity

Substance tested	Concentration added	Concentration obtained [ng/mL]	Conclusion
HDL (High-Density-Lipoprotein)	10 500 ng/mL	< 4.130	< LoB
LDL (Low-Density-Lipoprotein)	14 400 ng/mL	< 4.130	< LoB
Albumin	800 ng/mL	< 4.130	< LoB
AOPP (Advanced oxidation protein products)	100 μ mol/L	< 4.130	< LoB

12 PRECAUTIONS

- All reagents in the kit package are for *in vitro* diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or Thimerosal as bactericides. Sodium azide and Thimerosal are toxic. Substrates for the enzymatic colour reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- The stop solution consists of diluted sulphuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped up immediately with copious quantities of water. Do not breath vapour and avoid inhalation.

13 TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay. Furthermore we recommend not assembling wells of different microtiter plates for analysis, even if they are of the same batch.
- Control samples should be analysed with each run.
- Reagents should not be used beyond the expiration date stated on the kit label.
- Substrate solution should remain colourless until use.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- Do not mix plugs and caps from different reagents.
- The assay should always be performed according to the enclosed manual.

14 GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- The guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. DRG can therefore not be held responsible for any damage resulting from incorrect use.

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15 REFERENCES / LITERATURE / LITERATUR

1. Corsi MM, Dogliotti G, Pedroni F, Ermetici F, Malavazos A, Ambrosi B, (2007) ADMA: a possible role in obese patients. Poster P173 of the 6th World Congress on Hyperhomocysteinemia, Saarbrücken, Germany, June 5-9, 2007, erschienen in CCLM 45(5)
2. Koubaa N, Nakbi A, Smaoui M, Abid N, Chaaba R, Abid M, Hammami M (2007) Hyperhomocysteinemia and elevated ox-LDL in Tunisian type 2 diabetic patients: Role of genetic and dietary factors. Clin Biochem 40(13-14):1007-14. Epub 2007 Jun 14
3. Licastro F, Dogliotti G, Goi G, Malavazos AE, Chiappelli M, Corsi MM (2007) Oxidated low-density lipoproteins (oxLDL) and peroxides in plasma of Down syndrome patients. Arch Gerontol Geriatr 44 Suppl 1:225-32
4. Pfützner A, Kost I, Löbig M, Knesovic M, Armbruster FP, Forst T (2005) Clinical Evaluation of a New ELISA Method for Determination of Oxidized LDL Particles - a Potential Marker for Arteriosclerotic Risk in Diabetes Mellitus. Abstract of the 5th Diabetes Technology Meeting, San Francisco, 10.-12. November 2005

SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitäts-kennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico In vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	Numero di Catalogo	Número de catálogo	Référence de catalogue
	Batch code *	Chargencode *	Codice del lotto	Codigo de lote	Numéro de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis *	Utilizzare prima del	Establa hasta	Utiliser jusque
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Caution *	Achtung *			
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Contenu
<i>Volume/No.</i>	Volume / No.	Volumen/Anzahl	Volume/Quantità	Volumen/Número	Volume/Quantité