## Instructions for Use

# Sample Preparation System (Stool)

IVD

REF EIA-5674

$\Sigma$  100

**DRG** 

DRG Instruments GmbH, Germany  
Frauenbergstraße. 18, D-35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

Distributed by:

**DRG** 

DRG International, Inc., USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.  
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.  
Si prega di usare la versione valida dell'inserito del pacco a disposizione con il kit.  
Por favor, se usa solo la version valida de la metodico técnico incluido aqui en el kit.**

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## PREPARATION OF A STOOL SAMPLE SOLUTION

### 1 DESIGN OF THE STOOL PREPARATION SYSTEM

(REF EIA-5674, unfilled stool sample system)

Stool sample tube

**Yellow:** Dipstick with click connector

**Blue:** Cone-shaped insert to strip off excessive stool

### 2 PREPARATIONS FOR STOOL SAMPLING

#### 2.1 Processing of faecal material

- a. The stool sample has to be thawed.  
For particularly heterogeneous samples, we recommend a mechanical homogenisation using an applicator, inoculation loop or similar device.
- b. Dilute the extraction buffer concentrate according to the Instructions for Use (IFU) of the kit used.

Then add the volume of extraction buffer stated in the IFU to the empty sample tube. The dilution factor of the stool suspension depends on the amount of stool sampled and the volume of the buffer (see dilution table).

#### 2.2 Dilution table

##### Dilution with 1.5 mL buffer:

Applied amount of stool:	15 mg
Buffer volume:	1.5 mL
Dilution factor:	<b><u>1:100</u></b>


##### Dilution with 0.75 mL buffer:

Applied amount of stool:	15 mg
Buffer volume:	0.75 mL
Dilution factor:	<b><u>1:50</u></b>


Please note: for fluid stool samples: 15 mg = 15 µL

### 3 APPLICATION OF THE STOOL SAMPLE TUBE

#### Steps 1-2

#### Steps 3 - 4.

1. **Remove dipstick**  
Open the upper yellow part of the cap and remove the attached dipstick. The blue cone-shaped insert remains on the tube!
2. **Obtain stool sample**  
Insert the dipstick into the stool at 3 different sites. Make sure that the notches at the lower part of the stick are covered completely with stool.
3. **Place dipstick back into the tube**  
Stick the dipstick with the attached stool only once back into the sample tube. Excess material will be stripped off by the cone-shaped insert leaving 15 mg of sample in the tube.
4. **Prepare stool sample suspension**  
Close the tube properly by turning until you hear the yellow shutter click. Then shake well until no stool sample remains in the notches and the stool is completely suspended in the extraction buffer. The suspension now contains 15 mg of stool (if the stool sample is fluid, pipet 15  $\mu$ L into the tube). For details about dilutions and volumes please refer to the dilution table.  
**Attention:** The stool dilution in the sample tube hardly changes color. Repeated insertion of sample into the tube will falsify the results!
5. **Apply stool sample solution**  
Open the complete cap of the tube including the lower blue cone-shaped insert. The stool-buffer solution is now ready for use according to the corresponding ELISA manual.

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