## Instructions for Use

# SARS-CoV-2 (RBD) Total Ab ELISA

IVD

REF EIA-6154

 96

**DRG Instruments GmbH**, Germany  
Frauenbergstraße 18, 35039 Marburg  
Phone: +49 (0)6421-1700 0, Fax: +49 (0)6421-1700 50  
Website: [www.drg-diagnostics.de](http://www.drg-diagnostics.de)  
E-mail: [drg@drg-diagnostics.de](mailto:drg@drg-diagnostics.de)

Distributed by:

**DRG International, Inc.**, USA  
841 Mountain Ave., Springfield, NJ 07081  
Phone: (973) 564-7555, Fax: (973) 564-7556  
Website: [www.drg-international.com](http://www.drg-international.com)  
E-mail: [corp@drg-international.com](mailto:corp@drg-international.com)

**Please use only the valid version of the Instructions for Use provided with the kit.  
Verwenden Sie nur die jeweils gültige, im Testkit enthaltene, Gebrauchsanweisung.  
Si prega di usare la versione valida delle istruzioni per l'uso a disposizione con il kit.  
Por favor, use sólo la versión válida de las instrucciones de uso que se suministran con el kit.  
Utilisez seulement la version valide des Instructions d'utilisation fournies avec le kit.**

**Table of Contents / Inhaltsverzeichnis / Tabella dei Contenuti / Tabla de Contenidos**

1	INTENDED USE.....	2	1	FINALIDAD PREVISTA .....	27
2	PRINCIPLE OF THE TEST .....	2	2	FUNDAMENTO DEL ENSAYO .....	27
3	WARNINGS AND PRECAUTIONS .....	3	3	PRECAUCIONES .....	27
4	REAGENTS.....	4	4	COMPONENTES DEL KIT .....	28
5	SPECIMEN COLLECTION AND PREPARATION .....	5	5	MUESTRAS .....	29
6	ASSAY PROCEDURE.....	5	6	PROCEDIMIENTO DE ENSAYO .....	29
7	RESULTS .....	7	7	RESULTADOS.....	31
8	QUALITY CONTROL.....	7	8	CONTROL DE CALIDAD.....	31
9	PERFORMANCE CHARACTERISTICS .....	8	9	CARACTERÍSTICAS DEL ENSAYO .....	32
10	LIMITATIONS OF USE.....	10	10	LIMITACIONES DE USO.....	32
11	LEGAL ASPECTS .....	10	11	ASPECTOS LEGALES.....	33
1	ZWECKBESTIMMUNG .....	11	12	REFERENCES / LITERATURE.....	34
2	TESTPRINZIP .....	11			
3	HINWEISE UND VORSICHTSMAßNAHMEN .....	12			
4	BESTANDTEILE DES KITS .....	13			
5	PROBENVORBEREITUNG .....	14			
6	TESTDURCHFÜHRUNG .....	14			
7	ERGEBNISSE .....	16			
8	QUALITÄTSKONTROLLE.....	16			
9	ASSAY-CHARAKTERISTIKA.....	17			
10	GRENZEN DES TESTS.....	18			
11	RECHTLICHE GRUNDLAGEN .....	19			
				SYMBOLS USED .....	35
1	DESTINAZIONE D'USO.....	20			
2	PRINCIPIO DEL TEST .....	20			
3	PRECAUZIONI .....	20			
4	COMPONENTI DEL KIT.....	21			
5	CAMPIONI .....	22			
6	ATTUAZIONE DEL TEST.....	22			
7	RISULTATI .....	24			
8	CONTROLLO QUALITÀ.....	24			
9	CARATTERISTICHE DEL TEST.....	25			
10	LIMITAZIONE DEL TEST.....	25			
11	ASPETTI LEGALI .....	26			

## 1 INTENDED USE

The **DRG SARS-CoV-2 (RBD) Total Ab ELISA** is an enzyme immunoassay for the **qualitative *in vitro* diagnostic** measurement of IgG/IgM/IgA antibodies to SARS-CoV-2 in serum or plasma (EDTA, lithium heparin or citrate plasma).

The DRG SARS-CoV-2 (RBD) total Ab ELISA is intended for use as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection. It can support the diagnosis of COVID-19 disease and supplement direct pathogen detection via real-time reverse transcription polymerase chain reaction (RT-PCR), or pulmonary CT-Scans. In addition, serology can help to gather epidemiological information on the disease prevalence.

### 1.1 Summary and Explanation

Coronavirus disease 2019 (COVID-19) is an infectious disease caused by the newly discovered Severe Acute Respiratory Syndrome Corona Virus 2 (SARS-CoV-2), which was first identified during an outbreak of respiratory illness cases in Wuhan City, China in December 2019 and has been declared a global pandemic by the WHO in March 2020<sup>1,2</sup>. SARS-CoV-2 is a positive-sense single-stranded RNA virus and belongs to the Betacoronavirus Genus, which also includes SARS-CoV (2003) and MERS (2012), but also other human coronaviruses like strains HKU1, OC43, NL63, and 229E which cause common cold with normally mild symptoms typically throughout winter month<sup>3</sup>.

SARS-CoV-2 is predominantly transmitted by droplet infection via coughing or sneezing and through close contact with infected patients. Smear infection through contaminated surfaces seem to play a minor role. Infection mainly occurs in the respiratory tract, but may also be caused through conjunctiva of the eyes. The incubation of COVID-19 ranges between 1-14 days, with the majority of cases manifesting within 4 to 6 days. The clinical manifestations of COVID-19 vary widely from asymptomatic courses to severe pneumonia with lung failure and death. However, in around 80 % of cases the symptoms are mild to moderate including fever, cough, respiratory problems, fatigue or temporarily loss of taste. Reported fatality rates increase with age, comorbidities (such as diabetes, cardiovascular or pulmonary disease, chronic liver disease, cancer) and genetic predisposition<sup>4</sup>.

The coronavirus spike (S) glycoprotein forms a homotrimeric class I viral fusion complex on the outer envelope of the virion that mediates the binding of the virus to its host cells. Subunit S1 of the spike protein contains a receptor binding domain (RBD) which interacts with the human ACE-2 receptor on the surface of the host cell and initiates the entry of the virus<sup>5-7</sup>. Antibodies which are directed against the RBD protein of SARS-CoV-2 are discussed to have neutralizing effects and to protect host cells from viral infection<sup>8,9</sup>.

Antibody response to SARS-CoV-2 infection can be detected as early as 3 days after symptom onset and reaches a peak after 2-3 weeks. Seroconversion for IgG and IgM occurs simultaneously or sequentially<sup>10,11</sup>. In consequence, the sensitivity of serological assay increases from the first week of symptom onset (< 40 %) to 100 % > 15 days after onset<sup>12</sup>. Duration of time antibodies are present post-infection is not well characterized.

Diagnosis of COVID-19 mainly relies on real-time reverse transcription polymerase chain reaction (RT-PCR) testing of respiratory specimens. However, PCR can't identify people who went through an infection, recovered, and cleared the virus from their bodies. Furthermore, chest X-rays, pulmonary computed tomography (CT) scans, and lung ultrasounds are important tools in the early diagnosis of pneumonia in patients with COVID-19<sup>4</sup>.

## 2 PRINCIPLE OF THE TEST

The DRG SARS-CoV-2 (RBD) total Ab ELISA is a solid phase enzyme immunoassay in the one-step antigen capture format.

Microtiter wells as a solid phase are coated with recombinant Receptor-Binding-Domain (RBD) of the Spike protein antigen of SARS-CoV-2.

Samples and controls are incubated in the coated wells together with enzyme conjugate (recombinant RBD protein coupled with horseradish peroxidase).

If anti-RBD antibodies are present in the sample, immobilized immune complexes are formed.

After a washing step to remove all unbound substances, the solid phase is incubated with substrate solution. The colorimetric reaction is stopped by addition of stop solution, and optical density (OD) of the resulting yellow product is measured.

The intensity of this color is directly proportional to the amount of anti-SARS-CoV-2 RBD antibodies in the patient specimen. Optical density (OD) at 450 nm is read using an ELISA microplate reader.

The presence of anti-SARS-CoV-2 RBD antibodies in an individual specimen is determined by comparing the OD values of the specimen to the OD values of the Cut-off Control.

**3 WARNINGS AND PRECAUTIONS**

1. This kit is for in vitro diagnostic use only. For professional use only.
2. All reagents of this test kit which contain human serum or plasma have been tested and confirmed negative for HIV I/II, HBsAg and HCV by FDA approved procedures. All reagents, however, should be treated as potential biohazards in use and for disposal.
3. Before starting the assay, read the instructions completely and carefully. Use the valid version of instructions for use provided with the kit. Be sure that everything is understood.
4. The microplate contains snap-off strips. Unused wells must be stored at 2 °C to 8 °C in the sealed foil pouch and used in the frame provided.
5. Pipetting of samples and reagents must be done as quickly as possible and in the same sequence for each step.
6. Use reservoirs only for single reagents. This especially applies to the substrate reservoirs. Using a reservoir for dispensing a substrate solution that had previously been used for the conjugate solution may turn solution coloured. Do not pour reagents back into vials as reagent contamination may occur.
7. Mix the contents of the microplate wells thoroughly to ensure good test results. Do not reuse microwells.
8. Do not let wells dry during assay; add reagents immediately after completing the rinsing steps.
9. Allow the reagents to reach room temperature (20 °C to 26 °C) before starting the test. Temperature will affect the optical density readings of the assay.
10. Never pipet by mouth and avoid contact of reagents and specimens with skin and mucous membranes.
11. Do not smoke, eat, drink or apply cosmetics in areas where specimens or kit reagents are handled.
12. Wear disposable latex gloves when handling specimens and reagents. Microbial contamination of reagents or specimens may give false results.
13. Handling should be done in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation.
14. Do not use reagents beyond expiry date as shown on the kit labels.
15. All indicated volumes have to be performed according to the protocol. Optimal test results are only obtained when using calibrated pipettes and microtiter plate readers.
16. Do not mix or use components from kits with different lot numbers. It is advised not to exchange wells of different plates even of the same lot. The kits may have been shipped or stored under different conditions and the binding characteristics of the plates may result slightly different.
17. Avoid contact with Stop Solution containing 0.5 M H<sub>2</sub>SO<sub>4</sub>. It may cause skin irritation and burns.
18. Some reagents contain ProClin 300, BND and/or MIT as preservatives. In case of contact with eyes or skin, flush immediately with water.
19. TMB substrate has an irritant effect on skin and mucosa. In case of possible contact, wash eyes with an abundant volume of water and skin with soap and abundant water. Wash contaminated objects before reusing them. If inhaled, take the person to open air.
20. Chemicals and prepared or used reagents have to be treated as hazardous waste according to the national biohazard safety guideline or regulation.
21. For information on hazardous substances included in the kit please refer to Safety Data Sheets. Safety Data Sheets for this product are available upon request directly from DRG.

## 4 REAGENTS

### 4.1 Reagents provided

1. **Microtiterwells**, 12 x 8 (break apart) strips, 96 wells;  
Wells coated with SARS-CoV-2 RBD antigen.
2. **Neg. Control**; 1 vial, 1.0 mL; ready to use;  
Contains non-mercury preservative
3. **Pos. Control**; 1 vial, 0.6 mL, ready to use;  
Contains non-mercury preservative
4. **Cut-Off Control**, 1 vial, 0.6 mL, ready to use;  
Contains non-mercury preservative.
5. **Enzyme Conjugate**, 1 vial, 12 mL, ready to use;  
colored red;  
SARS-CoV-2 RBD conjugated to horseradish peroxidase;  
Contains non-mercury preservative.
6. **Substrate Solution**, 1 vial, 14 mL, ready to use;  
Tetramethylbenzidine (TMB).
7. **Stop Solution**, 1 vial, 14 mL, ready to use;  
Contains 0.5 M H<sub>2</sub>SO<sub>4</sub>.  
Avoid contact with the stop solution. It may cause skin irritations and burns.
8. **Wash Solution**, 1 vial, 30 mL (40X concentrated);  
See "Reagent Preparation".

### 4.2 Materials required but not provided

- A calibrated microplate spectrophotometer (450 nm, with reference wavelength at 620 nm to 630 nm)
- Calibrated variable precision micropipettes
- Incubator 37 °C
- Manual or automatic equipment for rinsing wells
- Absorbent paper
- Distilled water
- Timer

### 4.3 Storage Conditions

When stored at 2 °C to 8 °C unopened reagents will retain reactivity until expiration date. Do not use reagents beyond this date. Opened reagents must be stored at 2 °C to 8 °C. Microtiter wells must be stored at 2 °C to 8 °C. Once the foil bag has been opened, care should be taken to close it tightly again.

### 4.4 Reagent Preparation

Bring all reagents and required number of strips to room temperature (20 °C to 26 °C) prior to use.

#### **Wash Solution**

Add distilled water to the 40X concentrated Wash Solution.

Dilute 30 mL of concentrated *Wash Solution* with 1170 mL distilled water to a final volume of 1200 mL.

*The diluted Wash Solution is stable for 1 week at room temperature.*

### 4.5 Disposal of the Kit

The disposal of the kit and all used materials/reagents must be performed according to the national regulations. Special information for this product is given in the Safety Data Sheet, section 13.

### 4.6 Damaged Test Kits

In case of any damage to the test kit or components, DRG must be informed in writing, at the latest one week after receiving the kit. Damaged single components should not be used for a test run. They have to be stored until a final solution has been found. After this, they should be disposed of according to the official regulations.

## 5 SPECIMEN COLLECTION AND PREPARATION

Serum or plasma (EDTA, lithium heparin or citrate plasma) can be used in this assay.

*Note:* Samples containing sodium azide should not be used in the assay.

In general, it should be avoided to use haemolytic, icteric, or lipaemic specimens. For further information refer to chapter “*Interfering Substances*”.

### 5.1 Specimen Collection

#### Serum:

Collect blood by venipuncture (e.g. Sarstedt Monovette for serum), allow to clot, and separate serum by centrifugation at room temperature. Do not centrifuge before complete clotting has occurred. Patients receiving anticoagulant therapy may require increased clotting time.

#### Plasma:

Whole blood should be collected into centrifuge tubes containing anticoagulant (e.g. Sarstedt Monovette with the appropriate plasma preparation) and centrifuged immediately after collection.

### 5.2 Specimen Storage and Preparation

Specimens should be capped and may be stored for up to 7 days at 2 °C to 8 °C prior to assaying.

Specimens stored for a longer time should be frozen only once at -20 °C prior to assay. Thawed samples should be inverted several times prior to testing.

### 5.3 Specimen Dilution

If a sample exceeds the measuring range of the photometer used, this sample can be diluted using *Neg. Control*.

## 6 ASSAY PROCEDURE

### 6.1 General Remarks

- All reagents and specimens must be allowed to come to room temperature before use. All reagents must be mixed without foaming.
- Once the test has been started, all steps should be completed without interruption.
- Use new disposal plastic pipette tips for each standard, control or sample in order to avoid cross contamination.
- Optical density is a function of the incubation time and temperature. Before starting the assay, it is recommended that all reagents are ready, caps removed, all needed wells secured in holder, etc. This will ensure equal elapsed time for each pipetting step without interruption.
- As a general rule the enzymatic reaction is linearly proportional to time and temperature.
- Close reagent vials tightly immediately after use to avoid evaporation and microbial contamination.

## 6.2 Test Procedure

1. Select the required number of microtiter strips or wells and insert them into the holder.

Please allocate at least:

- |         |                |                                |
|---------|----------------|--------------------------------|
| 1 well  | (e.g. A1)      | for the <i>Neg. Control</i>    |
| 2 wells | (e.g. B1 + C1) | for the <i>Cut-Off Control</i> |
| 1 well  | (e.g. D1)      | for the <i>Pos. Control</i>    |

It is left to the user to determine controls and patient samples in duplicate.

2. Dispense  
**20 µL** of *Neg. Control* into well A1  
**20 µL** of *Cut-Off Control* into well B1 + C1  
**20 µL** of *Pos. Control* into well D1 and  
**20 µL** of each sample with new disposable tips into appropriate wells
3. Dispense **80 µL Enzyme Conjugate** into each well.  
Thoroughly mix for 10 seconds. It is important to have a complete mixing in this step.
4. Incubate for **60 minutes at 37 °C**.  
**Do not use a cover foil to cover the wells during incubation!**
5. Rinse the wells **3 times** with **400 µL** diluted *Wash Solution* per well, if a plate washer is used.  
**- OR -**  
Briskly shake out the contents of the wells.  
Rinse the wells **3 times** with **300 µL** diluted *Wash Solution* per well for manual washing. Strike the wells sharply on absorbent paper to remove residual droplets.  
**Important note:**  
The sensitivity and precision of this assay is markedly influenced by the correct performance of the washing procedure!
6. Add **100 µL** of **Substrate Solution** to each well.
7. Incubate for **15 minutes** at room temperature (20 °C to 26 °C).  
**Do not use a cover foil to cover the wells during incubation!**
8. Stop the enzymatic reaction by adding **50 µL** of **Stop Solution** to each well.  
**Note:** Highly positive patient samples can cause dark precipitates of the chromogen!
9. Measure the optical density of the solution in each well at **450 nm (reading) and at 620 nm to 630 nm (background subtraction, recommended)** with a microtiter plate reader.  
It is recommended that the wells be read **within 10 minutes** after adding the *Stop Solution*.

## 7 RESULTS

Where applicable calculate the mean OD values of all duplicates.

### 7.1 Validation of the Test Run

In order for an assay run to be considered valid, the following criteria must be met:

<b>Neg. Control in A1:</b>	OD < 0.15
<b>Cut-Off Control in B1 + C1:</b>	0.15 < OD < 0.50
<b>Pos. Control in D1:</b>	OD > 0.50

### 7.2 Interpretation of qualitative Results

$$\text{NEGATIVE} \quad \frac{\text{Mean OD}_{\text{Patient}}}{\text{Mean OD}_{\text{Cut-Off Control}}} < 0.9$$

$$\text{GREY ZONE} \quad 0.9 < \frac{\text{Mean OD}_{\text{Patient}}}{\text{Mean OD}_{\text{Cut-Off Control}}} < 1.1$$

Repeat test after e.g. 2 weeks with new patient samples.

Results in the second test again in the grey zone → NEGATIVE

$$\text{POSITIVE} \quad \frac{\text{Mean OD}_{\text{Patient}}}{\text{Mean OD}_{\text{Cut-Off Control}}} > 1.1$$

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## 8 QUALITY CONTROL

It is recommended to use control samples according to state and federal regulations. The use of control samples is advised to assure the day to day validity of results. Use controls at both normal and pathological levels.

If the results of the assay do not fit to the established acceptable ranges of control materials patient results should be considered invalid.

In this case, please check the following technical areas: Pipetting and timing devices; photometer, expiration dates of reagents, storage and incubation conditions, aspiration and washing methods.

After checking the above mentioned items without finding any error contact your distributor or DRG directly.

## 9 PERFORMANCE CHARACTERISTICS

### 9.1 Specificity of Antigen (Cross-Reactivity)

The following samples were tested for cross-reactivity of the assay:

Sample positive to	n
Human Coronavirus NL63 Nucleoprotein *	5
Human Coronavirus 229E Nucleoprotein IgG *	9
MERS IgG *	3
Coronavirus HKu1 IgG *	13
Coronavirus OC43 IgG *	12
Human Coronavirus 229E/NL63 Nucleoprotein IgG *	5
Human Coronavirus NL63 Nucleoprotein and Parainfluenza 1/2/3 IgG *	1
Human Coronavirus 229E/NL63 Nucleoprotein IgG and Parainfluenza 1/2/3 IgG *	1
Adenovirus IgG	10
Bordetella pertussis IgG	6
Chlamydia pneumonia IgG	14
Chlamydia trachomatis IgG	6

Sample positive to	n
EBV IgG	6
Enterovirus IgA	7
Enterovirus IgG	7
Enterovirus IgM	4
HCV IgG	28
Influenza A IgG	5
Influenza B IgG	5
Legionella pneumophila IgM	2
Mycoplasma pneumonia IgG	4
Parainfluenza 1/2/3 IgG	16
Parvovirus IgG	5
RSV IgG	10
Rubella IgG	3

\*marked samples were used as Coronavirus panel for evaluation of diagnostic specificity.

Sample positive to	n
Rheumatoid factor IgM	6
Rheumatoid factor IgA	4
different ANA	16
HAMA	6

Sample	n
Pregnant women	14

The DRG SARS-CoV-2 (RBD) Total Ab ELISA shows no cross-reactivity to the evaluated antibodies.

### 9.2 Analytical Sensitivity

The analytical sensitivity of the DRG ELISA was calculated by adding 2 standard deviations to the mean OD of 20 replicate analyses of the *Neg. Control* and was found to be 1.350 DU (OD = 0.039).

### 9.3 Diagnostic Specificity, Sensitivity and Cut-Off

ROC-Analysis is used for the determination.

Using the Cut-Off 0.33 OD, the following specificities were found.

Cohort	n	Non-reactive	Reactive	Specificity
Diagnostic Routine panel	183	183	0	100.0 %
Blood Donors panel	359	357	2	99.4 %
Coronavirus panel *	49	49	0	100.0 %
<b>Total</b>	<b>591</b>	<b>589</b>	<b>2</b>	<b>99.7 %</b>

\* See chapter 9.1 Specificity of Antigens

Using the Cut-Off 0.33 OD, the following sensitivity was found.

Days post PCR confirmation	n	Non-reactive	Reactive	Sensitivity
> 14	86	0	86	100 %

Please note that the sensitivity decreases for samples collected earlier than 14 days after infection with the coronavirus (<sup>9-12</sup>).

## 9.4 Comparison Studies

The DRG SARS-CoV-2 (RBD) Total Ab ELISA was compared with CE-certified Roche Elecsys Anti-SARS-CoV-2 (ECLIA).

n = 35		Roche Elecsys Anti-SARS-CoV-2 (ECLIA)	
		pos.	neg.
DRG ELISA	pos.	35	0
	neg.	0	0

Agreement: 100.0 %

## 9.5 Reproducibility

### 9.5.1 Intra-Assay

The intra-assay (within-run) precision of the DRG ELISA was determined by 20x measurements of 12 samples covering the measuring range of the ELISA. For each level of diagnosis, the mean CV of three samples was calculated (3 samples x 20 measurements = 60 determinations).

Sample	n	OD range	CV (%)
Negative samples	60	0.017 – 0.019	9.1
Low positive samples	60	0.388 – 0.703	4.6
Positive samples	60	0.935 – 1.524	4.1
High positive samples	60	2.517 – 3.627	3.9

### 9.5.2 Inter-Assay

Inter-assay precision for samples was determined for 4 patient samples covering the measuring range in 3 independent runs on 3 days with 10 determinations. CV was calculated from 30 determinations.


Sample	n	Mean OD	Mean Conc. [DU]	CV (%) Conc.
1	30	0.020	0.8	8.9
2	30	0.287	11.8	12.0
3	30	1.228	50.2	8.4
4	30	3.333	136.3	8.5

Inter-assay precision for controls was determined by measuring each control 24 times in 12 independent runs. CV was calculated from 20 determinations.

Control	n	Mean OD	CV (%)
Neg. Control	24	0.027	14.8
Cut-Off Control	24	0.339	6.5
Pos. Control	24	0.978	5.7

## 9.6 Linearity

Four samples (serum) containing different amounts of analyte were serially diluted 1:2 up to 1:16 with *Neg. Control* and assayed with the DRG ELISA. The percentage recovery was calculated by comparing the expected and measured values for the analyte.

## 10 LIMITATIONS OF USE

Reliable and reproducible results will be obtained when the assay procedure is performed with a complete understanding of the package insert instruction and with adherence to good laboratory practice.

Any improper handling of samples or modification of this test might influence the results.

Bacterial or fungal contamination of patient samples or reagents, or cross-contamination between reagents may cause erroneous results.

This test is only for qualitative detection. Test results should not be the sole basis for clinical diagnosis and treatment. The confirmation of infection with SARS-CoV 2 must be combined with the patient's clinical signs in conjunction to other tests.

In the first week of the onset of the infection with SARS-CoV 2, patients results may be negative for antibodies. The sensitivity increases over the first two weeks and reaches a plateau after two to three weeks.

In addition, patients with low immunity or other diseases that affect immune function, failure of important systemic organs, and use of drugs that suppress immune function can also lead to negative results of SARS-CoV 2 antibodies.

### 10.1 Interfering Substances

Haemoglobin (up to 4 mg/mL), Bilirubin (up to 0.5 mg/mL) and Triglyceride (up to 7.5 mg/mL) have no influence on the assay results.

### 10.2 High-Dose-Hook Effect

"High Dose Hook Effect" is not detected for results up to 150 DU.

## 11 LEGAL ASPECTS

### 11.1 Reliability of Results

The test must be performed exactly as per the manufacturer's instructions for use. Moreover the user must strictly adhere to the rules of GLP (Good Laboratory Practice) or other applicable national standards and/or laws. This is especially relevant for the use of control reagents. It is important to always include, within the test procedure, a sufficient number of controls for validating the accuracy and precision of the test.

The test results are valid only if all controls are within the specified ranges and if all other test parameters are also within the given assay specifications. In case of any doubt or concern please contact DRG.

### 11.2 Therapeutic Consequences

Therapeutic consequences should never be based on laboratory results alone even if all test results are in agreement with the items as stated under point 11.1. Any laboratory result is only a part of the total clinical picture of a patient.

Only in cases where the laboratory results are in acceptable agreement with the overall clinical picture of the patient should therapeutic consequences be derived.

The test result itself should never be the sole determinant for deriving any therapeutic consequences.

### 11.3 Liability

Any modification of the test kit and/or exchange or mixture of any components of different lots from one test kit to another could negatively affect the intended results and validity of the overall test. Such modification and/or exchanges invalidate any claim for replacement.

Claims submitted due to customer misinterpretation of laboratory results subject to point 11.2 are also invalid.

Regardless, in the event of any claim, the manufacturer's liability is not to exceed the value of the test kit. Any damage caused to the test kit during transportation is not subject to the liability of the manufacturer.

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
**12 REFERENCES / LITERATURE**

1. Ji T. et al. Detection of Covid-19: A review of the current literature and future perspectives. *Biosensors and Bioelectronics*. 2020 166; 112455.
2. World Health Organization. Coronavirus disease 2019 (COVID-19) situation report-85, Data as received by WHO from national authorities by 10:00 CET, 14 April 2020. Accessed April 14, 2020. <https://www.who.int/docs/default-source/coronaviruse/situation-reports/20200414-sitrep-85-covid-19>.
3. Zhong N et al. Epidemiology and cause of severe acute respiratory syndrome (SARS) in Guangdong People's Republic of China in February 2003. *Lancet*. 2003 362:1353-1358
4. da Silva SJR et al. Clinical and laboratory diagnosis of SARS-CoV-2, the virus causing COVID-19. *ACS Inf Dis*. August 2020; <https://dx.doi.org/10.1021/acsinfecdis.Oc00274>.
5. Brielle ES, Schneidman-Duhovny D, and Linial M. The SARS-CoV-2 exerts a distinctive strategy for interacting with the ACE2 human receptor. *Viruses*. 2020 12(5):497.
6. Gui M et al. Cryo-electron microscopy structures of the SARS-CoV spike glycoprotein reveal a prerequisite conformational state for receptor binding. *Cell Research*. 2017 27:119-129.
7. Turonova B et al. In situ structural analysis of SARS-CoV-2 spike reveals flexibility mediated by three hinges. *Science*. 2020 10.1126/science.abd.5223.
8. Premkumar L et al. The receptor binding domain of the viral spike protein is an immunodominant and highly specific target of antibodies in SARS-CoV-2 patients. *Science Immunology* 2020 10.1126/sciimmunol.abc8413.
9. Robbiani DF et al. Convergent antibody responses to SARS-CoV-2 in convalescent individuals. *Nature*. 2020 584:437-456.
10. Kontou PI et al. Antibody tests in detecting SARS-CoV-2 infections: a meta-analysis. *Diagnostics*. 2020 10:319-34.
11. Long QX et al. Antibody responses to SARS-CoV-2 in patients with COVID-19. *Nature Medicine*. 2020 26:845-848.
12. Zhao J et al. Antibody responses to SARS-CoV-2 in patients of novel coronavirus disease 2019. *Clin Inf Dis*. 2020 Accepted manuscript doi:10.1093/cid/ciaa344.

## SYMBOLS USED

Symbol	English	Deutsch	Italiano	Español	Français
	European Conformity	CE-Konformitätskennzeichnung	Conformità europea	Conformidad europea	Conformité normes européennes
	Consult instructions for use *	Gebrauchsanweisung beachten *	Consultare le istruzioni per l'uso	Consulte las instrucciones de uso	Consulter les instructions d'utilisation
	<i>In vitro</i> diagnostic medical device *	<i>In-vitro</i> -Diagnostikum *	Dispositivo medico-diagnostico in vitro	Producto sanitario para diagnóstico in vitro	Dispositif médical de diagnostic in vitro
	Catalogue number *	Artikelnummer *	No. di Cat.	No de catálogo	Référence
	Batch code *	Fertigungslosnummer, Charge *	Lotto no	Número de lote	No. de lot
	Contains sufficient for <n> tests *	Ausreichend für <n> Prüfungen *	Contenuto sufficiente per "n" saggi	Contenido suficiente para <n> ensayos	Contenu suffisant pour "n" tests
	Temperature limit *	Temperaturbegrenzung *	Temperatura di conservazione	Temperatura de conservación	Température de conservation
	Use-by date *	Verwendbar bis *	Data di scadenza	Fecha de caducidad	Date limite d'utilisation
	Manufacturer *	Hersteller *	Fabbricante	Fabricante	Fabricant
	Biological risks *	Biologische Risiken *	Rischi biologici	Riesgos biológicos	Risques biologiques
	Caution *	Achtung *	Attenzione	Precaución	Attention
	For research use only	Nur für Forschungszwecke	Solo a scopo di ricerca	Sólo para uso en investigación	Seulement dans le cadre de recherches
<i>Distributed by</i>	Distributed by	Vertreiber	Distributore	Distribuidor	Distributeur
<i>Content</i>	Content	Inhalt	Contenuto	Contenido	Conditionnement
<i>Volume/No.</i>	Volume / No.	Volumen / Anzahl	Volume / Quantità	Volumen / Número	Volume / Quantité
<i>Microtiterwells</i>	Microtiterwells	Mikrotiterwells	Micropozzetti	Placas multipocillo	Microplaques
<i>Enzyme Conjugate</i>	Enzyme Conjugate	Enzymkonjugat	Tracciante enzimatico	Conjugado enzimático	Conjugué enzymatique
<i>Substrate Solution</i>	Substrate Solution	Substratlösung	Soluzione di substrato	Solución de sustrato	Solution substrat
<i>Stop Solution</i>	Stop Solution	Stopplösung	Soluzione d' arresto	Solución de parada	Solution d'arrêt
<i>Standard</i>	Standard	Standard	Standard	Estándar	Standard
<i>Control</i>	Control	Kontrolle	Controllo	Control	Contrôle
<i>Pos. Control</i>	Positive Control	Positive Kontrolle	Controllo positivo	Control positivo	Positif Contrôle
<i>Neg. Control</i>	Negative Control	Negative Kontrolle	Controllo negativo	Control negativo	Négatif Contrôle
<i>Cut-off Control</i>	Cut-off Control	Grenzwert-Kontrolle	Controllo valore limite	Control valor limite	Valeur limite Contrôle
<i>Wash Solution</i>	Wash Solution	Waschlösung	Soluzione di lavaggio	Solución de lavado	Solution de lavage
<i>Sample Diluent</i>	Sample Diluent	Probenverdünnungsmedium	Diluyente dei campioni	Solución para dilución de la muestra	Solution pour dilution de l'échantillon