Lyophilized Serum Calibrator

HDL CALIBRATOR (HDL Cholesterol Calibrator)

Lyophilized human based calibrator for the use in tests for quantitative in vitro determination of high density lipoproteins cholesterol (HDL-C) on photometric systems.

| |
|------------|
| REF |
|------------|

| |
|--------------|
| Cont. |
|--------------|

| | | |
|------------------|-----------------|-----------------------------|
| F03710SV | 1 x 3 ml | single vial |
| F03710VSV | 1 x 3 ml | unlabeled single vial (OEM) |

Additionally offered:

| | | |
|--------|------------|--------------------------|
| F03100 | 5 x 100 ml | Chol. HDL Direct Reagent |
| F03115 | 5 x 50 ml | Chol. HDL Direct Reagent |

CALIBRATOR COMPOSITION

HDL-Cholesterol in human serum based matrix
 non reactive stabilizers
 preservatives

CALIBRATOR PREPARATION

1. Open the vial very carefully, avoiding any loss of the lyophilized material.
2. Add exactly 3 ml of dist. water (inaccurate reconstitution of the control and error in assay technique can cause erroneous results).
3. Close vial carefully and let stand for 30 min.
4. Dissolve contents completely by swirling gently, avoiding the formation of foam.
5. Do not shake !

CALIBRATOR STABILITY AND STORAGE

Storage: at 2 – 8°C
 Stability: until indicated date of expiration
 Stability after reconstitution:
 7 days at 2 – 10 °C
 30 days at –20 °C
CLOSE IMMEDIATELY AFTER USE
FREEZE ONLY ONCE!

Do not use the product if there is visible evidence of microbial growth in the vial.
 Improper handling and/or storage of the calibrator can affect results.

WARNINGS AND PRECAUTIONS

1. Warning.
 H317: May cause an allergic skin reaction.
 P261: Avoid breathing dust/fume/gas/mist/vapours/spray.
 P280: Wear protective gloves/protective clothing/eye protection/face protection.
 P303+P361+P353: In on skin (or hair): Remove/take off immediately all contaminated clothing. Rinse skin with water/shower.
 P333+P313: If skin irritation or rash occurs: Get medical advice/attention.
2. Each individual donation intended for use in manufacture of DIALAB HDL calibrator was tested for HbsAg, anti-HCV and anti-HIV 1 and HIV 2 by FDA required tests. Since no test method can assure that products derived from human blood do not contain HIV-1/2 and Hepatitis B and Hepatitis C virus, this material and all patient samples should be handled as though capable of transmitting infectious diseases.

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