## Lyophilized Serum Calibrator

### HDL/LDL CALIBRATOR (HDL/LDL-Cholesterol Calibrator)

Lyophilized human based calibrator for the use in tests for quantitative in vitro determination of high density lipoproteins cholesterol (HDL-C) and low density lipoproteins cholesterol (LDL-C) on photometric systems.

**REF**

F14715SV

**Cont.**

1 x 3 ml single vial

#### CALIBRATOR COMPOSITION

HDL- and LDL-Cholesterol in human serum

#### CALIBRATOR PREPARATION

1. Open the vial very carefully, avoiding any loss of the lyophilized material.
2. Add exactly 3 ml of deionised or distilled water (inaccurate reconstitution of the control and error in assay technique can cause erroneous results).
3. Close vial carefully and dissolve contents completely by swirling gently, avoiding the formation of foam. Do not shake!

#### CALIBRATOR STABILITY AND STORAGE

Storage: at 2 – 10 °C

Stability: until date of expiration stated on the label

Stability after reconstitution:

7 days at 2 – 10 °C

Close immediately after use.

Do not use the product if there is visible evidence of microbial growth in the vial.

Improper handling and/or storage of the calibrator can affect results.

#### WARNINGS AND PRECAUTIONS

1. If the calibrators come in contact with the mouth, eyes or skin, wash off immediately with a large amount of water. Consult a physician if necessary.
2. This product is prepared from human sera which was found to be negative for HBsAg, anti-HIV 1&2 antibodies and anti-HCV antibody. Because no test can offer complete assurance that products derived from human sera will not transmit infectious agents, it is recommended that this product should be handled with the same precautions used for patient specimens.
3. Be careful not to cut yourself with the aluminium cap when removing it from the vial.

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