Mythic 18-22 FLUSH - CLEANER

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Data wydania: 10. 2022.

Mythic 18-22 FLUSH - CLEANER

Reagent name	Variant	Cat. No
Mythic 18-22 FLUSH - CLEANER	250 ml	HM18-011-025

INTENDED USE

Mythic 18-22 FLUSH - CLEANER is intended for periodic and emergency cleaning of the hematology analyzers' measuring system. It is only for *in vitro* diagnostics, for healthcare professional users.

OPERATION PRINCIPLE

Mythic 18-22 FLUSH - CLEANER is a concentrated solution of hypochlorite intended for intensive oxidative cleaning, rinsing and washing of the hematology analyzers' measuring system. Mythic 18-22 FLUSH - CLEANER is used to remove sediments formed from tested samples.

CONCENTRATION OF THE ACTIVE INGREDIENTS IN THE REAGENT

active chlorine	3,5 %
sodium hydroxide	0,2 %
stabilizer	0,5 %
sodium chloride	

REAGENT STABILITY

The reagent must be used before the expiry date indicated on the label.
Once opened, the reagent is stable 9 weeks at room temperature.

WARNING AND NOTES

- Protect from direct sunlight!
- Do not use reagent beyond expiry date printed on the label.
- Do not combine the remaining reagent with the reagent from a new package.
- Use personal protective equipment to prevent contact with samples, reagents and controls.
- Avoid discharge to sewage system or to environment.
- Please refer to the MSDS for detailed information concerning safe storage and use of the product.
- Mythic 18-22 FLUSH - CLEANER meets the criteria for classification in accordance with Regulation (EC) No 1272/2008.

Ingredients:

Mythic 18-22 FLUSH - CLEANER contains sodium hypochlorite.

Danger

- H314 Causes severe skin burns and eye damage.
- H318 Causes serious eye damage.
- H400 Very toxic to aquatic life.
- H411 Toxic to aquatic life with long lasting effects.
- H290 May be corrosive to metals.

P280 Wear protective gloves, protective clothing, eye protection or face protection.
P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER or doctor.
P273 Avoid release to the environment.
P391 Collect spillage.
P234 Keep only in original packaging.

ASSAY PROCEDURE

Due to differences in cleaning procedures depending on the analyzer model, follow the instructions in the device user's manual and/or Service recommendations.

WASTE MANAGEMENT ¹

After use, the reagents should be handled as potentially infectious and disposed of in accordance with local legal requirements.

- Reagents residues: 18 01 06*
- Empty packages: 15 01 10*
- Wastewater from the analyzer: 18 01 03*

INCIDENTS ²

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer (website address: incidents@cormay.pl) and the competent authority of the Member State in which the user and/or the patient is established.

- Serious incident** means any incident that directly or indirectly led, might have led or might lead to any of the following:
- the death of a patient, user or other person,
 - the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
 - a serious public health threat.

LITERATURE

- European Commission notice on technical guidance on the classification of waste (2018/C 124/01) of 9 April 2018.
- Regulation (EU) 2017/746 of the European Parliament and of the Council of 5 April 2017 on *in vitro* diagnostic medical devices.

LIST OF CHANGES

Previous version: 05	Current version: 06
Sections updated: <i>WARNING AND NOTES</i>	

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