# “DIAQUICK” HSV 1/2 IgM Cassette

for serum and plasma samples

**REF**

**Content**

- J15040CE** - 25 cassettes individually packed (25 x Ref. No: J15040B)  
 - 25 disposable pipettes  
 - 1 vial buffer sufficient for 25 tests  
 - 1 package insert

For professional in vitro diagnostic use only

**GENERAL INFORMATION**

<b>Method</b>	sandwich type immunochromatographic assay
<b>Shelf life</b>	24 months from date of production
<b>Storage</b>	2 – 30 °C
<b>Sample</b>	human serum or plasma
<b>Results</b>	after 15 minutes

**INTENDED USE**

The “DIAQUICK” HSV 1/2 IgM Cassette (serum, plasma) is a lateral flow chromatographic immunoassay for the qualitative detection of IgM antibodies to HSV 1/2 in serum or plasma to aid in the diagnosis of HSV 1/2 infection.

**SUMMARY**

Herpes Simplex Virus (HSV) is a common pathogen of humans. There are two distinct types of HSV: type 1 and 2. HSV-1 is usually associated with infection in the oropharyngeal area and eyes, while HSV-2 causes mostly genital and neonatal infections.<sup>1,2</sup> However, the tissue specificity is not absolute.<sup>3</sup> Both HSV-1 and HSV-2 can be found at either site of infection without being specific for one region, i.e., either genital or oropharyngeal. Neonatal HSV infection may remain localized or become disseminated. Localized infection may involve skin, eyes, mouth or the central nervous system. Infants with neonatal HSV bear the risk of developing severe neurological or ocular sequelae. The detection of antibodies to HSV has been acknowledged to be a screening method for HSV infections. The methods for such antibody detection include complement fixation, indirect immunofluorescent antibody detection, plaque neutralization and ELISA. Antibodies of the IgM class are produced during the first 2-3 weeks of HSV infection and exist only transiently in most patients. The detection of IgM antibodies against HSV helps to distinguish between primary and recurrent infections.

The detection of anti-HSV 1/2 IgM antibodies enables the effective diagnosis of acute or recent HSV 1/2 infection. The “DIAQUICK” HSV 1/2 IgM Cassette is a rapid chromatographic immunoassay for the qualitative detection of IgM antibodies to HSV 1/2 in serum or plasma specimens.

**TEST PRINCIPLE**

The “DIAQUICK” HSV 1/2 IgM Cassette is a qualitative, lateral flow immunoassay for the detection of IgM antibodies to HSV-1 and/or HSV-2 in serum or plasma specimens. In this test, antigens of HSV 1/2 are coated in the test line region of the test. During testing, the serum or plasma specimen reacts with goat anti human IgM coated particles on the test strip. The mixture then migrates forward on the membrane by capillary action and reacts with the HSV 1/2 specific antigens on the membrane in the test line region. The presence of a coloured line in the test line region indicates a positive result for HSV 1/2 infection, while its absence indicates a negative result for the infection. To serve as a procedural control, a coloured line will always appear in the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

**REAGENTS**

The test contains goat anti-human IgM, and HSV 1/2 antigen. A streptavidin-IgG is employed in the control line system.

**STORAGE AND STABILITY**

Store as packaged in the sealed pouch at room temperature or refrigerated (2 – 30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

**PRECAUTIONS**

- For in vitro diagnostic use only. Do not use after the expiration date.
- Do not smoke, drink, or eat in areas where specimens or kits reagents are handled.
- Wear protective clothing such as laboratory coats, disposable gloves and eye protection, when specimens are being tested.
- Humidity and temperature can adversely affect results.
- The used test should be discarded according to local regulations.

**MATERIAL REQUIRED BUT NOT PROVIDED**

- Specimen collection container
- Centrifuge (for plasma only)
- Timer

**SPECIMEN COLLECTION AND PREPERATION**


- The “DIAQUICK” HSV 1/2 IgM Cassette can be performed using either serum or plasma specimens.
- Separate the serum or plasma from blood as soon as possible to avoid haemolysis. Only clear, non-haemolysed specimens can be used.
- Testing should be performed immediately after the specimens have been collected. Do not leave the specimens at room temperature for prolonged periods. Specimens may be stored at 2 – 8 °C for up to 3 days. For long-term storage, specimens should be kept below -20 °C.
- Bring specimens to room temperature prior to testing. Frozen specimens must

be completely thawed and mixed well prior to testing. Specimens should not be frozen and thawed repeatedly.

- If specimens are to be shipped, they should be packed in compliance with local regulations for the transportation of etiologic agents.

**ASSAY PROCEDURE**

- Bring the pouch to room temperature before opening it. Remove the test cassette from the sealed pouch and use it as soon as possible. Best results will be obtained if the assay is performed within one hour.
- Place the test cassette on a clean and level surface. Hold the dropper vertically and transfer 1 drop of serum or plasma (approx. 10 µL) and 2 drops of buffer (approx. 80 µL) to the specimen well of the test cassette. Avoid trapping air bubbles in the specimen well. See the illustration below.
- Wait for the coloured line(s) to appear. The result should be read at 15 minutes. Do not interpret results after 20 minutes.

**INTERPRETATION OF RESULTS**

(Please refer to the illustration above)

**POSITIVE:** Two coloured lines appear. One coloured line should always appear in the control line region (C) and another line should be in the test line region.

**\*NOTE:** The intensity of the colour in the test line regions may vary depending on the concentration of HSV1/2 antibodies present in the specimen. Therefore, any shade of colour in the test line region should be considered positive.

**NEGATIVE:** One coloured line appears in the control line region (C). No line appears in the test line regions.

**INVALID:** Control line fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

**QUALITY CONTROL**

Internal procedural controls are included in the test. A coloured line appearing in the control region (C) is an internal positive procedural control. It confirms sufficient specimen volume and correct procedural technique. Control standards are not supplied with this kit; however, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance.

**EXPECTED VALUES**

The “DIAQUICK” HSV 1/2 IgM Cassette has been compared with leading commercial EIA HSV 1/2 tests, demonstrating an overall accuracy of 98.9%.

**LIMITATIONS**

- The “DIAQUICK” HSV 1/2 IgM Cassette is for in vitro diagnostic use only. This test should be used for the detection of IgM antibodies to HSV-1 and/or HSV-2 in serum or plasma specimens. Neither the quantitative value nor the rate of increase in the concentration of IgM antibodies to HSV-1 and/or HSV-2 can be determined by this qualitative test.
- The “DIAQUICK” HSV 1/2 IgM Cassette will only indicate the presence of IgM antibodies to HSV-1 and/or HSV-2 in the specimen and should not be used as the sole criteria for the diagnosis of HSV 1/2 infection.
- As with all diagnostic tests, all results must be considered with other clinical information available to the physician.
- If the test result is negative and clinical symptoms persist, additional follow-up testing using other clinical methods is suggested. A negative result at any time does not preclude the possibility of HSV 1/2 infection.

**PERFORMANCE CHARACTERISTICS**

**Sensitivity and Specificity**

The “DIAQUICK” HSV 1/2 IgM Cassette was compared with leading commercial EIA HSV 1/2 tests; the results show that “DIAQUICK” HSV 1/2 IgM Cassette has a high sensitivity and specificity.

Method	HSV 1/2 EIA (IgM)		Total Results
	Positive	Negative	
„DIAQUICK“ HSV 1/2 IgM Cassette	Positive	1	21
	Negative	305	307
	Total Results	22	306

Sensitivity = 90.9% (70.8 – 98.9%)\*  
 Accuracy = 99.1% (97.4 – 99.8%)\*

Specificity = 99.7% (98.2 – 100%)\*  
 \* 95% Confidence Interval

**Precision**

**Intra-Assay**

Within-run precision has been determined by using 10 replicates of three specimens: a negative, a low positive, and a high positive. The negative, low positive, and high positive values were correctly identified >99% of the time.

**Inter-Assay**

Between-run precision has been determined by 10 independent assays on the same three specimens: a negative, a low positive, and a high positive. Three different lots of the “DIAQUICK” HSV 1/2 IgM Cassette have been tested over a 3-day period using negative, low positive, and high positive specimens. The specimens were correctly identified >99% of the time.

#### Cross-reactivity

The "DIAQUICK" HSV 1/2 IgM Cassette has been tested for HAV, HBV, HCV, HIV, RF, Syphilis, H. Pylori, Rubella, TOXO, CMV positive specimens. The results showed no cross-reactivity.

#### Interfering Substances

The "DIAQUICK" HSV 1/2 IgM Cassette has been tested and no interference was observed in specimens containing 110 mg/mL human albumin, 1 mg/mL bilirubin, 10 mg/mL haemoglobin, 0.2 mg/mL cholesterol and 15 mg/mL triglycerides. The following compounds have also been tested using the "DIAQUICK" HSV 1/2 IgM Cassette and no interference was observed.

Acetaminophen	Caffeine	EDTA
Acetylsalicylic Acid	Gentisic Acid	Ethanol
Ascorbic Acid	Phenylpropanolamine	Glucose
Bilirubin: 1g/dL	Salicylic Acid	Phenothiazine

#### REFERENCES

1. Bosma T.J. Corbett KM, O'Shea S. Banatvala .W. Best JM. PCR for detection of rubella virus RNA in clinical samples. J Clin Microbiol 1995; 33:1075-9 Thomas HIJ, Islogran-Capner P, Connor NS. Adaptation of commercial rubella specific
2. IgG kit to assess specific IgG avidity. Serodiagn Immunother Infect Disease 1993.1:13-6.
3. Mellinger AK, Cragan ID. Atkinson WL et al. High incidence of congenital rubella syndrome after a rubella outbreak. PediatrInfect Dis J 1995:14:573-5

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**BIBLIOGRAPHIE**

1. Bosma T.J, Corbett KM, O'Shea S, Banatvala .W, Best JM. PCR for detection of rubella virus RNA in clinical samples. J Clin Microbiol 1995; 33:1075-9  
 Thomas HJ, Islogran-Capner P, Connor NS. Adaptation of commercial rubella specific
2. IgG kit to assess specific IgG avidity. Serodiagn Immunother Infect Disease 1993.1:13-6.
3. Mellinger AK, Cragan ID, Atkinson WL et al. High incidence of congenital rubella syndrome after a rubella outbreak. PediatrInfect Dis J 1995;14:573-5